


## Cytochrome C Tryptic Digest – pH 2.5, XTerra RP<sub>18</sub>

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract

This application brief highlights the analysis of cytochrome C using XTerra RP<sub>18</sub> columns.

---

## Introduction

Cytochrome C tryptic digest has been studied in this application brief.

---

## Experimental

### HPLC Conditions

Column:	XTerra RP <sub>18</sub> 4.6 x 50 mm, 3.5 µm (p/n: 186000434)
Mobile phase A:	H <sub>2</sub> O
Mobile phase B:	ACN
Mobile phase C:	100 mM HCOOH, pH 2.5
Flow rate:	0.75 mL/min, 0.2 mL/min to MS
Injection volume:	50 µL (25 µg)
Temperature:	Ambient
Detection:	MS ESI+
Instrument:	Alliance 2790, Micromass Quattro Ultima

## Gradient

Time (min)	Profile		
	%A	%B	%C
0.0	90	0	10
45.0	50	40	10
45.0	10	80	10
50.0	10	80	10

## MS Conditions

MS sytem: Micromass ZQ

Ion source: ESI+

Capillary: 3.5 kV

Cone: 30 V

Extractor: 3.0

Source temp.: 150 °C


Desolvation temp.: 400 °C

Cone gas flow: 60 L/hr

Desolvation gas flow: 550 L/hr

---

## Results and Discussion



---

## Featured Products

Alliance HPLC <<https://www.waters.com/514248>>

WA20738.029, June 2002