


Nota applicativa

## Oxybutynin in Rat Plasma

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---


### Abstract

This application brief highlights the analysis of oxybutynin using XTerra MS C<sub>18</sub> columns.

---

## Introduction

Oxybutynin in rat plasma has been studied in this application brief.



Oxybutynin

---

## Experimental

### HPLC Conditions


Column:	XTerra MS C <sub>18</sub> 2.1 x 30 mm, 3.5 μm (p/n: 186000398)
Mobile phase A:	1.0% NH <sub>4</sub> OH
Mobile phase B:	ACN
Isocratic mobile phase composition:	45% A; 55% B

Flow rate:	0.2 mL/min
Injection volume:	30 $\mu$ L
Detection:	MS ESI+
Instrument:	Alliance 2790, Micromass Quattro Ultima

## MS Conditions

Ion source:	ESI+
Source temp.:	150 $^{\circ}$ C
Gas cell:	$1.5 \times 10^{-3}$ mbar, 25 eV
Desolvation temp.:	350 $^{\circ}$ C
Cone gas flow:	150 L/hr
Drying gas flow:	600 L/hr
Cone voltage:	30 V

Oasis® MCX Extraction Method  
Oasis® MCX Extraction Plate, 10 mg/96 well  
Part Number 186000259




---

## Results and Discussion

## CID mass spectra

cid25  
6 1 ( 1.005)

Daughters of 358ES+  
1.37e8



---

## Featured Products

Alliance HPLC <<https://www.waters.com/514248>>

WA20738.083, June 2002