


Note d'application

Verapamil - pH 2.5, LC-MS

Waters Corporation




This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief highlights the analysis of verapamil by LC-MS using XTerra MS C₁₈ columns.

Introduction

Verapamil has been analyzed in this application brief.



Experimental


HPLC Conditions

| | |
|-------------------------------------|--|
| Column: | XTerra MS C ₁₈ 2.1 x 30 mm, 3.5 µm (p/n: 186000398) |
| Mobile phase A: | 0.1% HCOOH in H ₂ O, pH 2.5 |
| Mobile phase B: | 0.1% HCOOH in ACN, pH 2.5 |
| Flow rate: | 0.2 mL/min to MS |
| Isocratic mobile phase composition: | 77% A; 23% B |
| Injection volume: | 20 µL of 100 pg/µL |
| Temperature: | Ambient |
| Detection: | MS ESI+, SIR 455.45 |
| Instrument: | Alliance 2795 HT, Micromass ZQ |

MS Conditions

| | |
|------------------------------|--------------|
| MS system: | Micromass ZQ |
| Source: | ESI+ |
| Capillary (KV): | 3.0 |
| Cone (V): | 35 |
| Extractor: | 3.0 |
| RF Lens: | 0.5 |
| Source temp.: | 150 |
| Desolvation temp.: | 350 |
| Cone gas flow (L/Hr): | 60 |
| Desolvation gas flow (L/Hr): | 500 |
| LM resolution: | 15 |
| HM resolution: | 15 |
| Ion energy: | 1.0 |
| Multiplier (V): | 650 |

Results and Discussion



Featured Products

Alliance HPLC <<https://www.waters.com/514248>>

WA20738.115, June 2002