

Simultaneous Analysis of Risperidone, 9-Hydroxyrisperidone, and Haloperidol in Plasma

Michelle Wood, Michael Morris

Waters Corporation



For research use only. Not for use in diagnostic procedures.

Abstract

In this application note, we describe a rapid and sensitive method for the simultaneous analysis of haloperidol, risperidone, and 9-hydroxyrisperidone in plasma.

Introduction

Risperidone (Risperdal) is an antipsychotic drug that is used in the management of schizophrenia. Therapeutic activity is thought to be mediated through the antagonism of both the dopamine type 2 (D_2 , see Figure 1) and serotonin type 2 (5-HT₂) receptors of the central nervous system. Risperidone is extensively metabolized in the liver to 9-hydroxyrisperidone, which is the predominant circulating species. Since the metabolite has the same activity as the parent drug, it is important that both compounds are quantified when monitoring plasma concentrations. Although no therapeutic or toxic range has been established, concentrations are generally monitored to ensure that total drug concentrations do not exceed 200 µg/L in plasma.




Figure 2. Simultaneous analysis of 9-hydroxyrisperidone (A), risperidone (B), haloperidol (C) and internal standard trazodone (D) using LC-MS/MS. The above responses were obtained with a 10 μ L injection of the 5 μ g/L plasma calibrator.

Compound name: Risperidone

Correlation coefficient: $r = 0.998492$, $r^2 = 0.996987$

Calibration curve: $1.18683 * x + 0.0443565$

Response type: Internal Std (Ref 1), Area * (IS Conc. / IS Area)

Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




Figure 3. Standard curve for risperidone extracted from human plasma (the inset figure has been expanded to show the response from 0-10 $\mu\text{g/L}$). Fifty microlitres of sample was extracted with 100 μL precipitant prior to analysis by LC-MS/MS.

Conclusion

We describe a rapid and sensitive method for the simultaneous analysis of haloperidol, risperidone and 9-hydroxyrisperidone in plasma. The method involves a simple protein precipitation step prior to analysis using LC-MS/MS and is sufficiently sensitive to enable the quantification of these drugs in plasma taken from patients who are undergoing therapy with these particular antipsychotic agents. The wide dynamic range of the developed LC-MS/MS technique also lends itself well to the monitoring of plasma concentrations where wide inter-individual variations are often observed.

References

1. Janssen Pharmaceutica Products. Package insert for Risperdal. December 2000.
2. McNeil Pharmaceutical. Package insert for Haldol. October 1998.

Featured Products

- [MassLynx MS Software <https://www.waters.com/513662>](https://www.waters.com/513662)
- [OpenLynx Open Access <https://www.waters.com/10008851>](https://www.waters.com/10008851)

720000703, August 2003



© 2021 Waters Corporation. All Rights Reserved.