

## Prochlorperazine - Isolation of Degradation Products, Transfer from Analytical to Prep

---

Waters Corporation



---

### Abstract


This application brief highlights the analysis of Prochlorperazine using Symmetry and SymmetryPrep columns.



---

## Introduction

This application brief highlights the isolation of degradation products of Prochlorperazine and transfer from analytical to preparatory scale.



1. Prochlorperazine Edisylate

---

## Experimental

### HPLC Method

Column:	Symmetry C <sub>18</sub> , 3.9 x 150 mm, 5 µm (p/n: WAT046980)  SymmetryPrep C <sub>18</sub> , 7.8 x 150 mm, 7 µm (p/n: WAT066288)
Mobile phase A:	0.1% TFA in water
Mobile phase B:	Acetonitrile
Flow rate:	0.7 mL/min  2.8 mL/min
Injection volume:	Prochlorperazine edisylate  A. 0.8 mg, B. 3.2mg



Detection:

UV @ 280 nm



Gradient

Time (min)	Profile	
	%A	%B
0	90	10
50	40	60

---

Results and Discussion





---

## Featured Products

WA31763.130, June 2003



