

Bupropion

Waters Corporation




This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief demonstrates analysis of bupropion.

Introduction

The compound analyzed in this study is bupropion (USP Tailing Factor - 1.05).




1. Bupropion

Experimental

HPLC Method

| | |
|-------------------|-----------------------------------------------------|
| Column: | Symmetry C ₈ , 3.9 x 150 mm, 5 µm |
| Guard column: | Sentry Guard Column, 3.9 x 20 mm, 5 µm |
| Part number: | Column - WAT046970, Guard - WAT054250 |
| Mobile phase: | 20 mM potassium phosphate, pH 7.0/methanol 35:65 |
| Flow rate: | 1.0 mL/min |
| Injection volume: | 10 µL of 20 µg/mL sample |
| Detection: | UV @ 254 nm |

Results and Discussion



Featured Products

WA31763.43, June 2003

