

Note d'application

## Propranolol in Rat Plasma by Mixed-Mode SPE and LC-MS/MS (Low Elution Volume)

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---


### Abstract

This application brief demonstrates the analysis of Propranolol in Rat Plasma by LC-MS/MS.

---

## Introduction

The compound analyzed in this study is Propranolol.



PROPRANOLOL

---

## Experimental

### HPLC Method

Column:	XTerra MS C <sub>18</sub> , 2.1 x 30 mm, 3.5 μm
Part number:	186000398
Mobile phase A:	Water + 0.5 % NH <sub>4</sub> OH
Mobile phase B:	ACN + 0.5 % NH <sub>4</sub> OH
Flow rate:	0.2 mL/min

Temperature: Ambient

LC instrument: Alliance 2795

## Gradient

Time (min)	Profile
------------	---------

	%A
--	----

0	5
---	---

1	95
---	----

## MS Conditions

MS instrument: Micromass Quattro Triple Quadrupole

Ion source: ESI+

Source temp.: 150 °C

Gas cell: 2.0 e<sup>-3</sup> bar Argon

Desolvation temp.: 350 °C

Drying gas flow: 500 L/hr

Cone gas flow: 50 L/hr

Cone voltage: 25 V

Collision energy: 20

Capillary voltage: 3.5 Kv

MRM transition:


Metoclopramide (IS)  $m/z$  299.8  $\rightarrow$  226.7

Propranolol  $m/z$  259.9  $\rightarrow$  154.9

Amitriptyline  $m/z$  278.1  $\rightarrow$  232.9

Nortriptyline  $m/z$  263.9  $\rightarrow$  190.8


# NEW! OASIS<sup>®</sup> $\mu$ ELUTION PLATE



## OASIS® MCX GENERIC EXTRACTION PROTOCOL

Conditions for Oasis® MCX  $\mu$ Elution Plate



Part Number 186001830



---

Results and Discussion

# CALIBRATION CURVE



---

## Featured Products

WA31764.132, June 2003

© 2021 Waters Corporation. All Rights Reserved.