

## Propranolol in Rat Plasma by Mixed-Mode SPE and LC-MS/MS (Low Elution Volume)

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.


### Abstract

This application brief demonstrates the analysis of Propranolol in Rat Plasma by LC-MS/MS.

---

## Introduction

The compound analyzed in this study is Propranolol.



## PROPRANOLOL

## Experimental

### HPLC Method

|                 |   |
|-----------------|---|
| Column:         | XTerra MS C <sub>18</sub> , 2.1 x 30 mm, 3.5 µm |
| Part number:    | 186000398                                       |
| Mobile phase A: | Water + 0.5 % NH <sub>4</sub> OH                |
| Mobile phase B: | ACN + 0.5 % NH <sub>4</sub> OH                  |
| Flow rate:      | 0.2 mL/min                                      |
| Temperature:    | Ambient   |
| LC instrument:  | Alliance 2795                                   |

## Gradient

| Time (min) | Profile |
|------------|---------|
|            | %A      |
| 0          | 5       |
| 1          | 95      |

## MS Conditions

|                    |  |
|--------------------|--|
| MS instrument:     | Micromass Quattro Triple Quadrupole  |
| Ion source:        | ESI+   |
| Source temp.:      | 150 °C   |
| Gas cell:          | 2.0 e <sup>-3</sup> bar Argon  |
| Desolvation temp.: | 350 °C   |
| Drying gas flow:   | 500 L/hr   |
| Cone gas flow:     | 50 L/hr  |
| Cone voltage:      | 25 V   |
| Collision energy:  | 20   |
| Capillary voltage: | 3.5 Kv   |
| MRM transition:    | Metoclopramide (IS) <i>m/z</i> 299.8 → 226.7<br>Propranolol <i>m/z</i> 259.9 → 154.9 |

---

Amitriptyline  $m/z$  278.1  $\rightarrow$  232.9

Nortriptyline  $m/z$  263.9  $\rightarrow$  190.8


# NEW! OASIS<sup>®</sup> $\mu$ ELUTION PLATE



# OASIS® MCX GENERIC EXTRACTION PROTOCOL



Conditions for Oasis® MCX  $\mu$ Elution Plate

Part Number 186001830



## Results and Discussion

# CALIBRATION CURVE



## Featured Products

WA31764.132, June 2003



© 2021 Waters Corporation. All Rights Reserved.