

## Protriptyline in Urine by Mixed-Mode Weak Cation Exchange and LC-MS/MS

---

Waters Corporation

This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application brief demonstrates analysis of protriptyline in urine by mixed-mode weak cation exchange and LC-MS/MS.

---

### Introduction

The compound analyzed in this study is Protriptyline.

---



**Protriptyline**

---

## Experimental

### LC Conditions

Column:	XTerra MS C <sub>18</sub> 2.1 x 20 mm <i>IS</i> , 3.5 μm
Part number:	186001923
Mobile phase A:	10 mM NH <sub>4</sub> HCO <sub>3</sub> , pH 10
Mobile phase B:	MeOH with 10 mM NH <sub>4</sub> HCO <sub>3</sub> , pH 10
Flow rate:	0.4 mL/min
Injection volume:	10 μL
Column temperature:	Ambient
Instrument:	Waters 2777 Sample Manager and Waters 1525μ

---

## Gradient

Time (min)	%A	%B
0.0	95	5
3.0	5	95
4.0	5	95
4.1	95	5
5.0	95	5

## MS Conditions

Waters Micromass Quattro Ultima


ESI+

Source temp.: 150 °C

Desolvation temp.:	350 °C
Cone gas flow:	50 L/Hr
Desolvation gas flow:	550 L/Hr
Collision cell:	2.2e <sup>-3</sup> bar (Argon gas)
Cone voltage:	60 volts
CID:	25eV
MRM transition:	<i>m/z</i> 264.0 → 191.1

## Oasis® WCX $\mu$ Elution Plate

Part Number: 186002499




---

## Results and Discussion

---

## 107% Recovery



---

## Featured Products

WA31764.198, June 2003

© 2022 Waters Corporation. All Rights Reserved.