

## Pharmaceutical Compounds in Plasma

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application brief demonstrates analysis of pharmaceutical compounds in plasma.

---


## Introduction

The compounds used in this study are –


1. Caffeine
2. Practolol
3. N-acetyl procainamide
4. Propranolol
5. Methoxyverapamil
6. Amitriptyline




Caffeine




Practolol




N-acetyl procainamide



Propranolol



Methoxyverapamil



Amitriptyline

---

Experimental

LC Conditions

|                       |                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------|
| Column:               | Xterra MS C <sub>18</sub> , 2.1 x 20 mm <i>IS</i> , 3.5 μm, (p/n: 186001923)                |
| Mobile phase A:       | 10 mM NH <sub>4</sub> HCO <sub>3</sub> , pH 10                                              |
| Mobile phase B:       | Methanol                                                                                    |
| Flow rate:            | 0.4 mL/min                                                                                  |
| Injection volume:     | 20 μL                                                                                       |
| Sample concentration: | 5 μg/mL                                                                                     |
| Temperature:          | Ambient                                                                                     |
| Detection:            | MS                                                                                          |
| Instrument:           | Waters 277 Sample Manager, Waters 1525 Binary HPLC Pump and Waters Micromass Quattro Ultima |

## Gradient

| Time (min) | Profile |    |
|------------|---------|----|
|            | %A      | %B |
| 0.0        | 100     | 0  |
| 5.0        | 5       | 95 |

## MS Conditions

## Quattro Ultima

|                               |      |
|-------------------------------|------|
| ES <sup>+</sup> MRM cone (V): | 5.0  |
| Capillary:                    | 3.5  |
| Source temp. (°C):            | 150  |
| Desolvation temp. (°C):       | 400  |
| Cone gas flow (L/Hr):         | 50   |
| Desolvation gas fLow(L/Hr):   | 550  |
| LM resolution 1 & 2:          | 13.5 |
| HM resolution 1 & 2:          | 13.5 |
| Ion energy 1:                 | 0.4  |
| Ion energy 2:                 | 0.8  |
| Multiplier (V):               | 650  |

## Sample Preparation

|              |                                                                                                  |
|--------------|--------------------------------------------------------------------------------------------------|
| SPE:         | Oasis HLB $\mu$ Elution Plate, (p/n:186001828BA)<br>p/n:186001828BA                              |
| Condition:   | 200 $\mu$ L Methanol                                                                             |
| Equilibrate: | 200 $\mu$ L Water                                                                                |
| Load:        | 150 $\mu$ L Spiked rat plasma with 2% H <sub>3</sub> PO <sub>4</sub> ,<br>diluted with 1:1 Water |

|         |                             |
|---------|-----------------------------|
| Wash:   | 200 µL 5% Methanol in Water |
| Elute:  | 50 µL Methanol              |
| Dilute: | 100 µL Water                |


---

## Results and Discussion

Compounds (MRM Transitions, Cone voltage (V), CID (eV))

1. Caffeine (195 > 138; 20; 20)
2. Practolol (267 > 190.1; 30; 20)
3. N-acetyl procainamide (278 > 205.1; 25;20)
4. Propranolol (260 > 183; 25;18)
5. Methoxyverapamil (485.2 > 165.1; 45; 30)
6. Amitriptyline (278 > 233.1; 25; 17)

The top figure is the total ion current, followed by the extracted ion signals for each of the three analytes.



---

## Featured Products

WA31787.19, June 2003

