


## Separation of a Tryptic Digest of Human Serum Albumin

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract

This application brief demonstrates separation of a tryptic digest of human serum albumin.

---

## Introduction

The compound used in this study is tryptic digest of human serum albumin.

---


## Experimental

Column:	BioSuite C <sub>18</sub> 3 µm PA-A (2.1 x 250 mm)
Eluent A:	0.1% Formic acid in water
Eluent B:	0.1% Formic in acetonitrile
Gradient:	0–40% B in 300 min
Flow rate:	0.16 mL/min
Injection volume:	0.2 mL/min
Temperature:	40 °C
Detection:	UV 220 nm

---

## Results and Discussion

Column: BioSuite™ C<sub>18</sub> 3µm PA-A 2.1 x 250 mm  
Eluent A: 0.1% Formic acid in water  
Eluent B: 0.1% Formic in acetonitrile  
Gradient: 0-40% B in 300 min.  
Flow Rate: 0.2 ml/min  
Temperature: 40° C  
Detection: UV 220 nm



## Featured Products

WA31793, June 2003