

## Camphorsulfonic Acid in Rat Plasma on Oasis WAX

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---


### Abstract

This application brief demonstrates analysis of camphorsulfonic acid in rat plasma on Oasis Wax.

---

## Introduction

Camphorsulfonic acid is a strong acid with a pKa of approximately 1.5. The best SPE recoveries for this type of acid are on Oasis WAX products.



## Experimental

### Test Conditions

## Oasis WAX 10 mg 96-Well Plate

Condition:	500 $\mu$ L MeOH
Equilibrate:	500 $\mu$ L H <sub>2</sub> O
Load:	500 $\mu$ L (250 $\mu$ L rat plasma, diluted 1:1 with 4% H <sub>3</sub> PO <sub>4</sub> )
Wash 1:	500 $\mu$ L 2% FA, pH 2.7
Wash 2:	500 $\mu$ L MeOH
Elute:	250 $\mu$ L (125 $\mu$ L x 2) 5% NH <sub>4</sub> OH in MeOH
Options:	<ol style="list-style-type: none"><li>1. Dilute 250 <math>\mu</math>L H<sub>2</sub>O with 2% FA</li><li>2. Evaporate/ Reconstitute</li><li>3. Direct inject</li></ol>
Inject:	10 $\mu$ L

## Oasis WAX 96-Well Plate $\mu$ Elution Plate

Condition:	200 $\mu$ L MeOH
Equilibrate:	200 $\mu$ L H <sub>2</sub> O
Load:	100 $\mu$ L (50 $\mu$ L rat plasma diluted 1:1 with 4% H <sub>3</sub> PO <sub>4</sub> )
Wash 1:	200 $\mu$ L 2% FA, pH 2.7
Wash 2:	200 $\mu$ L MeOH

Elute:	50 $\mu$ L (25 $\mu$ L x 2) 5% $\text{NH}_4\text{OH}$ in MeOH
Options:	<ol style="list-style-type: none"><li>1. Direct injection</li><li>2. Dilute with 50 <math>\mu</math>L <math>\text{H}_2\text{O}</math> with 2% FA</li><li>3. Evaporate/ Reconstitute</li></ol>
Inject:	10 $\mu$ L
Column:	SunFire $\text{C}_{18}$ 2.1 x 20 mm IS, 3.5 $\mu\text{m}$
Mobile phase A:	10 mM $\text{CH}_3\text{COO-NH}_4^+$ , pH 5.5
Mobile phase B:	MeOH with 10 mM $\text{CH}_3\text{COO-NH}_4^+$ , pH 5.5
Flow rate:	0.4 mL /min
Injection volume:	10 $\mu$ L
Column temp:	Ambient
Instrument:	2777 Sample Manager, 1525 $\mu$ Binary HPLC Pump and Quattro Premier

## Gradient


Time (min)	Profile	
	%A	%B
0.0	95	5
3.0	5	95
4.0	5	95
4.1	95	5
5.0	95	5

## Quattro Premier

ESI- source temp:	150 °C
Desolvation temp:	350 °C
Cone gas flow:	50 L /Hr
Desolvation gas flow:	600 L /Hr
Collision cell:	2.2e <sup>-3</sup> bar (Ar gas)

	MRM transition	Cone (V)	CID (eV)
Camphorsulfonic acid	<i>m/z</i> 231.1 → 79.8	60	30

## Results and Discussion



WA60084, June 2007