

응용 자료

Camphorsulfonic Acid in Rat Plasma on Oasis WAX

Waters Corporation




This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief demonstrates analysis of camphorsulfonic acid in rat plasma on Oasis Wax.

Introduction

Camphorsulfonic acid is a strong acid with a pKa of approximately 1.5. The best SPE recoveries for this type of acid are on Oasis WAX products.



Experimental

Test Conditions

Oasis WAX 10 mg 96-Well Plate

Condition:	500 µL MeOH
Equilibrate:	500 µL H ₂ O
Load:	500 µL (250 µL rat plasma, diluted 1:1 with 4% H ₃ PO ₄)
Wash 1:	500 µL 2% FA, pH 2.7
Wash 2:	500 µL MeOH
Elute:	250 µL (125 µL x 2) 5% NH ₄ OH in MeOH
Options:	<ol style="list-style-type: none">1. Dilute 250 µL H₂O with 2% FA2. Evaporate/ Reconstitute3. Direct inject
Inject:	10 µL

Oasis WAX 96-Well Plate µ Elution Plate

Condition:	200 µL MeOH
Equilibrate:	200 µL H ₂ O
Load:	100 µL (50 µL rat plasma diluted 1:1 with 4% H ₃ PO ₄)
Wash 1:	200 µL 2% FA, pH 2.7
Wash 2:	200 µL MeOH

Elute:	50 μ L (25 μ L x 2) 5% NH_4OH in MeOH
Options:	<ol style="list-style-type: none"> 1. Direct injection 2. Dilute with 50 μL H_2O with 2% FA 3. Evaporate/ Reconstitute
Inject:	10 μ L
Column:	SunFire C_{18} 2.1 x 20 mm IS, 3.5 μ m
Mobile phase A:	10 mM $\text{CH}_3\text{COO-NH}_4^+$, pH 5.5
Mobile phase B:	MeOH with 10 mM $\text{CH}_3\text{COO-NH}_4^+$, pH 5.5
Flow rate:	0.4 mL /min
Injection volume:	10 μ L
Column temp:	Ambient
Instrument:	2777 Sample Manager, 1525 μ Binary HPLC Pump and Quattro Premier

Gradient


Time (min)	Profile	
	%A	%B
0.0	95	5
3.0	5	95
4.0	5	95
4.1	95	5
5.0	95	5

Quattro Premier

ESI- source temp:	150 °C
Desolvation temp:	350 °C
Cone gas flow:	50 L /Hr
Desolvation gas flow:	600 L /Hr
Collision cell:	2.2e ⁻³ bar (Ar gas)

	MRM transition	Cone (V)	CID (eV)
Camphorsulfonic acid	<i>m/z</i> 231.1 → 79.8	60	30

Results and Discussion



WA60084, June 2007