

응용 자료

## Oasis 2x4 Method: Proof of Concept

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application brief highlights Oasis 2x4 Method which is a viable approach to SPE sorbent and protocol selection.

# Introduction

In order to prove that the Oasis 2x4 Method is a viable approach to SPE sorbent and protocol selection, a group of representative small molecules was spiked into rat plasma. The molecules include a base, a quaternary ammonium salt, a neutral, an acid and a strong acid. The spiked plasma samples were then extracted following the Oasis 2x4 Method. All molecules elute in the correct Elute 1 or Elute 2 fractions according to what theory predicted, proving that the method is viable.

Follow the simple steps outlined in this flow chart to achieve high recoveries and the cleanest extracts:

- Characterize your analyte [Neutral, Acid or Base, pKa].
- Select one of the four Oasis sorbents.
- Apply the indicated Protocol [1 or 2].
- Determine SPE recoveries by LC analysis.



---

## Experimental

10 mg Oasis 96-well Plates

### Protocol 1

Condition:	500 $\mu$ L MeOH
Equilibrate:	500 $\mu$ L H <sub>2</sub> O
Load:	500 $\mu$ L sample (250 $\mu$ L plasma diluted 1:1 with 4% H <sub>3</sub> PO <sub>4</sub> in H <sub>2</sub> O)
Wash 1:	500 $\mu$ L 2% HCOOH in H <sub>2</sub> O
Elution 1:	2 x 125 $\mu$ L MeOH
Elution 2:	2 x 125 $\mu$ L 5% NH <sub>4</sub> OH in MeOH
Dilution:	250 $\mu$ L water

### Protocol 2

Condition:	500 $\mu$ L MeOH
Equilibrate:	500 $\mu$ L H <sub>2</sub> O
Load:	500 $\mu$ L sample (250 $\mu$ L plasma diluted 1:1 with 4% H <sub>3</sub> PO <sub>4</sub> in H <sub>2</sub> O)
Wash 1:	5% NH <sub>4</sub> OH in H <sub>2</sub> O
Elution 1:	2 x 125 $\mu$ L MeOH


Elution 2:	2 x 125 $\mu$ L 2% HCOOH in MeOH
Dilution:	5% NH <sub>4</sub> OH in H <sub>2</sub> O (To neutralize acid for high pH LC)
Column:	XBridge C <sub>18</sub> 2.1 x 20 mm IS, 3.5 $\mu$ m
Mobile Phase A:	10 mM NH <sub>4</sub> HCO <sub>3</sub> , pH 10
Mobile Phase B:	10 mM NH <sub>4</sub> H
Injection Volume:	10.0 $\mu$ L
Column Temperature:	Ambient
Detection:	UV @ 254 nm (Prednisone)
Instrumentation:	2777 Sample Manager, 1525 $\mu$ Binary HPLC Pump, Quattro Premier and 2996 PDA

## Gradient


Time (min)	Profile
	%A
0	95
3.0	5
4.8	5
5.0	95

Time (min)	Profile
7.0	95
Quattro Premier	
Capillary:	3.4 kV
Source Temp.:	120 °C
Desolvation Temp.:	350 °C
Cone Gas Flow:	50 L /Hr
Desolvation Gas Flow:	700 L /Hr
Collision Cell Pressure:	2.59e <sup>-3</sup> mbar
MRM Transitions:	Imipramine 281.2 > 85.95 ESI+ Decanesulfonic Acid 220.9 > 79.7 ESI - Ibuprofen 205.2 > 161 ESI - Valethamate 306.3 > 162.8 ESI+


## Results and Discussion




Imipramine (B)  
pKa = 9.4  
100 ng/mL




Prednisone (N)  
20 µg/mL



Valethamate (QA)  
pKa >12  
100 ng/mL




Decanesulfonic Acid (SA)  
pKa <0.5  
200 ng/mL




Ibuprofen (A)  
pKa = 5.2  
100 ng/mL


## Elute 2, Primary Analyte Data



## Prednisone SPE Recovery Data



## Elute 1, Counter Analyte Data



## Featured Products

2998 Photodiode Array (PDA) Detector <<https://www.waters.com/1001362>>

WA60090, June 2007