

Multiresidue Analysis of Pesticides in Vegetables and Fruits

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief shows the Japan Ministry of Health, Labor and Welfare (JPMHLW) Method for multi-residue analysis of pesticides in vegetables and fruit.


Introduction

This application brief shows the Japan Ministry of Health, Labor and Welfare (JPMHLW) Method for multi-residue analysis of pesticides in vegetables and fruit. This sample preparation method calls for an extract from the commodity, followed by a SPE extract from a Sep-Pak Vac Carbon Black/Aminopropyl column.

Experimental

SPE Procedure

Sep-Pak® Vac Carbon Black/Aminopropyl 6 cc/500 mg/500 mg



LC Conditions

| | |
|-----------------|--------------------------------------------------|
| System: | Alliance HPLC 2695 |
| Column: | XTerra MS C ₁₈ , 3.5 µm, 2.1 x 150 mm |
| Flow rate: | 0.2 mL/min |
| Mobile phase A: | water |
| Mobile phase B: | methanol |

Mobile phase C: 100 mM ammonium acetate in water

Injection volume: 5 µL

Column temp.: 40 °C

Gradient

| Time (min) | %A | %B | %C |
|------------|----|----|----|
| 0 | 80 | 15 | 5 |
| 1 | 55 | 40 | 5 |
| 3.5 | 55 | 40 | 5 |
| 6 | 45 | 50 | 5 |
| 8 | 40 | 55 | 5 |
| 17.5 | 0 | 95 | 5 |
| 30 | 80 | 15 | 5 |
| 47 | 80 | 15 | 5 |

MS Conditions

MS System: Waters Quattro Premie XE

Ionization mode: Positive electrospray (ESI⁺)
Multiple reaction monitoring

Results and Discussion

| Pesticides* | Spike Conc. / $\mu\text{g/g}$ | Recovery (%) |
|---------------------------|-------------------------------|--------------|
| Abamectin | 0.1 | 102.0 |
| Anibfos | 0.1 | 111.7 |
| Azinphos-methyl | 0.1 | 107.6 |
| Benzofenap | 0.1 | 139.5 |
| Butafenacil | 0.1 | 104.5 |
| Chloridazon | 0.1 | 106.0 |
| Chromafenozide | 0.1 | 108.2 |
| Clomeprop | 0.1 | 104.4 |
| Cloquintocet-mexyl | 0.1 | 108.7 |
| Clothianidin | 0.1 | 101.5 |
| Cyazofamid | 0.1 | 108.3 |
| Cyflufenamid | 0.1 | 110.1 |
| Dimethirimol | 0.1 | 101.0 |
| Fenoxycarb | 0.1 | 108.7 |
| Ferimzone | 0.1 | 112.6 |
| Formetanate hydrochloride | 0.1 | 86.7 |
| Furathiocarb | 0.1 | 100.5 |
| Imidacloprid | 0.1 | 111.8 |
| Indoxacarb | 0.1 | 121.2 |
| Iprovalicarb | 0.1 | 106.2 |
| Isoxaflutole | 0.1 | 99.5 |
| Lactofen | 0.1 | 106.8 |
| Methoxyfenozide | 0.1 | 103.3 |
| Mibemectin A3 | 0.1 | 114.5 |
| Mibemectin A4 | 0.1 | 101.2 |
| Naproanilide | 0.1 | 115.9 |
| Oryzalin | 0.1 | 103.8 |
| Oxycarboxin | 0.1 | 85.1 |
| Oxydemeton-methyl | 0.1 | 108.0 |
| Phenmedipham | 0.1 | 102.2 |
| Pyrazolynate | 0.1 | 72.7 |
| Quizalofop-P-tefuryl | 0.1 | 145.3 |
| Simeconazole | 0.1 | 106.0 |
| Thiacloprid | 0.1 | 109.2 |
| Thiamethoxam | 0.1 | 108.3 |
| Tridemorph | 0.1 | 94.6 |
| Etriticonazole | 0.1 | 113.3 |

*Five replicate samples analyzed per level.

Featured Products

Alliance HPLC System <<https://www.waters.com/534293>>

720002586, April 2008

© 2021 Waters Corporation. All Rights Reserved.