

Nota applicativa

Determination of Polycyclic Aromatic Hydrocarbons in Municipal and Industrial Wastewater

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This Application note demonstrates the determination of Polycyclic Aromatic Hydrocarbons in municipal and industrial wastewater.

Introduction

Polycyclic aromatic hydrocarbons (PAHs) are one of the most widespread organic pollutants. PAHs are made up of fused aromatic rings and are formed during the combustion of carbon based fuels (wood, coal, diesel), as well as being present in crude oil. The United States Environmental Protection Agency (US EPA) has classified seven PAH compounds as being potentially carcinogenic including benz[a]anthracene, benzo[a]pyrene, benzo[b]fluoranthene, benzo[k]fluoranthene, chrysene, dibenz[a,h]anthracene, and indeno[1,2,3-cd]pyrene.

Experimental

HPLC conditions

| | |
|-------------|--|
| Instrument: | Waters Alliance HPLC system with PDA and fluorescence detectors |
| Eluent: | Water/acetonitrile |
| Column: | Waters PAH 4.6 x 250 mm @ 30 °C |
| Injection: | 20 µL of Supelco Standard EPA 610 (#48743) diluted 1:50 in 40:60 water/acetonitrile |
| Flow Rate: | 1.2 mL/min |

Detection: UV @ 254 nm and fluorescence using timed programmed wavelengths

Data: Waters Empower software

Sample preparation

Liquid/liquid extraction with MeCl₂

Eluent preparation

Filter and degas through a 0.45 µm filter.


A: Water

B: Acetonitrile


| Time | Flow | %A | %B | Curve |
|---------|------|----|-----|-------|
| Initial | 1.2 | 40 | 60 | - |
| 12.0 | 1.2 | 0 | 100 | 9 |
| 23.0 | 1.2 | 40 | 60 | 11 |

Eluent gradient.

Results and Discussion



Standard chromatogram, UV @ 254 nm, 1-20 ppm PAH analytes.



Standard chromatogram, fluorescence/programmed wavelengths, 1-20 ppm PAH analytes.

| | Analyte | UV max (nm) | EX (nm) | EM (nm) | Detection Limit (ppb) ¹ |
|----|-----------------------------|-------------|---------|---------|------------------------------------|
| 1 | Naphthalene | 220 | 277 | 330 | 0.14 |
| 2 | Acenaphthylene | 229 | NA | NA | NA |
| 3 | Acenaphthene | 227 | 270 | 323 | 0.01 |
| 4 | Fluorene | 261 | 265 | 310 | 0.03 |
| 5 | Phenanthrene | 251 | 252 | 365 | 0.02 |
| 6 | Anthracene | 252 | 250 | 402 | 0.01 |
| 7 | Fluoranthene | 236 | 284 | 467 | 0.02 |
| 8 | Pyrene | 240 | 332 | 378 | 0.01 |
| 9 | Benzo(a)anthracene | 287 | 284 | 390 | 0.01 |
| 10 | Chrysene | 267 | 270 | 367 | 0.04 |
| 11 | Benzo(b)fluoranthene | 256 | 298 | 436 | 0.09 |
| 12 | Benzo(k)fluoranthene | 307 | 303 | 432 | 0.01 |
| 13 | Benzo(a)pyrene ² | 296 | 280 | 410 | 0.03 |
| 14 | Dibenzo(a,h)anthracene | 297 | 294 | 398 | 0.01 |
| 15 | Benzo(g,h,i)perylene | 299 | 290 | 420 | 0.03 |
| 16 | Indeno(1,2,3-cd)pyrene | 250 | 305 | 480 | 0.49 |

1 Fluorescence mode used for detection limit determination, no pre concentration.

Seven replicates per 40 CFR pt. 136 App. B.

2 Regulated compound; action level 0.17 ppb.

PAH target analytes.

References

1. The Determination of Biodegradation Products of PAH Using LC-MS/MS: WA20747
2. PAHs in Drinking Water – Oasis Solution: WA31764.127
3. Waters PAH Columns Improve Analysis of PAH Compounds: 720000382EN

Featured Products

Alliance HPLC System <<https://www.waters.com/534293>>

Empower 3 Chromatography Data Software <<https://www.waters.com/513188>>

2998 Photodiode Array (PDA) Detector <<https://www.waters.com/1001362>>

2475 Fluorescence (FLR) Detector <<https://www.waters.com/514434>>

720002733, August 2008