

## ACQUITY UPLC Analysis of Thiourea

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---


### Abstract

This application brief demonstrates the analysis of Thiourea.

---

## Introduction

### Structure



**Thiourea**

---

## Experimental


### Test Conditions

|                         |   |
|-------------------------|---|
| Column:                 | ACQUITY UPLC BEH Amide, 2.1 x 150 mm, 1.7 $\mu$ m   |
| Part Number:            | 186004802   |
| Isocratic Mobile Phase: | 95/2.5/2.5 MeCN/IPA/H <sub>2</sub> O with 10 mM CH <sub>3</sub> COONH <sub>4</sub> and 0.01% NH <sub>4</sub> OH, pH 9.0 |
| Flow Rate:              | 0.2 mL/min  |
| Injection Volume:       | 5.0 $\mu$ L (PLNO)  |

|                       |   |
|-----------------------|---|
| Sample Concentration: | 10 µg/mL  |
| Sample Diluent:       | 75/25 MeCN/MeOH with 0.2% HCOOH                       |
| Column Temperature:   | 25 °C   |
| Weak Needle Wash:     | 95/5 MeCN/H <sub>2</sub> O                            |
| Detection:            | UV @ 245 nm   |
| Sampling Rate:        | 20 points/sec   |
| Filter Time Constant: | 0.2   |
| Instrument:           | Waters ACQUITY UPLC with ACQUITY UPLC<br>PDA Detector |

---

## Results and Discussion



---

## Featured Products

ACQUITY UPLC System <<https://www.waters.com/514207>>

ACQUITY UPLC PDA Detector <<https://www.waters.com/514225>>

WA60095, June 2009