

Note d'application

ACQUITY UPLC Analysis of Morphine

Waters Corporation



Abstract

This application brief highlights the analysis of morphine on ACQUITY UPLC BEH Amide Columns.

Introduction

The compounds analysed in this study are:

1. 6-acetylmorphine (100 ng/mL)
2. Morphine (100 ng/mL)
3. Morphine-3 β -D-glucuronide (5 g/mL)



Experimental

Chromatographic Conditions

Column:	ACQUITY UPLC BEH Amide, 2.1 x 50 mm, 1.7 μ m
Part Number:	186004800
Mobile Phase A:	50/50 MeCN/H ₂ O with 10 mM NH ₄ COOH and 0.125% HCOOH, pH 3

Mobile Phase B:	90/10 MeCN/H ₂ O with 10 mM NH ₄ COOH and 0.125% HCOOH, pH 3
Flow Rate:	0.6 mL/min
Injection Volume:	5 µL (PLNO)
Sample Diluent:	75/25 MeCN/MeOH with 0.2% HCOOH
Column Temperature:	30 °C
Weak Needle Wash:	95/5 MeCN/H ₂ O
Instrument:	Waters ACQUITY UPLC with ACQUITY SQD

Gradient


Time (min)	Profile	
	%A	%B
Initial	0.1	99.9
1.05	0.1	99.9
4.35	99.9	0.1
4.50	0.1	99.9
6.00	0.1	99.9

Mass Spectrometer Conditions

Ionization Mode:	ES+
Capillary:	3.0 KV

Cone:	30 V (6-Acetyl morphine and Morphine), 40 V (Morphine-3 β -D-glucuronide)
Source Temperature:	120 °C
Desolvation Temperature:	350 °C
Desolvation Gas Flow (L/Hr):	800
SIR <i>m/z</i> :	329.5 (6-Acetyl morphine); 287.5 (Morphine) 463.6 (Morphine-3 β -D-glucuronide)
Dwell Time:	0.1 s

Results and Discussion



Featured Products

SQ Detector 2 <<https://www.waters.com/134631584>>

WA60101, June 2009

© 2022 Waters Corporation. All Rights Reserved.