

## ACQUITY UPLC Analysis of Histidine Dipeptides

---

Waters Corporation

This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application highlights the analysis of Histidine Dipeptides on ACQUITY UPLC BEH Amide Columns.

---


### Introduction

The compounds used in this study are:


1. Creatinine (1 µg/mL)
  2. Creatine (5 µg/mL)
  3. Anserine (5 µg/mL)
  4. Canosine (5 µg/mL)
-




**Carnosine**



**Creatine**



**Creatinine**



**Anserine**

---

## Experimental

### Chromatographic Conditions

Column:	ACQUITY UPLC BEH HILIC, 2.1 x 50 mm, 1.7 $\mu\text{m}$
Part Number:	186004800
Mobile Phase A:	50/50 MeCN/H <sub>2</sub> O with 10 mM CH <sub>3</sub> COONH <sub>4</sub> and 0.04 % NH <sub>4</sub> OH, pH 9.0
Mobile Phase B:	95/5 MeCN/H <sub>2</sub> O with 10 mM CH <sub>3</sub> COONH <sub>4</sub> and 0.04 % NH <sub>4</sub> OH, pH 9.0

---

Flow Rate: 0.5 mL/min

Injection Volume: 5 µL

Sample Diluent: 75/25 MeCN/MeOH

Column Temperature: 30 °C

Weak Needle Wash: 95/5 MeCN/H<sub>2</sub>O


Instrument: Waters ACQUITY UPLC with ACQUITY SQD

## Gradient

Time (min)	Profile	
	%A	%B
Initial	0.1	99.9
5.00	65.0	35.0
5.01	0.1	99.9
6.00	0.1	99.9

---

## Results and Discussion



---

## Featured Products

SQ Detector 2 <<https://www.waters.com/134631584>>

WA60103, July 2009

© 2022 Waters Corporation. All Rights Reserved.