

## Separation of Nucleotide Phosphates on ACQUITY UPLC BEH Amide



Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

## Abstract


This application brief describes the separation of nucleotide phosphates on ACQUITY UPLC BEH Amide column.

---

## Introduction

The nucleotide phosphates used in this study are-

1. Adenosine monophosphate (AMP)
2. Uridine monophosphate (UMP)
3. Adenosine diphosphate (ADP)
4. Uridine diphosphate (UDP)
5. Adenosine triphosphate (ATP)
6. Uridine triphosphate (UTP)



---

Figure 1. Structures of the compounds used in this study.

---


# Experimental

## Test conditions

Column:	ACQUITY UPLC BEH Amide, 2.1 x 100 mm, 1.7 $\mu$ m
Part Number:	186004801
Isocratic Mobile Phase:	70/30 ACN/H <sub>2</sub> O with 27 mM KH <sub>2</sub> PO <sub>4</sub> , pH 4.5
Flow Rate:	0.5 mL/min
Injection Volume:	5 $\mu$ L (PLNO)
Sample Concentration:	shown on chromatogram
Sample Diluent:	80/20 ACN/H <sub>2</sub> O
Column Temperature:	25 °C
Weak Needle Wash:	95/5 ACN/H <sub>2</sub> O
Instrument:	Waters ACQUITY UPLC with ACQUITY PDA
Detection:	UV 260 nm
Sampling Rate:	20 Hz
Time Constant:	0.1 s

---

## Results and Discussion



*Figure 2. Sample Concentration*

1- AMP (50  $\mu\text{g/mL}$ ), 2- UMP (50  $\mu\text{g/mL}$ ), 3- ADP (100  $\mu\text{g/mL}$ ), 4- UDP (100  $\mu\text{g/mL}$ ), 5- ATP (100  $\mu\text{g/mL}$ ), 6- UTP (100  $\mu\text{g/mL}$ )

---

## Featured Products

[ACQUITY UPLC System](#)

[ACQUITY UPLC PDA Detector](#)

WA64069, August 2009

©2019 Waters Corporation. All Rights Reserved.