


LC-MS Isocratic Separation of Neurotransmitters on Atlantis HILIC Silica

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

Abstract


This application brief demonstrates the LC-MS isocratic separation of neurotransmitters on Atlantis HILIC Silica

Columns.


Introduction

The compounds used in this study are:

1. Acetylcholine
2. Choline



Acetylcholine (Ach)



Choline (Ch)

Experimental

LC Conditions

Column:	Atlantis HILIC Silica, 2.1 x 50 mm, 3 μ m
Part Number:	186002011
Mobile Phase A:	10 mM NH ₄ COOH with 0.125% HCOOH in H ₂ O
Mobile Phase B:	10 mM NH ₄ COOH with 0.125% HCOOH in 90:5:5 ACN:MeOH:H ₂ O

Isocratic Mobile Phase Composition:	10% A; 90% B
Flow Rate:	0.5 mL/min
Injection Volume:	10.0 µL (full loop)
Sample Concentration:	5 ng/mL each
Sample Diluent:	75:25 ACN:MeOH with 0.2% HCOOH
Column Temp.:	30 °C
Weak and Strong Needle Wash:	95:5 ACN:H ₂ O
Detection:	MS
Sampling Rate:	5 points/sec
Instrument:	Waters ACQUITY UPLC with TQD

MS Conditions

Ionization Mode:	ES+
Capillary:	0.5 kV
Cone:	25 V (acetylcholine), 40 V (choline)
Source Temperature:	120 °C
Desolvation Temperature:	350 °C
Desolvation Gas Flow:	800 L/Hr


SIR: 146.1 *m/z* (acetylcholine); 104.0 *m/z* (choline)

Dwell Time: 150 msec

ISD: 10 msec

ICD: 10 msec

Results and Discussion



Featured Products

· [ACQUITY UPLC System <https://www.waters.com/514207>](https://www.waters.com/514207)

WA64071, August 2009



© 2021 Waters Corporation. All Rights Reserved.