

Multi-Residue Analysis of Pesticides in Grapes Using AOAC QuEChERS Method by UPLC-MS/MS

日本ウォーターズ株式会社



This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief demonstrates multi-residue analysis of pesticides in grapes using AOAC QuEChERS method by UPLC-MS/MS.

Experimental

Test Conditions

LC Conditions

| | |
|-------------------|---|
| LC System: | Waters ACQUITY UPLC System |
| Column: | ACQUITY UPLC BEH C ₁₈ , 2.1 x 100 mm, 1.7 μm |
| Column Temp: | 40 °C |
| Sample Temp: | 4 °C |
| Flow Rate: | 0.3 mL/min. |
| Mobile Phase A: | Water + 0.1% formic acid |
| Mobile Phase B: | Methanol + 0.1% formic acid |
| Injection Volume: | 15 μL, Partial loop injection |

Gradient:


| Time | Flow Rate | A% | B% |
|------|-----------|----|----|
| 0.00 | 0.3 | 75 | 25 |

| Time | Flow Rate | A% | B% |
|-------|-----------|----|-----|
| 0.25 | 0.3 | 75 | 25 |
| 7.75 | 0.3 | 5 | 100 |
| 8.50 | 0.3 | 0 | 100 |
| 8.51 | 0.5 | 75 | 25 |
| 10.50 | 0.5 | 75 | 25 |
| 11.0 | 0.3 | 75 | 25 |

MS Conditions

| | |
|--------------|------------------------------------|
| Instrument: | Waters ACQUITY TQ Detector |
| Ionization: | Positive electrospray (ESI+) |
| Acquisition: | Multiple reaction monitoring (MRM) |

Results and Discussion



Featured Products

ACQUITY UPLC System <<https://www.waters.com/514207>>

720003369, March 2010