

Multi-Residue Analysis of Pesticides in Grapes Using AOAC QuEChERS Method by GC-MS

日本ウォーターズ株式会社



This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief demonstrates multi-residue analysis of pesticides in grapes using AOAC QuEChERS method by GC-MS.

Experimental

Extraction Procedure

1. Add 15 mL 1% acetic acid in acetonitrile into the 50 mL DisQuE extraction tube 1.
2. Add 15 g of homogenized sample into the 50 mL tube.
3. Add any internal standards and standard mixture.
4. Shake vigorously for 1 minute and centrifuge > 1500 rcf for 5 minute.
5. Transfer 1 mL of the acetonitrile extract into the clean-up tube 2.
6. Shake for 30 seconds and centrifuge >1500 rcf for 1 minute.
7. Transfer 0.5 mL extract into a tube.
8. Add any post-extraction internal standards.
9. Add 0.25 mL toluene.
10. Evaporate at 50 °C with N₂ to < 0.1 mL.
11. Bring volume up to 0.2 mL with toluene.
12. Transfer to vial with insert for analysis.

Test Conditions

GC Conditions

Instrument: Agilent 6890N GC

Column: RTX-5MS, 30 x 0.25 mm, (0.25 µm film)


GC Conditions

| | |
|-------------------|---|
| Carrier Gas: | Helium |
| Flow Rate: | 1.0 mL/min |
| Temp. Program: | Initial 100 °C, hold 1 min, then 10 °C/min to 320 °C, hold for 7 minute |
| Injection Volume: | 2 µL splitless |

MS Conditions

| | |
|--------------|---------------------------------|
| Instrument: | Waters Quattro micro GC-MS |
| Ionization: | Electron Impact (70 eV) |
| Acquisition: | Single Ion Recording (SIR) Mode |

Results and Discussion



Featured Products

720003370, March 2010

©2019 Waters Corporation. All Rights Reserved.