



---

## Abstract

Today's bioanalytical scientists face several challenges in their workflow, regardless of the regulated or non-regulated environment. Some of these challenges include:

- The ability to develop highly sensitive LC-MS methods for any molecule
- Ensuring high reproducibility and robustness while maintaining the shortest possible run time
- Addressing regulatory guidelines
- Ensuring high throughput value

This application note focuses on quantitation of clopidogrel by following the MRM of the active pro-drug with an LLOQ of less than 1 pg/mL in plasma.

### Benefits

Waters Regulated Bioanalysis System Solution is capable of addressing sensitivity challenges in the world of regulated bioanalysis.

---

## Introduction

Clopidogrel is a prodrug, which takes action on an adenosine diphosphate (ADP) receptor on platelet cell membranes. The drug specifically and irreversibly inhibits the P2Y<sub>12</sub> subtype of an ADP receptor, which is important in aggregation of platelets and cross-linking by the protein fibrin.



followed by an organo-aqueous solution, then eluted in an organo-aqueous elution solvent. The eluted samples were directly injected onto the system. For micro-elution plates, the elution resulted in a higher concentration of the sample, and hence higher a signal-to-noise (S/N) ratio compared to that obtained from the standard SPE. Atorvastatin was used as the internal standard (IS) for the estimation of clopidogrel.

---

## Experimental

### LC Conditions

LC system:	ACQUITY UPLC
Column:	ACQUITY UPLC HSS C <sub>18</sub> 1.8 $\mu$ m, 2.1 x 50 mm
LC column:	99% aqueous buffer over 0.4 min followed by a 90% organic elution until 2 min; then change back to initial conditions.
Column temp.:	40 °C
Flow rate:	0.400 mL/min
Injection volume:	5 $\mu$ L

### MS Conditions

MS system:	Xevo TQ-S
MS mode:	ESI positive
MRM transition:	322.1 $\rightarrow$ 212.0








Figure 4.1 Chromatogram of the blank internal standard.













©2019 Waters Corporation. All Rights Reserved.