# Waters<sup>™</sup>

#### アプリケーションノート

Improving the Productivity in Isolating a Naturally Occurring Bioactive Compound Using Supercritical Fluid Extraction and Preparative Supercritical Fluid Chromatography

John McCauley, Jo-Ann M. Jablonski, Andrew J. Aubin, Rui Chen

日本ウォーターズ株式会社



### Abstract

In this application note, we have demonstrated employing SFE and prep SFC to holistically improve the productivity in isolating a low-polarity, bioactive compound from a complex natural product extract. The SFE alleviated the sample complexity prior to analysis and purification, thereby improving sample loading and reducing solvent use in the ensuing chromatography. The SFE extract also lends itself well for SFC analysis and purification.

#### Benefits

- SFE alleviates the sample complexity in natural product extracts prior to chromatographic analysis and purification, enabling a more efficient purification downstream.
- SFC offers complementary separation to RPLC. In addition, there is a wide range of column chemistries available in SFC with vastly different separation mechanisms. The combinations of SFC/RPLC and SFC/SFC provide unmatched resolving power to meet the challenges, primarily arising from sample complexity, in natural product isolation.
- Both SFE and SFC reduce the use of organic solvents and provide an easy sample recovery under mild conditions, thereby increasing the overall purification productivity and cost-effectiveness.

### Introduction

Natural products are a productive source of leads for new drugs due to their high chemical diversity, biochemical specificity, and many "drug-likeness" molecular properties.<sup>1-4</sup> A large portion of today's existing drugs on the market are either directly derived from naturally occurring compounds or inspired by a natural product. In addition, natural products are also used in the forms of food supplements, nutraceuticals, alternative medicines, and as active ingredients in cosmetics.<sup>5</sup>

Isolation and purification of bioactive compounds play an important role in natural product research. The most commonly used process often involves extraction of target compounds from the cellular matrix, prepurification by various chromatographic techniques including flash chromatography (FC), low pressure liquid chromatography (LPLC), and medium pressure liquid chromatography (MPLC), followed by preparative high pressure liquid chromatography (prep HPLC).<sup>6</sup> However, this process is not without its challenges. For example, conventional extraction methods for natural products include Soxhlet extraction, maceration, percolation, and sonication. These methods are often time- and labor-intensive, consume large amounts of organic solvents, and can lead to the degradation of thermally labile compounds. Furthermore, prep chromatography is largely dominated by reversed-phase liquid chromatography (RPLC), whereby the separation is driven by the differentiating polarity of the analytes. While a generally applicable chromatographic technique for a variety of compound classes, RPLC does not necessarily guarantee an adequate resolution for all analytes, especially for the structural analogs and isomers of similar polarities often found in natural products. As a result, the purification step is perceived by many as a rate-limiting step and a major bottleneck for natural product drug discovery, as well as in the development of differentiated nutraceutical and cosmetic products<sup>.7</sup>

To that end, supercritical fluid (SF) based techniques, including supercritical fluid extraction (SFE) and supercritical fluid chromatography (SFC), can offer viable additions to the natural product isolation toolbox by leveraging the unique properties of supercritical CO<sub>2</sub>: high diffusivity, low viscosity, and superb solvation power. SFE has been successfully applied to the extraction of many bioactive compounds from medicinal plants, including steroids, terpenes, alkaloids, and phenolic compounds.<sup>6</sup> Preparative SFC has been widely adopted by the pharmaceutical industry for active pharmaceutical ingredient (API) purification. Its applications in natural product isolation, however, remain scarce.<sup>8</sup>

In this application note, we describe a systematic effort to holistically improve the productivity in isolating a naturally occurring terpene derivative with proven anti-cancer bioactivity from a raw plant sample. The process involves an extraction by SFE followed by three different, two-step purification routes, including MPLC+HPLC, MPLC+SFC, and SFC+SFC. The overall productivity and solvent consumption for each purification route are compared.

## Experimental

#### Materials and reagents

HPLC-grade methanol and isopropanol (IPA) were purchased from Thermo Fisher Scientific (Fair Lawn, NJ, USA). Denatured ethanol (reagent grade) was purchased from Sigma (St. Louis, MO, USA). The fine ground plant material was used as received.

#### Sample preparation

#### Solvent extraction

A total of 0.3 g of ground plant material and 6 mL methanol were placed into a 10 mL test tube. After sonication at 40 °C for 1 hour, the suspension was centrifuged for 5 min. The supernatant was transferred to a clean vial for further analysis.

#### Supercritical fluid extraction

The extraction experiments were performed on a Waters MV-10 ASFE System controlled by ChromScope Sample Prep Software. A total of 3 g of ground plant material was weighed into a 5 mL extraction vessel. The extraction was performed for 60 minutes with 8 mL/min CO<sub>2</sub>. The effluent was carried into a 100 mL collection vessel with a makeup flow of 1 min/mL of methanol/ isopropanol/hexane (1:1:1).

#### Chromatography

Analytical LC-MS experiments were performed on a Waters ACQUITY UPLC H-Class System/SQ Detector 2 and a Waters AutoPurification LC System. The analytical UPC<sup>2</sup>-MS experiments were performed on a Waters ACQUITY UPC<sup>2</sup>-MS System. All systems were controlled by MassLynx software. The MS-directed SFC preparative experiments were performed on a Waters Prep 100q SFC MS-Directed System controlled by MassLynx/FractionLynx Software. All UV-directed preparative experiments were performed on a Waters SFC 80 Preparative System controlled by ChromScope software. Detailed experimental parameters are summarized in Tables 1-3.

	Figure 2	2A	Figure	2B	Figure 6A	
Instrument	ACQUITY UPLC H-Class		AutoPurification		ACQUITY UPLC H-Class	
	System/SQD2 MS		LC MS System		System/SQD2 MS	
Flow rate (mL/min)	0.60		1.46		0.75	
Mobile phase A	Water	•	Water		Water	
Mobile phase B	Methanol		Methanol		Methanol	
Backpressure (psi)	N/A		N/A		N/A	
MS detection	ESI+		ESI+		ESI+	
Column	ACQUITY HSS T3		Atlantis T3		ACQUITY BEH C18	
	(1.8 µm, 3.0 x 150 mm)		(5 µm, 4.6 x 150 mm)		(1.7 µm, 2.1 x 50 mm)	
Temperature (°C)	60		Ambient		60	
Injection volume (µL)	1		Varying		0.5	
Gradient	Time (min)	%B	Time (min)	%B	Time (min)	%B
	0	92	0	88	0	80
	5	96	3.08	88	4	80
	5.25	92	8.21	94		
	6	92	8.61	100		
			9.22	88		
			20.90	88		

Table 1. Key experimental parameters for analytical LC.

	Figure 1	I	Figure 3	A	Figure 4	A	Figure	4B	Figure 6	6B
Instrument					ACQUITY UPC <sup>2</sup> System/TQD MS					
Flow rate (mL/min)				1.5						
Backpressure (psi)					1740					
MS Detection					APCI+					
Temperature (°C)					45					
Injection volume (µL)					1					
Mobile phase A					CO <sub>2</sub>					
Mobile phase B	Methano	Methanol Isoprop		ol	Isopropanol		Isopropanol		Isopropanol	
Caluma	ACQUITY UPC <sup>2</sup> 2-EP		ACQUITY UPC <sup>2</sup> 2-EP		ACQUITY UPC <sup>2</sup> 2-EP		GreenSep Nitro		ACQUITY UPC <sup>2</sup> 2-EP	
Column	(1.7 µm, 3.0 x 100 mm)		(1.7 µm, 3.0 x 100 mm)		(1.7 µm, 3.0 x 100 mm)		(1.8 µm, 3.0 x 100 mm)		(1.7 µm, 3.0 x 100 mm)	
Gradient	Time (min)	%B	Time (min)	%B	Time (min)	%B	Time (min)	%B	Time (min)	%B
	0	5	0	5	0	5	0	20	0	5
	4.5	45	2.50	25	2.50	25	2.5	50	2	25
	5	45	2.75	40	2.75	40	3	50	2.75	40
	5.25	5	3.25	40	3.25	40	3.25	20	3.25	40
	6	5	3.50	5	3.50	5	4	20	3.50	5
			4	5	4	5			4	5

Table 2. Key experimental parameters for UPC<sup>.2</sup>

	Figure 3E	3	Figure	5A	Figure 5B		
Instrument	Prep 100q SFC		SFC 80		SFC 80		
	MS-Directed System		Preparative System		Preparative System		
Flow rate (mL/min)	80		80		80		
Mobile phase A	CO <sub>2</sub>		CO <sub>2</sub>		CO <sub>2</sub>		
Mobile phase B	Isopropano	ol	Isopropa	anol	Ethanol		
Backpressure (psi)	1740		1740	)	1740		
Column	Viridis Silica 2-EP (5 µm, 19 x 150 mm)		Viridis Silica 2-EP (5 µm, 19 x 150 mm)		Nitro (5 µm, 21 x 150 mm)		
Temperature (°C)	40		40		40		
Sample diluent	Isopropanol		Isopropanol		Ethanol		
Injection volume (mL)	0.6		3		1		
Collection trigger	MS		UV		UV		
	Time (min)	%B	Time (min)	%B	Time (min)	%B	
	0	5	0	5	0	8	
Gradient	1	5	5	5	9	8	
	6.5	9	7	30			
	7	9	10	30			
	7.25	5	11	5			
	8	5	12	5			

Table 3. Key experimental parameters for preparative chromatography.

# **Results and Discussion**

#### Target compound extraction using SFE

Any solid-liquid extraction process, such as solvent extraction and SFE, is predominantly a solubility driven process. The process involves diffusion of the extracting solvent into the matrix, solubilization of the target analytes in the extracting solvent, diffusion of the target analytes in the extraction solvent, and transport of

the extracted analytes into a collection vessel.<sup>6</sup> Conventional polar extraction solvents, such as alcohols, often produce extracts comprised of mixtures of many polar and non-polar compounds. Supercritical CO<sub>2</sub>, on the other hand, is a highly lipophilic solvent. As a result, only relatively non-polar compounds are typically extracted by SFE using neat CO<sub>2</sub>. In the current study, the target compound is a terpene derivative with a nominal mass of 390.28 Da and a LogP of 3.0. The low molecular weight and the relatively low polarity make it an ideal candidate for extraction by SFE.

Figure 1 shows the UPC<sup>2</sup>-MS chromatograms of two extracts obtained by SFE (Figure 1A) and methanol extraction (Figure 1B) using a BEH 2-EP column. Since 2-EP is a polar stationary phase, the elution order of the compounds generally tracks their polarities; the later the elution, the more polar the compounds. While both extracts contain similar amount of the target compound, it is evident that SFE yielded a much simpler extract compared to methanol extraction. For the SFE extract, the peaks immediately after the target compound (1.20–1.75 min, blue rectangle) are much lower in intensity than those in the methanol extract. The peaks between 1.75–3.50 min (red rectangle) are only present in the methanol extract. Overall, the SFE extract is a much simpler mixture consisting of fewer polar components. The target compound was therefore enriched by SFE prior to chromatography. This makes the SFE extract ideal for large mass loading in prep chromatography and requires relatively low organic co-solvent (mobile phase B) composition to completely elute off the components in the extract; thereby shortening the total run time, reducing the solvent consumption, and increasing purification productivity. Detailed prep SFC experiments are described in a later section.



Figure 1. UPC<sup>2</sup>-MS chromatograms of the mixtures obtained by (A) SFE and (B) methanol extraction, using an ACQUITY UPC<sup>2</sup> 2-EP Column.

#### Conventional purification approach: MPLC + HPLC

One of most commonly used approaches in natural product isolation involves MPLC followed by HPLC. In the current case, the SFE extract first underwent a purification step by MPLC (results not shown), attaining the target compound of >97% purity (referred to as the MPLC fraction hereafter). The main remaining impurity has a nominal mass of 360.27 Da, and results from the demethoxylation of the target compound. The structural similarity between the target and impurity presented a challenge in RPLC purification. Figure 2A shows the UPLC-MS and UV chromatograms of the MPLC fraction. A baseline resolution between the target and the impurity was achieved using a 3.0 x 150 mm UPLC column, where the impurity was present as a sodium adduct with an m/z=383. The close elution of the two peaks, however, severely hampered the sample loadability in the ensuing RPLC purification. Figure 2B summarizes a loading study of the MPLC fraction on an analytical column (5 µm, 4.6 x 150 mm). The baseline resolution was only preserved with a 10-µL injection. With an 80-µL injection, the impurity peak completely merged into the target peak. In addition to the limited resolution, the elution order of the compounds also contributed to the low

purification productivity. With RPLC, the impurity elutes before the target compound. In the case where target and impurity are partially separated, such as the one with 40- $\mu$ L injection in Figure 2B, though it is still possible to obtain pure target compound by excluding the front of the target peak where the impurity co-elutes, such practice is generally inadvisable in prep chromatography as the front of a peak often accounts for a high percentage of the total peak. Based on the loading study performed on the analytical column, the maximum loading on a 19 x 150 mm semi-prep column without compromising yield or purity was projected to be 170  $\mu$ L. At ~20 mg/mL, this translates into a maximum loading of 3.4 mg/injection.





# Leveraging the orthogonality between RPLC and SFC for improved loading capacity: an MPLC + SFC approach

SFC offers an attractive alternative. SFC is generally considered a normal-phase chromatographic technique when a polar stationary phase, such as 2-EP, is used. As a result, the elution order often reverses that in RPLC using a non-polar  $C_{18}$  column. Figure 3A shows the UPC<sup>2</sup>-MS and UV chromatograms of the MPLC fraction using a BEH 2-EP column. Compared to Figure 2A, not only did the UPC<sup>2</sup> method provide a better resolution, the elution order of the target and the impurity also reversed. The chromatography was then scaled up to a 19 x 100 mm semi-prep column, and the resulting chromatogram is shown in Figure 3B. The resolution was well maintained with a 600-µL injection at 20 mg/mL. The total run time using SFC was 8 min compared to the 20-min run time using RPLC. By using prep SFC to replace prep RPLC, the overall productivity was increased by 9-fold: 2.5-fold from the reduced run time and 3.5-fold from the increased

sample loading.



Figure 3. (A) UPC<sup>2</sup>-MS and UV chromatograms of the MPLC fraction at 1 mg/mL; and (B) prep SFC-MS and UV chromatogram of the MPLC fraction at 20 mg/mL.

# Leveraging the orthogonality between different column chemistries in SFC for improved purification productivity: an SFC + SFC approach

Though the approach demonstrated in Figure 3 led to a notable improvement in productivity, the overall process still suffers from large solvent consumption, mainly due to the initial MPLC step. The target compound in the current study has a relatively low polarity. For this sample, a high percentage of organic solvent is required to elute the target compound in LC; hence, the large solvent consumption. In SFC, however, the lipophilic CO<sub>2</sub> is the main mobile phase that elutes the target compound, thus minimizing the use of organic solvents (mobile phase B). Moreover, the raw sample was extracted with neat CO<sub>2</sub> and is, therefore, inherently compatible with SFC.

There is a wide range of column chemistries available in SFC, with retention mechanisms encompassing polar interactions, hydrophobic interactions,  $\pi$ - $\pi$  interactions, and steric recognitions. With proper selection of column chemistries, SFC can offer orthogonal selectivity necessitated by the sample complexity intrinsic to natural product isolation. Figure 4 shows the UPC<sup>2</sup>-MS and UV chromatograms of the SFE extract using a BEH 2-EP (Figure 4A) and a nitro column (Figure 4B), respectively. While 2-EP columns typically render polar interactions between analytes and stationary phase, nitro columns often retain and separate analytes based on  $\pi$ - $\pi$  interactions. This kind of combination provides complementary separation around the target compound. As can be seen in Figure 4, using a 2-EP column, the target compound at m/z =391 is well separated from the impurity at m/z=361, but less separated from another later eluting impurity at m/z=239. In contrast, using a nitro column, the impurity at m/z=239 became an earlier eluting peak and was well separated from the target compound, but the impurity at m/z=361 co-eluted with the target compound.



Figure 4. UPC<sup>2</sup>-MS and UV chromatograms of the SFE extract using (A) a 2-EP; and (B) a nitro column.

# Leveraging the orthogonality between different column chemistries in SFC for improved purification productivity: an SFC + SFC approach

Though the approach demonstrated in Figure 3 led to a notable improvement in productivity, the overall process still suffers from large solvent consumption, mainly due to the initial MPLC step. The target compound in the current study has a relatively low polarity. For this sample, a high percentage of organic solvent is required to elute the target compound in LC; hence, the large solvent consumption. In SFC, however, the lipophilic  $CO_2$  is the main mobile phase that elutes the target compound, thus minimizing the use of organic solvents (mobile phase B). Moreover, the raw sample was extracted with neat  $CO_2$  and is, therefore, inherently compatible with SFC.

Based on the retention behavior illustrated in Figure 4, a two-step SFC purification strategy was implemented: using a 2-EP column to remove the main impurity with an m/z=361 followed by using a nitro

column to remove any remaining impurities after the first step, such as the one with an *m/z*=239. The resulting chromatograms are shown in Figure 5. The overall yield, defined as the weight of the purified pure target compound/the total weight of SFE extract taken for purification, was similar to those from the other two approaches: MPLC+HPLC and MPLC+SFC.



Figure 5. (A) SFC/UV chromatogram of the SFE extract at 133 mg/mL using a Viridis 2-EP column; and (B) SFC/UV chromatogram of the collected fraction from the Viridis 2-EP step on a nitro column.

Aliquots of the purified final product were analyzed by both UPC<sup>2</sup>-MS and UPLC-MS to ensure a true representation of the sample profile. The resulting chromatograms are shown Figure 6. Both impurities at m/z=361 and m/z=239 illustrated in Figure 4 were successfully removed. The results indicate that the final product has a purity >99% by UV.



Figure 6. Purity analysis of the final product by (A) UPLC-MS and (B) UPC<sup>2</sup>-MS.

The SFC purification process resulted in smaller fraction volumes compared to MPLC and HPLC. The SFC fractions were quickly dried under mild conditions, minimizing the possible compound loss due to thermal degradation associated with the post-purification dry-down process. Compared to LC, SFC offered an easier and faster compound recovery.

Process	Rate-limiting step	Productivity (g/24 hr)	Solvent	Solvent consumption (L/24 hr)	CO <sub>2</sub> use (kg/24hr)
MPLC+ HPLC*	HPLC	0.25	MeOH	95	N/A
MPLC+SFC**	SFC	2.25	MeOH/IPA	75	105
SFC+SFC	First step SFC	3.50	IPA/Ethanol	11	105

\*HPLC calculations were based on a 19 x 150 mm column.

\*\*SFC calculations were based on 19 x 150 mm columns.

Table 4. Comparison on productivity and solvent consumption of different purification processes.

## Conclusion

In this application note, we have demonstrated employing SFE and prep SFC to holistically improve the productivity in isolating a low-polarity, bioactive compound from a complex natural product extract. The SFE alleviated the sample complexity prior to analysis and purification, thereby improving sample loading and reducing solvent use in the ensuing chromatography. The SFE extract also lends itself well for SFC analysis and purification.

For the MPLC+HPLC purification route, the target compound and its demethoxylated derivative formed a critical pair in HPLC that limited the column loading and overall purification productivity. The same critical pair was better separated on a 2-EP column using SFC. The elution order of the pair was also altered, enabling an increased column loading. Overall, the MPLC+SFC route offered a 9-fold improvement in productivity. However, both routes still suffered from large solvent consumption because of the MPLC step. Finally, an SFC+SFC purification process was developed, leveraging the orthogonal selectivity between different column chemistries available in SFC. The SFC+SFC route not only led to a 16-fold improvement in productivity, but also a 90% reduction in solvent consumption. In addition, both SFE and SFC also provided an easy sample recovery under mild conditions that minimized potential compound loss due to thermal degradation associated with post-purification dry-down.

The supercritical fluid-based techniques, SFE and SFC, augment the conventional toolbox for natural product research by offering unique selectivity in both extraction and chromatography; and empower laboratories and manufacturers in pharmaceutical, traditional medicine, nutraceutical, dietary supplement, and cosmetic industries for more efficient and more costeffective natural product isolation and purification.

### References

- 1. Harvey AL, Strategies for discovering drugs from previously unexplored natural products. *Drug Discovery Today*, 2000; 5 (7):294–300.
- 2. Harvey AL, Natural products in drug discovery. *Drug Discovery Today*, 2008; 13 (19/20):894–901.
- 3. Li JWH, Vederas JC, Drug Discovery and natural products: end of an era or endless frontier? *Science*, 2009; 325 (10):161–165.
- 4. Harvey AL. Natural Products as a screening source, Curr. Opin. Chem. Biol., 2007, 11: 480–484.
- 5. Sarker SD, Latif Z, Gary AI, Natural product isolation: an overview, Natural Product Isolation,

2nd ed. Eds. Sarker SD, Latif Z, Gary AI, Humana Press Inc. Totowa, NJ. 2006, P1-25.

- 6. Sticher O, Natural product isolation, Nat. Prod. Rep., 2008, 25, 517–554.
- 7. Koehn FE, Carter G, The evolving role of natural products in drug discovery, *Nat. Rev. Drug Discov*. 2005, 4: 206–220.
- 8. Ramirez P, Garcia-Risco MR, Santoyo S, Senorans FJ, Ibanez E, Reglero G, *J. Pharmaceu. Biomed*. 2006, 41:1606–1613.

### **Featured Products**

ACQUITY UPC2 System <https://www.waters.com/134658367> ACQUITY UPLC H-Class PLUS System <https://www.waters.com/10138533> AutoPurification System <https://www.waters.com/10007147> Prep 100q SFC System <https://www.waters.com/10145748> ChromScope Software <https://www.waters.com/134647658> MassLynx MS Software <https://www.waters.com/513662> SFE Bio-Botanical Extraction System <https://www.waters.com/134826287> FractionLynx <https://www.waters.com/513795> MV-10 ASFE System <https://www.waters.com/134661549> SQ Detector 2 <https://www.waters.com/134631584> SFC 80q/200q and 350 Preparative Systems <https://www.waters.com/10145745>

720005168EN, October 2015

©2019 Waters Corporation. All Rights Reserved.