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1. Shion, H. *et al.* Enabling Routine and Reproducible Intact Mass Analysis When Data Integrity Matters. Waters Application Note, [720006472](#), 2019.
2. Ranbaduge, N. *et al.* Routine Peptide Mapping Analysis using the BioAccord System. Waters Tech Brief, [720006466](#), 2019.
3. Zhang, X. *et al.* Released N-linked Glycan Analysis Using the BioAccord System. Waters Tech Brief, [720006474](#), 2019.
4. Shion, H. *et al.* A Platform Method for the Molecular Mass Analysis of the Light Chains and Heavy Chains of Monoclonal Antibodies using the BioAccord System. Waters Tech Brief, [720006529](#), 2019.
5. Zhang, X. *et al.* Increasing Productivity and Confidence for N-linked Glycan Analysis of Biosimilars Using the BioAccord System. Waters Tech Brief, [720006545](#), 2019.
6. Shion, H. *et al.* Analysis of Antibody Drug Conjugates (ADCs) by Native Mass Spectrometry on the BioAccord System. Waters Tech Brief, [720006570](#), 2019.
7. Ranbaduge, N. *et al.* Establishing an Integrated Peptide Attribute Profiling and Monitoring Workflow for Improved Productivity and Confidence. Waters Application Note, [720006602](#), 2019.
8. Ippoliti, S. *et al.* Online IEX-MS of mAb Charge Variants Using a BioResolve SCX mAb Column, IonHance CX-MS pH Concentrates, and BioAccord System. Waters Application Note, [720006672](#), 2019.
9. Shion, H. *et al.* Evaluating BioAccord LC-MS System Performance for Deployment of a Compact, SmartMS-Enabled Biopharma Solution in Late Stage Development and QC Environment – Intact Mass Analysis. Waters Tech Brief, [720006722](#), 2020.

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©2019 Waters Corporation. All Rights Reserved.

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