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1. Shion, H.等人, *Enabling Routine and Reproducible Intact Mass Analysis When Data Integrity Matters* (《在数据完整性环境下进行可重现的完整蛋白质质量数例行分析》), 沃特世应用纪要, [720006472](#), 2019.
2. Ranbaduge, N.等人, *Routine Peptide Mapping Analysis using the BioAccord System* (《使用BioAccord系统进行常规肽图分析》), 沃特世技术简报, [720006466](#), 2019.
3. Zhang, X.等人, *Released N-linked Glycan Analysis Using the BioAccord System* (《使用BioAccord系统进行游离N-糖分析》), 沃特世技术简报, [720006474](#), 2019.
4. Shion, H.等人, *A Platform Method for the Molecular Mass Analysis of the Light Chains and Heavy Chains of Monoclonal Antibodies using the BioAccord System* (《使用BioAccord系统实现单克隆抗体轻链和重链质谱分析的平台方法》), 沃特世技术简报, [720006529](#), 2019.
5. Zhang, X.等人, *Increasing Productivity and Confidence for N-linked Glycan Analysis of Biosimilars Using the BioAccord System* (《使用BioAccord系统提高生物类似药N-糖分析的工作效率和可信度》), 沃特世技术简报, [720006545](#), 2019.
6. Shion, H.等人, *Analysis of Antibody Drug Conjugates (ADCs) by Native Mass Spectrometry on the BioAccord System* (《利用BioAccord系统通过非变性质谱法分析抗体偶联药物(ADC)》), 沃特世技术简报, [720006570](#), 2019.
7. Ranbaduge, N.等人, *Establishing an Integrated Peptide Attribute Profiling and Monitoring Workflow for Improved Productivity and Confidence* (《建立集成化的肽属性分析和监测工作流程以提高分析效率和结果可信度》), 沃特世应用纪要, [720006602](#), 2019.
8. Ippoliti, S.等人, *Online IEX-MS of mAb Charge Variants Using a BioResolve SCX mAb Column, IonHance CX-MS pH Concentrates, and BioAccord System* (《使用BioResolve SCX色谱柱、IonHance CX-MS pH浓缩液和BioAccord系统对mAb电荷异构体进行在线IEX-MS分析》), 沃特世应用纪要, [720006672](#), 2019.
9. Shion, H.等人, *Evaluating BioAccord LC-MS System Performance for Deployment of a Compact, SmartMS-Enabled Biopharma Solution in Late Stage Development and QC Environment - Intact Mass Analysis* (《评估SmartMS赋能的紧凑型生物制药解决方案BioAccord LC-MS系统部署于开发后期阶段和QC环境中用于完整蛋白质质量数分析的性能》), 沃特世技术简报, [720006722](#), 2020.

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[UNIFI科学信息系统 <https://www.waters.com/134801648>](https://www.waters.com/134801648)

[专为生物制药应用打造的BioAccord LC-MS系统 <https://www.waters.com/waters/nav.htm?cid=135005818>](https://www.waters.com/waters/nav.htm?cid=135005818)

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©2019 Waters Corporation. All Rights Reserved.

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