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
1. Johns Hopkins Resource Center.COVID-19 Case Tracker.<https://gisanddata.maps.arcgis.com/apps/opsdashboard/index.html#/bda7594740fd40299423467b48e>
2. Andersen K.G., Rambaut A., Lipkin W.I. *et al.* The proximal origin of SARS-CoV-2. *Nat Med* 2020; 26:450–452.
3. Tortorici AM, Veasler D. Chapter Four–Structural insights into coronavirus entry, Editor(s): Félix A. Rey, Advances in Virus Research, Academic Press, 2019, 105:93-116, ISSN 0065-3527, ISBN 9780128184561.
4. Walls, AC., Park YJ, Tortorici A, Wall A, McGuire AT, Veasler D. Structure, Function, and Antigenicity of the SARS-CoV-2 Spike Glycoprotein, *Cell*, 2020;181(2):281-292.e6, ISSN 0092–8674.
5. Clark GF. The Role of Glycans in Immune Evasion: The Human Fetoembryonic Defence System Hypothesis revisited, *Molecular Human Reproduction*. 2014 Mar; 20(3):185–199.
6. Chang D, Zaia J. Why Glycosylation Matters in Building a Better Flu Vaccine. *Molecular & Cellular Proteomics* 2019 Dec; 18 (12):2348–2358.
7. Khatri K, Klein JA, White MR, Grant OC, Leymarie N, Woods RJ, Hartshorn KL, Zaia J. Integrated Omics and Computational Glycobiology Reveal Structural Basis for Influenza A Virus Glycan Microheterogeneity and Host Interactions. *Molecular & Cellular Proteomics* 2016 Jun; 15 (6):1895–1912.
8. Kumar S, Maurya VK, Prasad AK, Bhatt MLB, Saxena SK. Structural, Glycosylation and Antigenic Variation between 2019 Novel Coronavirus (2019-nCoV) and SARS Coronavirus (SARS-CoV). *Virusdisease*. 2020;31(1):13–21.
9. Watanabe Y, Allen JD, Wrapp D, McLellan JS, Crispin M. Site-Specific Glycan Analysis of the SARS-CoV-2 Spike. *Science*. 2020;eabb9983.
10. Shajahan A, Supekar NT, Gleinich AS, Azadi P. Deducing the N- and O- Glycosylation Profile of the Spike Protein of Novel Coronavirus SARS-CoV-2. *Glycobiology*. 2020; cwaa042. <https://doi.org/10.1093/glycob/cwaa042>.
11. Sanda M, Morrison L, and Goldman R. N- and O-Glycosylation of the SARS-CoV-2 Spike Protein.

Anal. Chem. 2021.

12. Guttman M and Lee KK. Site-Specific Mapping of Sialic Acid Linkage Isomers by Ion Mobility Spectrometry. *Anal. Chem.* 2016; 88(10):5212–5217.
13. Jin C, Harvey DJ, Struwe WB, and Karlsson NG. Separation of Isomeric O-Glycans by Ion Mobility and Liquid Chromatography-Mass Spectrometry. *Anal. Chem.* 2019. 91(16): 10604–10613.
14. Harvey DJ, Seabright GE, Vasiljevic S, Crispin M, and Struwe WB. Isomer Information from Ion Mobility Separation of High-Mannos Glycan Fragments. *J. Am. Soc. Mass. Spectrom.* 2018; 29(5): 972–988.
15. Barroso A, Giménez A, Konijnenberg A, Sancho J, Sanz-Nebot V, Sobott F. Evaluation of Ion Mobility for the Separation of Glycoconjugate Isomers Due to Different Types of Sialic Acid Linkage, at the Intact Glycoprotein, Glycopeptide and Glycan Level. *J. Proteom.* 2018. 173:22–31.
16. Hofmann J, Hahm HS, Seeberger PH, and Pagel K. Identification of Carbohydrate Anomers using Ion Mobility-Mass Spectrometry. *Nature.* 2015. 526:241-244.
17. Depland AD, Renois-Predelus G, Schindler B, Compagnon I. Identification of Sialic Acid Linkage Isomers in Glycans using Coupled InfraRed Multiple Photon Dissociation (IRMPD) Spectroscopy and Mass Spectrometry. *Int.J. Mass Spectrom.* 2018; 434:64–69.

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SELECT SERIES Cyclic IMS <<https://www.waters.com/135021297>>

ACQUITY UPLC M-Class系统 <<https://www.waters.com/134776759>>

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© 2021 Waters Corporation. All Rights Reserved.