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1. Senanayake SPJN, Shahidi F. Modification of Fats and Oils via Chemical and Enzymatic Methods. Bailey's Industrial Oil and Fat Products. 2005.
2. Srivastava Y, Semwal AD, Sajeevkumar VA, Sharma GK. Melting, Crystallization and Storage Stability of Virgin Coconut Oil and Its Blends by Differential Scanning Calorimetry (DSC) And Fourier Transform Infrared Spectroscopy (FTIR). *J Food Sci Technol*. 2017 Jan; 54(1):45–54.
3. Bro R. Multivariate Calibration: What Is in Chemometrics for the Analytical Chemist? *Anal Chim Acta*. 2003; 500(1):185–94.
4. Forina M, Lanteri S, Casale M. Multivariate Calibration. *J Chromatogr A*. 2007; 1158(1–2):61–93.
5. Ferreira JA, Santos JM, Breitzkreitz MC, Ferreira JMS, Lins PMP, Farias SC, *et al*. Characterization of the Lipid Profile From Coconut (*Cocos nucifera* L.) Oil of Different Varieties by Electrospray Ionization Mass Spectrometry Associated With Principal Component Analysis and Independent Component Analysis. *Food Res Int*. 2019; 123:189–97.
6. Bataglion GA, da Silva FMA, Santos JM, Dos Santos FN, Barcia MT, de Lourenço CC, *et al*. Comprehensive Characterization of Lipids From Amazonian Vegetable Oils by Mass Spectrometry Techniques. *Food Res Int*. 2014 Oct; 64:472–81.
7. Pizzo JS, Galuch MB, Manin LP, Santos PDS, Zappiello CD, Santos Junior O, *et al*. Direct Infusion Electrospray Ionisation Mass Spectrometry Applied in the Detection Of Adulteration of Coconut Oil With Palm Kernel Oil. *Food Addit Contam Part A, Chem Anal Control Expo risk Assess*. 2019 Nov; 36(11):1597–604.

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- [RADIAN ASAP直接质谱检测器 <https://www.waters.com/waters/nav.htm?cid=135073413>](https://www.waters.com/waters/nav.htm?cid=135073413)
- [LiveID软件 <https://www.waters.com/134939519>](https://www.waters.com/134939519)

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© 2022 Waters Corporation. All Rights Reserved.