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---

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1. Holz G, Chepurny O. Glucagon-Like Peptide-1 Synthetic Analogs: New Therapeutic Agents for Use in the Treatment of Diabetes Mellitus. *Current Medicinal Chemistry*. 2003 Nov 1;10(22):2471–83. <<http://online.uspnf.com>>
2. Magda Sara Wojtara, Mazumder A, Syeda Y, Nikodem Mozgala. Glucagon-Like Peptide-1 Receptor Agonists for Chronic Weight Management. *Advances in medicine*. 2023 Sep 20;2023(Article ID 9946924):1–7.
3. United States Pharmacopeia. Exenatide Injection USP Monograph. [online.uspnf.com](http://online.uspnf.com) <<http://online.uspnf.com>> .2019.
4. United States Pharmacopeia. Glucagon USP Monograph. [online.uspnf.com](http://online.uspnf.com) <<http://online.uspnf.com>> .2022.
5. Kumar Kuna A, Ganapathi S, Radha G. A Novel RP- HPLC Method Development and Forced Degradation Studies for Semaglutide in Active Pharmaceutical Ingredients and Pharmaceutical Dosage Form. *International Journal of Research in Pharmaceutical Sciences*. 2019 Apr 15;10(2):865–73.
6. Pedaprolu JN, Bonthu M, Vatchavai B, Kamatham S, Kolli S, Kapuganti AN. A New Stability-Indicating and Validated RP-HPLC Method for the Estimation of Liraglutide in Bulk and Pharmaceutical Dosage Forms. *Eurasian Journal of Analytical Chemistry*. 2016 Dec 8;12(2):31–44.
7. Pang E. Assessing Immunogenicity Risk of Peptides: the Synthetic Peptide Guidance and PSGs SBIA 2022: Advancing Generic Drug Development: Translating Science to Approval Day (1), Session 1A: (Peptide Immunogenicity Risk and Impurity Assessment Considerations) [Internet]. U.S. Food and Drug Administration; 2020 [cited 2023 Nov 6]. Available from: <https://www.fda.gov/media/166571/download> <<https://www.fda.gov/media/166571/download>> .
8. Birdsall R, Kellet J, Ippoliti S, Qing Yu Y. Increasing Recovery and Chromatographic Performance of “Acidic” Peptides Using Waters ACQUITY Premier Solution [Internet]. [www.waters.com](http://www.waters.com). Waters Corporation; 2021. Available from: [720007173](https://www.waters.com/720007173). March 2021.
9. Bigos P, Birdsall R, Walter T. MaxPeak Premier Solutions: Improving Consumer Safety Through Innovative

Science [Internet].www.waters.com.Waters Corporation; 2021.Available from: [720008054 <https://www.waters.com/nextgen/global/library/library-details.html?documentid=720008054>](https://www.waters.com/nextgen/global/library/library-details.html?documentid=720008054) .

10. MaxPeak Premier Peptide Reversed-Phase Column and Method Screening Kit: Practical Steps in Developing Robust Peptide Separations Available from [720008131 <https://www.waters.com/webassets/cms/support/docs/720008131en.pdf>](https://www.waters.com/webassets/cms/support/docs/720008131en.pdf) .
11. Bao Z, Cheng YC, Luo MZ, Zhang JY.Comparison of the Purity and Impurity of Glucagon-for-Injection Products under Various Stability Conditions.Scientia Pharmaceutica [Internet].2022 Jun 1 [cited 2023 Nov 10];90(2):32.Available from: [https://www.mdpi.com/2218-0532/90/2/32 <https://www.mdpi.com/2218-0532/90/2/32>](https://www.mdpi.com/2218-0532/90/2/32) .

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© 2025 Waters Corporation. All Rights Reserved.

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