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1. Hansen K J, Clemen L A, Ellefson M E., Johnson H O. Compound-Specific Quantitative Characterisation of Organic Fluorochemicals in Biological Matrices.*Environ.Sci.Technol.* 35, 766–770, 2001.
2. Giesy J P and Kannan K Global Distribution of Perfluorooctane Sulfonate in Wildlife. *Environ.Sci.Technol.* 35: 1339–1342, 2001.
3. Kelly B C, Ikononou M G, Blair J D, Surridge B, Hoover D, Grace R, Gobas F A P C. Perfluoroalkyl Contaminants in an Arctic Marine Food Web: Trophic Magnification and Wildlife Exposure. *Environ.Sci.Technol.* 43 (11), 4037–4043, 2009.
4. EFSA report.Results of the Monitoring of Perfluoroalkylated Substances in Food in the Period 2000–2009.February 2011.
5. United States Environmental Protection Agency.Analysis of Per and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS.4th Draft method 1633, July 2023.
6. Organtini K, Rosnack K, Hancock P. 根据EPA 1633分析全氟烷基和多氟烷基化合物(PFAS) 第1部分：建立和评估方法.沃特世应用纪要. [720008117ZH](#), 2023年.
7. Khoury-Hollins H., Reid L., Adams S. 使用高分辨率飞行时间质谱仪Xevo™ G3 QTof直接进样方法筛查饮用水中的全氟烷基和多氟烷基化合物(PFAS).沃特世应用纪要.[720008198ZH](#), 2024年.
8. Commission Recommendation (EU) 2022/1431, on the Monitoring of Perfluoroalkyl Substances in Food.L 221/105, 26.8.2022.
9. Travis D Evidence That There are Two Forms of Fluoride in Human Serum.*Nature* 214: 1050–1051, 1968.
10. Lau C, Perfluorinated Compounds: An Overview.In: DeWitt, J.C. (Ed.), Toxicological Effects of Perfluoroalkyl and Polyfluoroalkyl Substances.Springer, pp.1–21, 2015.
11. Dickman R A, Aga D S. A Review of Recent Studies on Toxicity, Sequestration, and Degradation of Per- and Polyfluoroalkyl Substances (PFAS) *Journal of Hazardous Material.* 436: 129120, 2022.
12. Trudel D, Horowitz L, Wormuth M, Scheringer M, Cousins IT, Hungerbühler K. Estimating Consumer Exposure to PFOS and PFOA.*Risk Anal.* 28:251–69, 2008.

13. De Silva A O, Armitage J M, Bruton T A, Dassuncao C, Heiger-Bernays W, Hu X C, Kärrman A, Kelly B, Ng C, Robuck A, Sun M, Webster T F, Sunderland E M. PFAS Exposure Pathways for Humans and Wildlife: A Synthesis of Current Knowledge and Key Gaps in Understanding Environ *Toxicol Chem.* 40: 631–657, 2021.
14. Schrenk D, Bignami M, et al, EFSA Panel on Contaminants in the Food Chain (EFSA CONTAM Panel), Risk to Human Health Related to the Presence of Perfluoroalkyl Substances in Food.EFS2;18(9) 2020.
15. EURL for Halogenated POPs in Feed and Food (2022): Guidance Document on Analytical Parameters for the Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Food and Feed, version 1.2 of 11 May 2022.在线获取(https://eurl-pops.eu/core-working-groups#_pfas <https://eurl-pops.eu/core-working-groups#_pfas>)
16. Adams S, Dreolin N, Organtini K L, Hancock P. 以鱼类参比样品作为案例研究, 使用高灵敏度LC-MS/MS方法分析动物产品中的全氟烷基和多氟烷基化合物(PFAS).沃特世应用纪要.720008108ZH, 2023年.
17. Twohig M, Fujimoto G, Mohan A, Organtini K L, Rosnack K J, Hird S. 环境样品中全氟烷基和多氟烷基化合物(PFAS)的非靶向分析方法.沃特世应用纪要.720007184ZH, 2021年.
18. U.S. Environmental Protection Agency.CompTox Chemicals Dashboard.
<https://comptox.epa.gov/dashboard/chemical-lists/pfasmaster> <
<https://comptox.epa.gov/dashboard/chemical-lists/pfasmaster>> accessed November 2023.
19. Suspect List of Possible Per- and Polyfluoroalkyl Substances (PFAS).
<https://data.nist.gov/od/id/mds2-2387> <<https://data.nist.gov/od/id/mds2-2387>> Accessed March 2023.
20. Bugsel B, Zweigle J, Zwiener C. Nontarget Screening Strategies for PFAS Prioritization and Identification by High Resolution Mass Spectrometry: A Review.*Trends in Environmental Analytical Chemistry.* 40: e00216, 2023.
21. Charbonnet J A, McDonough C A, Xiao F, Schwichtenberg T, Cao D, Kaserzon S, Thomas K V, Dewapriya P, Place B J, Schymanski E L, Field J A, Helbling D E, Higgins C P. Communicating Confidence of Per- and Polyfluoroalkyl Substance Identification via High-Resolution Mass Spectrometry.*Environ Sci Technol Lett.* 9: 473–481, 2022.

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© 2025 Waters Corporation. All Rights Reserved.

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