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1. Shion H, Boyce P, Berger SJ, Yu YQ. INTACT Mass- a Versatile waters_connect Application for Rapid Mass Confirmation and Purity Assessment of Biotherapeutics. Waters Application Note [720007547](https://www.waters.com/nextgen/global/library/application-notes/2022/intact-mass-a-versatile-waters-connect-application-for-rapid-mass-confirmation-and-purity-assessment-of-biotherapeutics.html) <<https://www.waters.com/nextgen/global/library/application-notes/2022/intact-mass-a-versatile-waters-connect-application-for-rapid-mass-confirmation-and-purity-assessment-of-biotherapeutics.html>> . February 2022.
 2. Doneanu CE, Boyce P, Shion H, Fredette J, Berger SJ, Gastall H, Yu YQ. LC-MS Analysis of siRNA, Single Guide RNA and Impurities Using the BioAccord System with ACQUITY Premier and New Automated INTACT Mass Application. Waters Application Note [720007546](https://www.waters.com/nextgen/global/library/application-notes/2022/lc-ms-analysis-of-sirna-single-guide-rna-and-impurities-using-the-bioaccord-system-with-acquity-premier-and-new-automated-intact-mass-application.html) <<https://www.waters.com/nextgen/global/library/application-notes/2022/lc-ms-analysis-of-sirna-single-guide-rna-and-impurities-using-the-bioaccord-system-with-acquity-premier-and-new-automated-intact-mass-application.html>> . April 2022.
 3. D' Esposito RJ, Doneanu CE, Gastall H, Berger SJ, Yu YQ. RNA CQA Analysis using the BioAccord LC-MS System and INTACT Mass waters_connect Application. Waters Application Note [720008130](https://www.waters.com/nextgen/global/library/application-notes/2023/rna-cqa-analysis-using-the-bioaccord-lc-ms-system-and-intact-mass-waters-connect-application.html) <<https://www.waters.com/nextgen/global/library/application-notes/2023/rna-cqa-analysis-using-the-bioaccord-lc-ms-system-and-intact-mass-waters-connect-application.html>> . November 2023.
 4. Ranbaduge N, Shion H, Yu YQ. Streamlined LC-MS Analysis of Stress Induced Impurities of a Synthetic Peptide using the BioAccord System and the waters_connect INTACT Mass Application. Waters Application Note [720007752](https://www.waters.com/nextgen/global/library/application-notes/2022/streamlined-lc-ms-analysis-of-stress-induced-impurities-of-a-synthetic-peptide-using-the-bioaccord-system-and-the-waters-connect-intact-mass-application.html) <<https://www.waters.com/nextgen/global/library/application-notes/2022/streamlined-lc-ms-analysis-of-stress-induced-impurities-of-a-synthetic-peptide-using-the-bioaccord-system-and-the-waters-connect-intact-mass-application.html>> . October 2022.
-

5. Riccardi F, Bo MD, Macor P, Toffoli G. A comprehensive overview on antibody-drug conjugates: from the conceptualization to cancer therapy. *Front. Pharmacol.* 14. 2023.
6. Shion H, Yu YQ, Chen W. Analysis of Antibody Drug Conjugates (ADCs) by Native Mass Spectrometry on the BioAccord System. Waters Application Note 720006570 <<https://www.waters.com/nextgen/global/library/application-notes/2019/analysis-of-antibody-drug-conjugates-adcs-by-native-mass-spectrometry-on-the-bioaccord-system.html>> . May 2019.
7. Martín M, Pandiaella A, Vargas-Castrillón E, Díaz-Rodríguez E, Iglesias-Hernangómez T, Cano CM, Fernández-Cuesta I, Winkow E, Perelló MF. Trastuzumab deruxtecan in breast cancer. *Critical Reviews in Oncology/Hematology*. 2024; 198: 104355. (<http://doi.org/10.1015/j.critrevonc.2025.104355> <<https://www.sciencedirect.com/science/article/pii/S1040842824000982?via%3Dihub>>)
8. Chen L, Wang L, Shion H, Yu C, Yu YQ, Zhu L, Li M, Chen W, Gao K. In-depth structural characterization of Kadcyla (ado-trastuzumab emtansine) and its biosimilar candidate. *MAbs*. 2016; 8(7): 1210–1223.

UNIFI、BioAccord和waters_connect是沃特世科技公司的商标。ENHERTU (fam-trastuzumab-deruxtecan-nxki)是Daiichi-Sankyo Company, Limited的商标。KADCYLA (ado-trastuzumab-emtansine)是Genentech, Inc.的商标。

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BioAccord LC-MS系统 <<https://www.waters.com/nextgen/global/products/mass-spectrometry/mass-spectrometry-systems/bioaccord-lc-ms-system.html>>

waters_connect信息学平台 <https://www.waters.com/nextgen/global/products/informatics-and-software/waters_connect/waters-connect-informatics-platform.html>

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© 2025 Waters Corporation. All Rights Reserved.

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