

Purifying Oligonucleotides

High Efficiency Preparative Chromatography
to Improve Yields and Turnaround Time

Understanding Oligonucleotide Purification

Oligonucleotide Structures and Modifications

Synthetic oligonucleotides are short nucleic acid sequences. They are chemically synthesized and typically 20–30 residues long but can sometimes be up to 100 nucleotides in length. They can be used both therapeutically and diagnostically. For accurate base pairing and recognition, it is critical that they are made with the correct sequence and free of impurities.

These so-called “oligos” are composed of deoxyribose or ribose sugars connected to nucleobases, such as adenine (A), thymine (T), cytosine (C) and guanine (G) or uracil (U). Many therapeutic oligos and oligo reagents are modified. These modifications are fundamental for improvements in nuclease stability, efficacy of delivery and detection sensitivity. Among several chemical modifications and structural adjustments to oligonucleotides, it is Phosphorothioate (PS) backbone modification, Locked Nucleic Acids (LNA), 2'-O-methyl and 2'-fluoro, capping, PEGylation etc. that are common.

Need for Purification

Identification

- Rapid identification of synthetic modified oligos
- Chemical modifications like PEGylation, conjugation to other chemical moieties etc.

Purification

- Synthesis- and process-related impurities
- Challenging base modifications

Purity Requirements

- Purity range for synthetic oligonucleotides depends on the goal of purification
- In vivo screening (>95%)
- In vitro screening (≥ 85%)

Essential Strategies

Reversed Phase Chromatography (RP) is effective for separating oligonucleotides with different length and sequences, and conjugates such as GalNAc & lipid conjugated moieties. **Ion-pairing (IP-RP)** is the go-to choice for high resolution separations of modified and unmodified oligos especially in complex mixtures.

Anion-Exchange Chromatography (AEX) is useful for separation of oligonucleotides based on their charge differences, such as single-stranded vs. double-stranded DNA. Oligonucleotides with high GC content are strong candidates for purification using AEX chromatography.

Size Exclusion chromatography (SEC) can also be used as a desalting step after IP-RP. It can also be used to size select different lengths and conformers. Impurities generated during synthesis of larger Oligos, usually are shorter than the target oligo thus can be separated by size exclusion chromatography.


Oligonucleotide Structures



A single stranded oligonucleotide, like an antisense oligonucleotide (ASO), PCR primer, or fluorescent probe



Double stranded oligonucleotide representative of a small interfering RNA (siRNA)




A long conformationally complex oligonucleotide denoting a CRISPR single guide RNA.

A conformationally example of an aptamer

Oligonucleotide Modifications

The following provides select modifications as reported in *Pharmaceutics*. 2022 Feb; 14(2): 260. CC-BY.


Backbone Modifications



Sugar Modifications





Conjugates



Tools and Strategies for Successful Purification

Successfully Applying IP-RP Purifications

High resolution purification of oligos from their process and product related impurities is most easily established using ion pair reversed phase LC. Since oligos are hydrophilic, they will not inherently retain to the hydrophobic alkyl chains or polymeric sorbents used as stationary phases for reversed phase columns. Consequently, retention is achieved through the use of an ion pairing reagent, typically an alkyl amine. The protonated amine will ion pair to the negatively charged backbone of the oligonucleotide. Preparative flow HPLC with flow rate capabilities up to 150 mL/min are used to readily apply this form of chromatography to columns with inner diameters ranging from 10 to 100 mm wide. Mass directed fraction collection can also be applied with specialized AutoPurification Systems.



Effect of column packing procedures on bed density, as reported in Waters Application Note 720001939.



The LC Prep AutoPurification System provides the flexibility of high-throughput parallel runs for selective mass-directed fraction collection from hundreds of samples.

Stationary Phase Considerations

Waters batch certified Bridged Ethyl Hybrid (BEH™) Technology for oligo separations ensures extra performance, reproducibility and scale-up for purification. High temperature and high pH are prevalent strategies for minimizing secondary interactions in oligonucleotide separation and purification. BEH Particle Technology withstands harsh conditions such as high pH, temp and additives that usually are required for Oligonucleotide work.

Oligonucleotides and nucleic acids vary in size and structure, ranging from just a few bases to thousands. Selecting the appropriate pore size is essential for efficient mass transfer of the oligonucleotide into the pore structure. The chart to the left helps visualize the selection of pore size for oligonucleotides with different lengths.

Optimal Bed Density (OBD) Packing

For prep LC columns with smaller aspect ratio (length/diameter), high pressure slurry packing with small particles often fails to achieve the ideal bed density in well designed analytical columns. Excessive axial compression at the inlet can cause particle breakage, channeling and reduced local bed permeability. Carefully designed OBD™ Columns deliver improvements in column lifetime, efficiency, peak shape, and low back pressure.

Challenges to Address

- Loading capacity
- Low resolution
- Unspecific binding
- Solubility
- Low purity
- Irreproducibility batch to batch and column to column
- Long elution time
- Need for a different selectivities
- Struggles with column secondary interactions and passivation
- Few starting point methods for newly emerging conjugates and modifications


Method Considerations

Mobile Phase Options

Ion-pairing reagents in the mobile phase increases the hydrophobicity of charged components. The reagent interacts with the charged groups in the oligonucleotide, reducing their impact on overall hydrophobicity. The type and concentration of the ion-pairing reagent influence retention mechanism and selectivity. Acetonitrile and methanol are most frequently applied eluents, though ethanol could also be used in some situations.

Depending on the type of oligonucleotides a range of alkylamines with different hydrophobicity can be used. To selectively separate oligos based on their length, more hydrophobic alkylamines such as TPA (tripropylamine), HA (hexylamine) or DBA (dibutylamine) can be considered. TEAA (triethylammonium acetate) or HAA (hexylammonium acetate) are volatile and can be removed from a collected fraction by evaporation or lyophilization. Triethylammonium bicarbonate presents yet another interesting volatile mobile phase option.


Example Modifiers and Ion Pairing Reagents



Scalability

Scalability while sustaining the chromatography profile is crucial for separation of critical impurities from a target compound. During oligonucleotide synthesis, deletion sequences often occur, resulting in the presence of full-length products alongside impurities, such as n-1,

Analytical Run




Analytical	
Column	4.6 x 50 mm, 2.5 μ m XBridge Oligo C18
Flow Rate	0.4 mL/min
Injection Vol.	1 μ L
Mass Load	10 μ g

42X

Preparatory	
Column	30 x 50 mm, 2.5 μ m XBridge Oligo C18
Flow Rate	17 mL/min
Injection Vol.	42 μ L
Mass Load	420 μ g

Gradient Calculator	Method Optimization	Dead Volume
Waters Gradient Calculator		
Please Select Language		
Method 1	Method 2	
System (Dead Volume (mL))	0.000	0.000
Column ID (mm)	4.6	30.0
Column Length (mm)	50	50

Preparative Run



HPLC w/ TEAA modifier
XBridge Oligonucleotide
BEH 130Å C₁₈ 2.5 μ m 4.6 x 50 mm Column
Flow rate: 0.8 ml/min
Column temp: Ambient

Mass Load: 420 μ g
XBridge Oligonucleotide BEH C₁₈
130Å 2.5 μ m OBD Prep 30 x 50 mm Column

Time-Based Fractionation
Productivity = ~1mg of purified material / hour

Loading Capacity

There are several important parameters to consider when it comes to optimizing your loading capacity. Column efficiency, pore size and particle size play significant roles. Larger pore size columns allow for better mass transfer and sharper peak shapes of larger analytes but having reduced surface area means they can exhibit lower loading capacity. Most importantly, a wide range of column IDs are available to optimize your purification runs.

System Considerations

Your choice of preparative LC system depends on the goal, and quantity of sample to be purified. Options include a semi-prep, lab scale prep, pilot or process scale instrument. Waters provides semi-prep and lab-scale preparative LC systems that allow auto purification configurations and ensure analytical scale up to prep scale chromatography with optimized fluidic designs, minimized dispersion and low carryover.

Detector Directed Fractionation and Impurity Analysis

Mass spectrometry (MS) directed fraction collection in preparative runs, enables maximization of targeted compound recovery. Diode Array/ultraviolet-visible directed fraction collection can also be applied, albeit for less selective triggering. A streamlined and compatible fraction collection strategy improves overall lab productivity.

Column Hardware and MaxPeak Premier High Performance Surfaces


Oligonucleotides phosphodiester groups can adsorb to the metals found in column bodies and frits. Using low adsorption column hardware can improve peak shape by minimizing chelation, leading to more reliable retention and more consistent chromatographic parameters. XBridge Premier Oligonucleotide BEH C₁₈ OBD Prep columns are manufactured with MaxPeak™ HPS vapor deposited components to enhance peak shape and improve recovery.

Cost-effective purification

Most older generation prep columns are made with particles ranging from 7 to 50 µm in diameter primarily due to limitations in column technology as well and instrument performance at the time. Prep columns with smaller particle size and optimized packing technology enhance the separation of impurities from unmodified, modified and conjugated oligonucleotides, providing higher reproducibility.

XBridge Oligonucleotide BEH C₁₈ 5 µm Prep Columns combined with OBD Technology offer extended column lifetimes and reduced the need for multiple purification passes.

ID (mm)	Loading (mg)
4.6	≤4
10	≥5–15
19	≥15–<200
30	≥200–<350
50	≥350–<600




Mass loading studies performed on a 4.6 mm ID column, as reported in Waters Application Note 720002602EN.




Up to 5 mL/min




Up to 10 mL/min




Application Notes




[Increasing the Productivity of Oligonucleotide Purification through Column Scaling and Method Optimization](#)




[Analysis of siRNA Duplexes at Non-Denaturing UPLC Conditions Using MaxPeak Premier Column Technology](#)




[CRISPR Single Guide RNA Characterization by IP- RP-LC-MS with a Premier Oligonucleotide BEH 300 Å C₁₈ Column](#)




[Improved Peak Recovery and Peak Shape of Oligonucleotides Using Waters Premier Columns](#)




[Best Practices for Oligonucleotide Analysis Using Ion-Pair Reversed-Phase \(IP-RP\) Liquid Chromatography – Columns and Chemistries](#)




[Improved Chromatographic Analysis of Oligonucleotides with ACQUITY Premier Oligonucleotide BEH C₁₈ Columns](#)




[Enhanced Resolution for Longer Oligonucleotide Analytes With a MaxPeak™ Premier Oligonucleotide BEH C₁₈ 300 Å](#)




[Analysis of Oligonucleotide Impurities on the BioAccord System with ACQUITY Premier](#)




[Evaluating the Tools for Improving Purification Throughput](#)



[Mass Confirmation and Purity Analysis of PCR Primer and Probe Reagents Using High Sensitivity Ion-pairing LC-MS](#)




[Optimum Bed Density \[OBD™\] columns: Enabling Technology for Laboratory-Scale](#)




[Analysis of siRNA Drugs at Denaturing UPLC Conditions Using MaxPeak Premier Column Technology](#)

Waters Journal Articles and Other Key Resources



[The impact of low adsorption surfaces for the analysis of DNA and RNA oligonucleotides.](#)

J Chrom A 2022, 1677, 463324




[Using Hybrid Organic-Inorganic Surface Technology to Mitigate Analyte Interactions with Metal Surfaces in UHPLC](#)

A. Chem. 2021, 93, 14, 5773–5781



www.waters.com/prepcalculator



[Analysis of siRNA with Denaturing and Non-Denaturing Ion-Pair Reversed-Phase Liquid Chromatography Methods](#)

LCGC North America, 2023, 41, 2, 60–66

Ordering Information

The columns, standards, and reagents that can help you purify your oligonucleotides.

XBridge™ Oligonucleotide BEH C₁₈ OBD Prep 5 µm Columns

Diameter		C. Guard	50 mm	100 mm	150 mm	C. Guard	50 mm	100 mm	150 mm
	Pore Size	130 Å				300 Å			
Standard Column	10 mm	186011150 ¹	186011151	186011152	186011153	186011168 ¹	186011169	186011170	186011171
	19 mm	186011154 ²	186011155	186011156	186011157	186011172 ²	186011173	186011174	186011175
	30 mm	186011158 ³	186011159	186011160	186011161	186011176 ³	186011177	186011178	186011179
	50 mm	–	186011162	186011163	186011164	–	186011180	186011181	186011182

¹ Requires Prep Guard Holder P/N 289000779

² Requires Prep Guard Holder P/N 186008745

³ Requires Prep Guard Holder P/N 186006912



XBridge Premier Oligonucleotide BEH C₁₈ OBD Prep 5 µm Columns

Diameter		C. Guard	50 mm	100 mm	150 mm	C. Guard	50 mm	100 mm	150 mm
	Pore Size	130 Å				300 Å			
Standard Column	10 mm	186011205 ¹	186011206	186011207	186011208	186011216 ¹	186011217	186011218	186011219

¹ Requires MaxPeak™ Prep Guard Holder P/N 186011548

XBridge Oligonucleotide BEH C₁₈ OBD Prep 130 Å 2.5 µm Columns

Diameter		50 mm
Standard Column	10 mm	186008212
	19 mm	186008962
	30 mm	186008963
	50 mm	186008964

For custom made products, contact a Waters representative.

XBridge Oligonucleotide BEH C₁₈ Analytical Columns for Method Development Screening

Diameter		50 mm	100 mm
		130 Å	300 Å
Standard Column 2.5 µm	4.6 mm	186003953	
Standard Column 5 µm	4.6 mm	186011165	186011183

For more information about oligonucleotide analytical testing, please visit [waters.com/Oligo](https://www.waters.com/Oligo)



waters.com/GTx

waters.com/Oligo

waters.com/OligoRPLC

waters.com/OligoWorks

waters.com/GTxResolveHILIC

waters.com/GTxStandards

Waters™

Waters, MaxPeak, OligoWorks, XBridge, ACQUITY, BEH, and XSelect are trademarks of Waters Technologies Corporation.
All other trademarks are the property of their respective owners.

©2025 Waters Corporation. Aug 25-14347 720008634EN Rev. C

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
waters.com