



HPLC Troubleshooting: A System Monitoring Approach



Step 5

Eliminate Causes that are Not Consistent
If there is more than one failure indicator, eliminate the causes that are not present in all lists, because they do not explain all the failure indicators. **Only causes that are consistent across applicable lists should remain.**

Decrease in Peak Area
The new peak area is stable but does not match the benchmark data.

- Evaporative loss of analytes
- Degradation of analytes
- Leak in system (between sample injector and detector)
- Change in split factor (for systems with a splitter only)
- Sample preparation error
- Sample volume low in the vial
- Weak needle wash empty or low
- Leak in sample fluidics
- Bubble in sample fluidics line
- Weak needle wash not compatible with sample
- Loss of detector sensitivity
- Injection Needle damaged
- Injection needle not drawing enough liquid

Step 6

Eliminate Causes that Do Not Explain Other Known Information
Eliminate additional causes that don't make sense. In this example, we can eliminate "Loss of detector sensitivity" because it would not explain why only two of the analytes are affected. We can also eliminate "sample preparation error" because we are using a reference standard sample.

Missing Peak(s)
An analyte does not elute inside the gradient or is not detected within the run time of the sample.

- Evaporative loss of analyte
- Degradation of analyte
- Sample preparation error
- Loss of detector sensitivity
- Co-elution with another analyte
- Wrong sample injected

Step 1

Benchmark Performance
Use a Quality Control Reference Material (QCRM) to **establish a benchmark** on your LC system when the system is in a **known good state**, such as after installation or performance maintenance and testing by a service engineer.

Step 3

List Failure Indicators
Identify which failure indicators are present in your qualification check chromatogram from the menu below. In the example to the left we have two failure indicators:

- Missing Peak (Peak 5)
- Decreased Peak Area (Peak 7)

Step 2

Qualification Check
Perform periodic qualification checks to confirm **continued good system performance**. A failed qualification check should trigger troubleshooting to identify the cause.

For best troubleshooting results, use a qualification standard that contains a variety of analyte types, such as Waters Reversed Phase QCRM, P/N 186006363.

Step 4

List Analyte Properties and Pass/Fail Category
Identify the properties of your peak analytes and list whether they pass or fail the qualification check:


- Uracil – polar, neutral, PASS
- Propranolol – basic, PASS
- Amitriptyline – basic, PASS
- Butylparaben – neutral, PASS
- Naphthalene – neutral, volatile, **FAIL**
- Dipropyl phthalate – neutral, PASS
- Acenaphthene – neutral, volatile, **FAIL**

Step 7

Test Possible Failure Causes
Once you have narrowed down the list of possible failure causes, test the **most likely** or **easiest to test** causes first. If the problem is unresolved, move on to the next most likely cause, and continue until your problem is resolved. In this example, it is very easy to test the two remaining possible failure causes at the same time by using a **new vial** of our reference standard.

Step 8

Re-Run Qualification Check
Confirm that the problem has been resolved. In this example, we can conclude that the failure cause was "Evaporative loss of analytes" because we can see that both of our failed peaks are volatile compounds, while the non-volatile compounds remained at a similar concentration to benchmark, and the problem was resolved when we injected from a new vial.



Failure Indicators Menu

Failure Indicators represent how the chromatogram failed and appear in your chromatograms as changes from benchmark values in six different categories. Each failure indicator has a list of possible Failure Causes, which represent why the chromatogram failed. The failure cause is the source of the failure indicator, and is typically instrument, column, method, or analyst related.

Change in Peak Retention Time	Increase in Retention Time <i>The new retention time is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Leak in systemOrganic solvent not being pumped properly (pump plungers, seals or check valves degraded)Column not properly equilibrated (either to starting mobile phase conditions or column temperature)Change in mobile phase pHWrong solvent/method being usedBubble in fluidic linesColumn degradation	Decrease in Retention Time <i>The new retention time is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Organic solvent not being pumped properly (pump plungers, seals or check valves degraded)Column not properly equilibrated (either to starting mobile phase conditions or column temperature)Change in mobile phase pHSample diluent effectWrong solvent/method being usedColumn overloadColumn degradation	Erratic Retention Time <i>Unstable retention time across multiple injections, with no trend.</i> <ul style="list-style-type: none">Solvent delivery failureGradient mixing failureColumn temperature variabilityMobile phase additive variationsMobile phase line partially blockedColumn degradation	Drifting Retention Time <i>A trend of increasing or decreasing retention time across multiple injections.</i> <ul style="list-style-type: none">Column not properly equilibrated (either to starting mobile phase conditions or column temperature)Solvent delivery failureGradient mixing failureDegradation of mobile phase additivesMobile phase composition changeMobile phase not prepared properlyColumn degradation
Change in Peak Area	Drifting Peak Area <i>A trend of increasing or decreasing peak area across multiple injections.</i> <ul style="list-style-type: none">Drifting decreaseEvaporative loss of analyteDegradation of analyteSample temperature not equilibratedDrifting increaseEvaporative loss of sample diluentSample temperature not equilibratedNon-specific adsorption	Erratic Peak Area <i>The new retention time is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Leak in sample fluidicsSample syringe plunger degradedSyringe damagedBubble in sample fluidicsFluidic volumes not characterized/calibratedLeak in systemLoss of detector sensitivityColumn degradation/failure	Decrease in Peak Area <i>The new peak area is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Evaporative loss of analytesDegradation of analytesLeak in system (between sample injector and detector)Change in split factor (for systems with a splitter only)Sample preparation errorSample volume low in the vialWeak needle wash empty or lowLeak in sample fluidicsBubble in sample fluidics lineWeak needle wash not compatible with sampleLoss of detector sensitivityInjection needle damagedInjection needle not drawing enough liquidNon-specific adsorption	Increase in Peak Area <i>The new peak area is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Concentration of analytesEvaporative loss of sample diluentWrong method/injection volume usedChange in split factor (for systems with a splitter only)Sample preparation errorCo-elution of peaks
Change in Peak Shape	Peak Tailing <i>Peak symmetry value is increased compared to the benchmark data.</i> <ul style="list-style-type: none">Secondary interaction with stationary phaseChange in mobile phase pHSystem tubing not seated properlySlight overload of analyteSample diluent effectsColumn degradationNon-specific adsorption	Peak Fronting <i>Peak symmetry value is decreased compared to the benchmark data.</i> <ul style="list-style-type: none">Secondary interaction with stationary phaseChange in mobile phase pHSystem tubing not seated properlyColumn overloadSample diluent effectsColumn degradationChange in detector settings	Peak Splitting <i>A single analyte yields two or more peaks that share a base.</i> <ul style="list-style-type: none">Sample diluent effectsParticulates accumulating on inlet fritSystem tubing not seated properlySample volume too largeChange in mobile phase pH (mild splitting)Column degradation	Peak Shouldering <i>A single analyte yields two partially separated peaks that are not as well separated as a split peak.</i> <ul style="list-style-type: none">Sample diluent effectsParticulates accumulating on inlet fritColumn overloadWeak needle wash not compatible with separationSystem tubing not seated properlyChange in mobile phase pHColumn degradationNon-specific adsorption
Change in System Pressure	Increase in System Pressure <i>The new pressure is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Particulates accumulating in column fritsSample precipitation onto columnClog in system tubingWrong method being usedChange in column temperatureColumn degradation	Drifting System Pressure <i>A steady increase or decrease in system pressure across replicate injections with a consistent trend.</i> <ul style="list-style-type: none">Particulates accumulating in column inlet fritSample precipitation onto columnColumn degradation	Decrease in System Pressure <i>The new pressure is stable but does not match the benchmark data.</i> <ul style="list-style-type: none">Leak in systemMobile phase lines obstructedCheck valve degradationAir in solvent linesChange in column temperatureColumn degradation	Erratic System Pressure <i>The system does not maintain steady pressure across replicate injections/within a single injection.</i> <ul style="list-style-type: none">Plunger/seal degradationCheck valve problemSolvent inlet filter cloggedSolvent mixer failureColumn degradationAir in solvent lines
Change in Peak Width	Peak Narrowing <i>Peaks are symmetrical, but narrower compared to the benchmark.</i> <ul style="list-style-type: none">Change in system tubingDecrease in system bandspreadingChange in mobile phase pHChange in detector settings	Peak Widening <i>Peaks are symmetrical, but wider compared to the benchmark.</i> <ul style="list-style-type: none">Change in system tubingIncrease in system bandspreadingSystem tubing not seated properlyChange in mobile phase pHSample diluent effectsVolume overloadChange in detector settingsColumn degradation	Extra Peak <i>There is an extra peak in the separation.</i> <ul style="list-style-type: none">Carryover from previous injectionSample preparation failureAnalyte separating as two peaks (see Peak Splitting)Contaminated mobile phaseContaminated system	Missing Peak(s) <i>An analyte does not elute inside the gradient or is not detected within the run time of the sample.</i> <ul style="list-style-type: none">Evaporative loss of analyteDegradation of analyteSample preparation failureLoss of detector sensitivityCo-elution with another analyteWrong sample injectedNon-Specific adsorption