


WATERS OASIS μ ELUTION PLATES

- Elution volume as low as 25 μ L
- No evaporation and reconstitution
- Ideal for small sample volumes
- Up to a 15x increase in concentration



PATENTED INNOVATION

Now you can confidently perform SPE cleanup and analyte enrichment of very small sample volumes (10-25 μ L) up to a maximum of 375 μ L. The Waters Oasis® μ Elution plate combines patented* plate design, proven Oasis chemistries**, and straightforward protocols that deliver high analyte recovery and clean extracts in elution volumes as low as 25 μ L. Using the Oasis μ Elution plate achieves superior results compared to other conventional SPE formats in less time. This plate format produces concentrated extracts that can be directly injected into your LC/MS/MS, eliminating the need for the time-consuming evaporation step.

Eluting in 25 μ L without evaporation provides up to a 15-fold increase in analyte concentration, enabling sensitive, robust, and reproducible SPE results. Scientists in both pharmaceutical drug discovery and drug development can prepare the cleanest biological sample extracts for more sensitive LC/MS/MS analysis.

*U.S. Patent 6,723,236


** U.S. patents 5,882,521 (1996), 5,976,376 (1998), 6,106,721 (1999), 6,254,780 (2001), 6,322,695 (2001), 6,468,422 (2002), 6,726,842 (2004), 6,773,583 (2004), 6,723,236 (2004), additional patents pending.

RAPID AND SENSITIVE METHOD FOR THE DETERMINATION OF ROPINIROLE IN HUMAN PLASMA

Ropinirole is a potent drug marketed under the name Requip® in the U.S. and Adartrel® in Europe. It is prescribed for both restless leg syndrome and Parkinson's disease. A bioanalytical method was developed and validated using Oasis MCX μ Elution plates for sample preparation and concentration, and an ACQUITY UPLC®/MS/MS for highest sensitivity. The required LLOQ (lower limit of quantitation) of 0.005 ng/mL (5 pg/mL) for ropinirole in human plasma was easily achieved. The response at 0.005 ng/mL is more than 10x the level found in the blank sample, surpassing the FDA guidelines for determining the limit of quantitation in bioanalytical methods.

This example dramatically demonstrates both the powerful extraction capability and capacity of Oasis MCX mixed-mode cation-exchange sorbent in the μ Elution plate format and the high sensitivity, speed, and specificity capabilities of the combination of Waters ACQUITY UPLC and Quattro Premier™ XE MS systems.


Oasis MCX μ Elution 96-well Plate Procedure




ACQUITY UPLC Conditions

Column:	ACQUITY UPLC BEH C ₁₈ ⁺ 2.1 x 50 mm, 1.7 μ m			
Column part number:	186002350			
Mobile phase A:	10 mM NH ₄ COOH, pH 9			
Mobile phase B:	CH ₃ OH			
Flow rate:	0.5 mL/min			
Gradient:	Time (min)	Profile A (%)	Curve B (%)	
	0.0	95	5	6
	2.0	2	98	6
	2.5	2	98	6
	2.6	95	5	6
	3.0	95	5	6
Injection volume:	8 μ L			
Column temperature:	45 °C			
Sample temperature:	15 °C			
Sample diluent:	CH ₃ OH + 5% NH ₄ OH			
Strong needle wash:	60:40 CH ₃ CN:IPA + 0.5% HCOOH (1200 μ L)			
Weak needle wash:	95:5 H ₂ O:CH ₃ OH (500 μ L)			
96-well collection plate part number:	186002481			

UPLC®/MS/MS Chromatogram of the XIC for a Sample of Ropinirole in Human Plasma at LLOQ of 0.005 ng/mL




Ropinirole Standard Curve and Linearity Data



Ropinirole standard linear curve fit with $1/x^2$ weighting, from 0.005 ng/mL to 10 ng/mL over 3.5 orders of magnitude, Correlation Coefficient: $r = 0.999054$, $r^2 = 0.998108$, and Calibration Curve Equation: $2.20969 \cdot x + 0.00746661$

SAMPLE ENRICHMENT: UP TO A 15X CONCENTRATION FACTOR




STRAIGHTFORWARD METHODOLOGY

The Oasis 2x4 method is a simple, logical approach to the selection of an SPE sorbent and protocol. Two protocols and four sorbents provide the flexibility to extract acids, bases, and neutrals with high SPE recoveries while removing matrix components that may interfere with the analysis.


This flow chart outlines the simple steps required to achieve high recovery and the cleanest extracts:

- Characterize your analyte [neutral, acid or base; pK_a].
- Select one of the four Oasis sorbents.
- Apply the indicated protocol [1 or 2].
- Determine SPE recoveries by LC analysis.


Oasis 2x4 Method Optimized for μ Elution Plate



Oasis MCX μ Elution Plate Protocol



Oasis MCX 10-mg Plate Protocol




Those who have used the Oasis 2x4 method with cartridges will note that the elutropic strength of the final elution solvent in the above protocols has been increased to optimize recovery in minimal elution volumes. We recommend as a starting point, using a 60:40 ratio of acetonitrile to methanol, not 100% methanol; this provides the optimal elution strength, viscosity, and solubility properties appropriate for high analyte recovery of a diverse set of analytes in only 25 μ L.

OASIS 2X4 METHOD PROOF OF CONCEPT

Recovery Study

To demonstrate the logic, simplicity, and effectiveness of the Oasis 2x4 method, five samples of rat plasma were prepared, each spiked with one of the previously characterized test analytes shown below:

Proof of Concept: Analytes Spiked into Rat Plasma



Each plasma sample was diluted [1:1, v:v] and acidified with phosphoric acid [4% in water]. Respective aliquots were then processed using the protocol and the Oasis mixed-mode ion-exchange sorbent designated by the Oasis 2x4 method for the corresponding sample type. SPE recovery was determined by LC/MS/MS analysis. The neutral analyte was processed on all four sorbents used.


Matrix Effects Study

Matrix effects* are an alteration in MS response caused by interfering components in a sample. They may cause a loss of signal, ion suppression or gain in signal, ion enhancement. Phospholipids or lysophospholipids have been identified and reported as a major contributor to matrix effects in plasma samples.** Matrix effects decrease analytical method robustness and reproducibility, raise limits of detection/quantitation, and may lead to spurious results.

SPE cleanup using Oasis 2x4 method with Oasis μ Elution plates effectively reduces plasma interferences and their corresponding matrix effects.

Quantitative Assessment of Matrix Effects:

Matrix Effects Data for the Oasis 2x4 Method Proof of Concept



* Neue and McDonald, Waters Whitepaper


** Chambers et al., J. Chromatogr. B, 2007, Jun 1;852 (1-2): 22-34

μELUTION PLATE LOADING CAPACITY




SPE device capacity is defined as the total mass of analytes and endogenous sample components retained by the sorbent bed under loading conditions. Breakthrough will occur when the capacity of the sorbent bed is exceeded. The physicochemical properties of Oasis sorbents are designed to provide exceptionally high loading capacity, even though each well in a Waters Oasis μElution plate contains only 2 mg of Oasis sorbent.

To determine the Oasis μElution plate capacity, increasing volumes of plasma and urine samples (from 50 μL to 350 μL in 50 μL increments) were spiked with 200 ng/mL imipramine (non polar base) and 200 ng/mL atenolol (polar base). The plasma aliquots were diluted 1:1 with 4% aqueous H₃PO₄ and the urine aliquots were diluted 1:1 with H₂O and then loaded onto the μElution plate. SPE recovery was calculated and plotted for each loading level.

SPE Protocol for Oasis MCX μElution 96-well Plate



SPE Recovery for 200 ng/mL Imipramine and 200 ng/mL Atenolol on Oasis MCX μElution Plate



OASIS 96-WELL μ ELUTION PLATES (1/pkg)

Description	Part Number
Oasis HLB 30 μ m	186001828BA
Oasis MCX 30 μ m	186001830BA
Oasis MAX 30 μ m	186001829
Oasis WCX 30 μ m	186002499
Oasis WAX 30 μ m	186002500
Oasis Sorbent Selection Plate, 3 rows each MCX, WAX, MAX and WCX	186004475



Oasis μ Elution Plate

MANIFOLD FOR EXTRACTION PLATES

Description	Qty	Part Number
Extraction plate manifold for Oasis 96-well plates	1/box	186001831
Extraction plate manifold kit A (includes extraction plate manifold, reservoir tray, sealing cap and 350 μ L sample collection plate)		WAT097944
Disposable reservoir tray	25/box	WAT058942
Sample collection plate, 350 μ L	50/box	WAT058943
Sealing cap for 96-well collection plate	50/pkg	WAT058959
SPE vacuum pump 115 V, 60 Hz		725000417
SPE vacuum pump 240 V, 50 Hz		725000418
Vacuum box gasket kit Kit includes: 2 foam top gaskets 2 orange O-rings		186003522



SPE Vacuum Pump
(Includes two gauges and pressure regulator)



Vacuum Box Gasket Kit

Austria and European Export (Central South Eastern Europe, CIS and Middle East) 43 1 877 18 07, **Australia** 61 2 9933 1777, **Belgium** 32 2 726 1000, **Brazil** 55 11 4134 3788, **Canada** 1 800 252 4752 x2205, **China** 86 21 6879 5888, **CIS/Russia** 7 095 336 7000, **Czech Republic** 420 2 617 1 1384, **Denmark** 45 46 59 8080, **Finland** 358 9 5659 6288, **France** 33 1 30 48 72 00, **Germany** 49 6196 400600, **Hong Kong** 852 29 64 1800, **Hungary** 36 1 350 5086, **India and India Subcontinent** 91 80 2837 1900, **Ireland** 353 1 448 1500, **Italy** 39 02 265 0983, **Japan** 81 3 3471 7191, **Korea** 82 2 820 2700, **Mexico** 52 55 5200 1860, **The Netherlands** 31 76 508 7200, **Norway** 47 6 384 60 50, **Poland** 48 22 833 4400, **Puerto Rico** 1 787 747 8445, **Singapore** 65 6273 1221, **Spain** 34 93 600 9300, **Sweden** 46 8 555 11 500, **Switzerland** 41 56 676 70 00, **Taiwan** 886 2 2543 1898, **United Kingdom** 44 208 238 6100

All other countries: Waters Corporation U.S.A. 1 508 478 2000/1 800 252 4752



Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, The Science of What's Possible, Oasis, ACQUITY UPLC, UPLC and Quattro Premier are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2008 Waters Corporation. Printed in the U.S.A.
July 2008 720000467EN IH-DS

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com

