

HPLC to UPLC® Method Migration: An Overview of Key Considerations and Available Tools

Dr. Michael Swartz, Ph. D.
Principal Consulting Scientist
Worldwide Pharmaceutical Business Operations
Waters Corporation

PittCon 2007

We Appreciate Your Opinion

- Please Complete Evaluation Form
- Qualify for drawing to win a \$25 AMEX gift card (drawn at end of this session – must be present to win)

AND

- Qualify to win a “Gold Package” registration to Inform 2007 (\$1595 value; one winner per day; winner to be notified post-Pittcon)

Thank You!


- UPLC® Overview
- Principles of Method Migration
- Examples
- Software Tools Available
- Summary



- A New Class of Separation Science
 - Based on chromatography columns with very small particles
 - Based on instruments designed to take advantage of the small particles
- Provides Improved Resolution, Speed, and Sensitivity with no Compromises
- Suitable for Chromatographic Applications in General
 - Appropriate for developing new methods
 - Appropriate for improving existing methods

Smaller Particles The Enabler of Productivity

HPLC




Flow Rate [mL/min]:

ID = 1.0 mm	0.04	0.07	0.10	0.13	0.17	0.20	0.24
ID = 2.1 mm	0.15	0.3	0.45	0.6	0.75	0.9	1.05
ID = 4.6 mm	0.7	1.4	2.1	2.8	3.5	4.2	4.9



HPLC vs. UPLC[®]: Speed, Sensitivity and Resolution

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Faster, More Sensitive Methods


Faster, More Sensitive,
Higher Resolution Methods




- Why Convert HPLC Methods to UPLC?
 - Get faster results with more resolution
 - More information
 - More robust methods
 - Better situational response time (stat samples faster, research decisions with more information, process monitoring, product release)
 - More samples analyzed per system, per scientist

Increased Productivity

- The New Method Must Preserve
 - Complete resolution of all relevant analytes
 - Peak homogeneity/purity
 - Certainty of peak identification
 - Quantitative accuracy and precision




Method Migration Success?



Method Migration Success!



Waters

THE SCIENCE OF WHAT'S POSSIBLE.™





Caffeic Acid Derivatives in Echinacea Purpurea

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



8 Diuretics + impurity



- It is possible to transfer most HPLC methods to ACQUITY systems with attention to details
 - Instrument design and performance differences could result in result differences
 - Some adjustments may be required
- ***But:*** This is not taking advantage of ***UPLC!!***

- Method Transfer Between HPLC Systems
 - First step towards method conversion
 - HPLC to ACQUITY as HPLC
- Method Migration or Conversion
 - HPLC separation to an ACQUITY UPLC[®] separation
- Method Optimization
 - Separation goals
 - Method verification
- Method Development
 - Systematic approach

Keys for HPLC to UPLC[®] Method Migration Success


- Proper Column Choice
 - Chemistry
 - Dimensions
- Keep in Mind Instrument Differences
 - Gradient delay volume
 - Detector data rates
- Proper Geometric Scaling
- Optimize for UPLC
- Don't Forget LC Theory



Determination of Proper Column Chemistry

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



Determination of Proper Column Dimensions

- Internal diameter
 - Generally prefer 2.1 mm
 - Only use 1 mm for specific reason
 - Severely sample limited
 - Direct flow to mass spectrometer
- Length
 - If primary goal is SPEED
 - 50 mm length to start
 - If primary goal is RESOLUTION
 - 100 mm length to start



Instrument Differences: Compensating for System Volumes

- Compare system volumes
 - This volume should be converted to column volumes for the best comparison
- If target system gives **smaller** isocratic segment
 - ADD an initial hold to the gradient table to give the identical hold.
- If target system gives **larger** isocratic segment
 - No exact compensation is possible
 - Chromatographic effect of extra isocratic hold usually small

Method Migration Example

Original HPLC Method


Objective: Maintain resolution while increasing speed

Resolution (1,2) = 12

Resolution (2,3) = 28


Original Column: 4.6 x 150 mm, 5 μ m
@ 1.5 mL/min

Ambient temperature: 21- 22°C
Flow rate: 1.50 mL/min
Sample analytes: 1. Caffeine (100mg/mL),
2. Hydroquinidine (33mg/mL),
3. 3-Aminobenzophenone (39mg/mL)
Molecular weight(s): Less than 500
Sample diluent: DMSO
Injection: 10 μ L
Detection: 254nm
Mobile phase: A: 0.05% TFA in water
B: 0.05% TFA in acetonitrile



Method Migration


Column Comparison: Injection Volumes



$$20\mu\text{L injection}/2.49\text{mL} = 0.8\%$$

UPLC

2.1 x 50mm




$$20\mu\text{L injection}/0.19\text{mL} = 11\%$$

**Sample volume too large
for smaller column volume**

Method Migration

Column Comparison: Injection Volumes

Column transfer from 4.6mm to 2.1mm i.d.
No injection volume scaling



Method Migration Example

Injection Volume Considerations

- Geometrically scale injection volume to volume of column
- Capacity proportional to surface area and internal solvent volume
- Suggested minimum injection volume on the instrument is 0.5 – 1 μL
 - If calculated volume too small for injection, dilute 5 - 10x with initial strength mobile phase
 - Typically 5 μL maximum injection on 2.1 x 50 mm

Method Migration

Calculate Injection Volume

Target injection volume =

$$\text{Original injection volume} \times \frac{\text{Target Column Volume}}{\text{Original Column Volume}}$$

Scaling a **10 μ L** injection on 4.6 x 150mm to 2.1 x 50mm

$$10\mu\text{L} \times \frac{3.14 \times 1.1^2 \times 50}{3.14 \times 2.3^2 \times 150} =$$

$$10\mu\text{L} \times \frac{0.19}{2.49} = 10\mu\text{L} \times 0.076$$

$$= \mathbf{0.8\mu\text{L}}$$

- For Migration:
 - First, **adjust flow rate** proportional to column diameter squared for constant linear velocity (geometrically scaled)
 - Second, **adjust Gradient Table** to maintain the same number of column volumes of solvent through the target column
 - Finally, **adjust** flow rate (linear velocity) **for smaller particle**
 - Analyte molecular weight must be considered

Method Migration

Scale the Flow Rate to Column Geometry

Scaling a **1.5mL/min** flow rate on 4.6x150mm to 2.1x50mm*

$$\text{Target Flow Rate} = \text{Original Flow Rate} \times \frac{\pi \times r^2 \text{ of Target}}{\pi \times r^2 \text{ of Original}}$$

This reduces to:

$$\text{Target Flow Rate} = \text{Original Flow Rate} \times \frac{d^2_{\text{Target}}}{d^2_{\text{Original}}}$$

So:

$$1.5\text{mL/min.} \times \frac{2.1^2}{4.6^2} = 0.31\text{mL/min.}$$

*Note: this assumes same particle size

Method Migration Example

Gradient Profile

- Express gradient duration in percent change per column volume (cv) units
- Calculate each segment as a number of column volumes
- Calculate time required to deliver the same number of column volumes to the target column at the chosen flow rate

Original Gradient Profile

Gradient Step	Time Since Injection	Flow Rate	%A	%B	Curve
Initial	0	1.5	95	5	*
2	15	1.5	5	95	6
3	20	1.5	5	95	1
4	30	1.5	95	5	1

Gradient Segments

Express as Column Volumes

For 15 min at 1.5mL/min on a 4.6 x 150mm column

$$\text{Gradient Volume} = \text{Flow Rate} \times \text{Time} = 1.5\text{mL/min} \times 15\text{min} = 22.5\text{mL}$$

$$\text{Column Volume} = \pi \times r^2 \times L = 3.14 \times 0.23^2 \times 15.0 = 2.49\text{mL}$$

$$\text{Gradient Duration (cv)} = \frac{\text{Gradient Volume}}{\text{Column Volume}}$$

$$\text{Gradient Duration} = \frac{22.5\text{mL}}{2.49\text{mL}} = 9.03 \text{ cv}$$

Original Gradient Profile for Scaling

Step	Time Since Injection	Flow Rate	%A	% B	Curve	Segment Duration (min)	Segment Duration (cv)
Initial	0	1.5	95	5	*	0	0
2	15	1.5	5	95	6	15	9.03
3	20	1.5	5	95	1	5	3.01
4	30	1.5	95	5	1	10	6.02

Scaling Gradient Step Time

Maintain Duration (cv)

Original Step 2: 15 min @ 1.5 mL/min with duration of **9.03cv**

Calculate Target Step 2: (keeping duration @ **9.03cv**)

$$\text{Column Volume} = \pi \times r^2 \times L = 3.14 \times 0.105^2 \times 5.0 = 0.17\text{mL}$$

$$\text{Gradient Step Volume} = \text{Duration (cv)} \times \text{Target Column Volume}$$

$$= \mathbf{9.03cv} \times 0.17\text{mL} = 1.54\text{mL}$$

$$\text{Gradient Step Time} = \text{Gradient Step Volume} / \text{Flow Rate}$$

$$= 1.54\text{mL} / 0.31 \text{ mL/min} =$$

5 min

Scaled Gradient Profile


2.1x50mm Column

Adjust time for same number of column volumes
per gradient segment

Gradient Step	Time Since Injection	Flow Rate	% A	% B	Curve	Segment Duration (min)	Segment Duration (cv)
Initial	0	0.31	95	5	*	0	0
2	5	0.31	5	95	6	5.0	9.03
3	6.67	0.31	5	95	1	1.67	3.01
4	10	0.31	95	5	1	3.33	6.02

Smaller Particles The Enabler of Productivity

HPLC



Flow Rate [mL/min]:

ID = 1.0 mm	0.04	0.07	0.10	0.13	0.17	0.20	0.24
ID = 2.1 mm	0.15	0.3	0.45	0.6	0.75	0.9	1.05
ID = 4.6 mm	0.7	1.4	2.1	2.8	3.5	4.2	4.9

- Consider 1.7 μ m target particle (2.1mm ID column)
- Assume temperature and viscosity transferred
- Adjust flow rate based on van Deemter curve and approximate molecular weight
 - ~0.6 mL/min for smaller molecules
 - average 500 dalton (molecular weight) molecules
 - ~0.1 mL/min for larger molecules because diffusion is slower
 - e.g., ~2,000 dalton peptides

Scaling for UPLC® Flow Rate

Step Time to Maintain Duration (cv)

Original Step 2: 15 min. @ 1.5 mL/min with Duration of **9.03cv**

Calculate Target Step 2: (keeping duration @ **9.03cv**)

Target Column Volume (2.1 x 50) = 0.17mL

Gradient Step Volume = Duration (cv) x Target Column Volume

= **9.03cv** x 0.17mL = 1.54mL

Gradient Step Time = Gradient Step Volume / **UPLC™ Flow Rate**

= 1.54mL / **0.60 mL/min.** =

2.6min

Select UPLC® flow rate and adjust time to maintain same number of column volumes per segment

Gradient Step	Time Since Injection	Flow Rate	% A	% B	Curve	Segment Duration (min)	Segment Duration (cv)
Initial	0	0.6	95	5	*	0	0
2	2.61	0.6	5	95	6	2.61	9.03
3	3.48	0.6	5	95	1	0.87	3.01
4	5.22	0.6	95	5	1	1.74	6.02


Review steps

- Refer to current chromatography
 - Observe system volumes, solvents and detection technique
 - Define objectives/room for improvement
- Select column dimensions – scale flow for linear velocity
 - 50 mm length for speed
 - 100mm length for complex samples/resolution
- Scale injection volume to column dimensions
- Use a gradient and flow scaled from current method
 - Adjust gradient to accommodate system differences, column differences and particle size differences
- Use a gradient and flow rate scaled for UPLC®

Method Conversion: Original HPLC Method

Resolution (1,2) = 12
Resolution (2,3) = 28

Original Column: 4.6 x 150 mm, 5 μ m
@ 1.5 mL/min, 30 minute cycle time




Ambient temperature: 21- 22°C
Flow rate: 1.50 mL/min
Sample analytes: 1. Caffeine (100mg/mL),
2. Hydroquinidine (33mg/mL),
3. 3-Aminobenzophenone (39mg/mL)

Molecular weight(s): Less than 500
Sample diluent: DMSO
Injection: 10 μ L
Detection: 254nm
Mobile phase: A: Water
B: acetonitrile
C: 1.0% TFA in water

Gradient: (A/B/C) 90/5/5 to 0/95/5 in 15 minutes (Curve 6)
hold for 5 minutes then to initial conditions (curve 11)
hold for 10 minutes

Resolution 1,2 = 11
Resolution 2,3 = 26

ACQUITY UPLC BEH C18, 1.7 μm ,
2.1 x 50 mm @ 0.6mL/min




Method Conversion

Scaled UPLC[®] Magnified

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Resolution 1,2 = 11
Resolution 2,3 = 26

UPLC™: 2.1 x 50 mm column, 1.7 μm
@ 0.6 mL/min, 5.2 minute cycle time




Original HPLC Method

Critical resolution

Resolution 1,2=12

Resolution 2,3=28




UPLC® Magnified View

Critical resolution

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Resolution 1,2=11
Resolution 2,3=26



ACQUITY UPLC® Calculator

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Gradient Separations

File Run Help

Gradient Separations

Select your existing HPLC conditions here:

Column Length: 25 cm, 15 cm, 12.5 cm
 Column Diameter: 4.6 mm, 4.0 mm, 3.9 mm
 Particle Size: 10 µm, 8 µm, 5 µm
 MW of Sample: 100, 200, 300
 Flow Rate: 1.00 mL/min
 Temperature: 30 Celsius

Organic Modifier in B: MeCN 0%, MeOH 100%
 Pressure Units: bar, PSI, MPa

Results With Existing Column

Peak Capacity: 126
 Column Length: 15 cm
 Particle Size: 5 µm
 Pressure* [PSI]: 1227
 Column Diameter: 4.6 mm
 MW of Sample: 300

Geometrically Scaled UPLC(TM) Method Choices | **Optimally Scaled UPLC(TM) Method Choices** | Print

Scaled Gradient 2.1 mm Column
(HPLC Linear Velocity)

Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure* [PSI]	Injection Volume [µL]	Detailed Gradient Profile
5 cm	2.1	0.208	138	10.1	3158	1.4	View
10 cm	2.1	0.208	197	20.1	6315	2.8	View

Change Flow: Calculate HPLC Linear Velocity UPLC(TM) Linear Velocity

Scaled Gradient 1.0 mm Column
(HPLC Linear Velocity)

Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure* [PSI]	Injection Volume [µL]	Detailed Gradient Profile
5 cm	1.0	0.047	92	10.1	3158	0.3	View
10 cm	1.0	0.047	154	20.1	6315	0.6	View

Change Flow: Calculate HPLC Linear Velocity UPLC(TM) Linear Velocity

Geometrically Scaled UPLC(TM) Method Choices

Step	Time	Flow	%A	%B	Time Segment (min)
Init Cond.	0.00	1.00	100	0	0.00
Init Hold	10.00	1.00	84	16	10.00
3	20.00	1.00	64	36	10.00
4	30.00	1.00	15	85	10.00
5					
6					
7					
8					
9					
10					

Select Maximum Column Pressure 10000 PSI [Change](#)

[Click here for your ACQUITY UPLC™ column choices.](#) Calculate

Here are your UPLC™ methods. Select View to see method details.

start | Mi... | 2.1 | 2 M | A... | Gr... | 2:04 PM

Method Conversion



Significantly Different Column Chemistry

At what point does method conversion become method development?

- When column chemistries are different, try a modified method development approach
- Proposed Strategy 1:
 - Understand general column selectivity differences
 - Follow chemistry and system scaling protocols
 - Generic gradients : 10 to 90 % organic solvent
 - Run a short gradient (3 min) and a long gradient (6 min)
 - Optimization tools (simulation software) are useful
- Proposed Strategy 2:
 - Follow method development protocols

Automated Method Development

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Automated Method Development

- ACQUITY UPLC Column Manager, 4 column selection device
- ACQUITY UPLC Binary Solvent Manager, solvent select valves

- Four ACQUITY UPLC Chemistries 2.1 x 50 mm, 1.7 μm :
 - ACQUITY UPLC™ BEH C₁₈
 - ACQUITY UPLC™ BEH Shield RP₁₈
 - ACQUITY UPLC™ BEH C₈
 - ACQUITY UPLC™ BEH Phenyl
- Solvents:
 - Acetonitrile
 - Methanol
- Buffers:
 - pH 3 ammonium formate
 - pH 10 ammonium bicarbonate

Experimental Matrix





pH 3, Acetonitrile

pH 3, Methanol





pH 10, Acetonitrile

pH 10, Methanol





C18







Shield
RP18



C8



Phenyl



Automated Column and Solvent Scouting

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

BSMB006N_SMS105N_PDAD948M in McConville_Projects\Column Manager on CORE-SRV-09 as mcconvip/PowerUser - Editing SS Method: column scouting with ...

File Edit View Inject Actions Customize Help

Run Only Continue on Fault

Sample Set Method: column scouting with ACN

PlateWell	SampleName	Inj Vol (uL)	# of Inj	Function	Method Set / Report Method	Run Time (Minutes)	Next Inj Delay (Minutes)	An_Column	Gradient_Time	Column_Position
1				Equilibrate	2 to 35%ACN p750mL pH3 60V1	5.30	0.00			
2				Condition Column	2 to 35%ACN p750mL pH3 60V1	5.30	0.00			
3	1.A.2	70% excedrine std	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
4	1.A.4	salicylic acid	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
5	1.A.5	ASA	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
6	1.A.6	caffeine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
7	1.A.7	acetaminophen	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
8	1.A.8	excedrine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V1	5.30	0.00	ACQUITY 2.1 x 50 C18 1.7uM	4.000 column 1
9				Equilibrate	2 to 35%ACN p750mL pH3 60V2	5.30	0.00			
10				Condition Column	2 to 35%ACN p750mL pH3 60V2	5.30	0.00			
11	1.A.1	blk	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
12	1.A.2	70% excedrine std	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
13	1.A.4	salicylic acid	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
14	1.A.5	ASA	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
15	1.A.6	caffeine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
16	1.A.7	acetaminophen	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
17	1.A.8	excedrine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V2	5.30	0.00	ACQUITY 2.1 x 50 C8 1.7uM	4.000 column 2
18				Equilibrate	2 to 35%ACN p750mL pH3 60V3	5.30	0.00			
19				Condition Column	2 to 35%ACN p750mL pH3 60V3	5.30	0.00			
20	1.A.1	blk	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
21	1.A.2	70% excedrine std	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
22	1.A.4	salicylic acid	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
23	1.A.5	ASA	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
24	1.A.6	caffeine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
25	1.A.7	acetaminophen	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
26	1.A.8	excedrine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V3	5.30	0.00	ACQUITY 2.1 x 50 Phenyl 1.7uM	4.000 column 3
27				Equilibrate	2 to 35%ACN p750mL pH3 60V4	5.30	0.00			
28				Condition Column	2 to 35%ACN p750mL pH3 60V4	5.30	0.00			
29	1.A.1	blk	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
30	1.A.2	70% excedrine std	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
31	1.A.4	salicylic acid	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
32	1.A.5	ASA	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
33	1.A.6	caffeine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
34	1.A.7	acetaminophen	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4
35	1.A.8	excedrine	1.5	1	Inject Standards	2 to 35%ACN p750mL pH3 60V4	5.30	0.00	ACQUITY 2.1 x 50 Shield RP18 1.7uM	4.000 column 4

Single Samples / Sample Sets / Running /

ACQ B5MB1
ACQ B5MA2
Sample Set 0.00 0.01 0.02 0.03 0.04 0.05 Liters

Sample Set Time Remaining: 00:13:42
Total Samples Time Remaining: 04:16:33
New Sample Set Time: 04:02:50

ACQUITY PDA Detector
Running
Shutter: Open

ACQUITY Binary Solvent Manager
Running
7128 psi A2: 65.3%
0.750 mL/min B1: 34.7%

Instrument Method
2 to 35%MEOH_750mL_pH3_89V3_DL
Edit Monitor Setup

dry lab methanol #1:A,6 1 1.20

For Help, press F1

Sample Set - Injection Running

dry lab methanol

#1:A,6 1 1.20

Total Time and Solvents Calculated

BSMB006N_SMS105N_PDAD948M in McConville_Projects\Column Manager on CORE-SRV-09 as mcconvip/PowerUser - Run Samples

File Edit View Inject Actions Customize Help

Run Only Continue on Fault

Sample Set Name	# of Injections	Duration	Wait for User	Run Mode	Interactive Sys Suit	SS Method Name	Acquired by
dry lab methanol	8	8:28	<input type="checkbox"/>	Run Only	Continue on Fault	dry lab methanol	mcconvip
column scouting with ACN	27	242:55	<input type="checkbox"/>	Run Only	Continue on Fault	column scouting with ACN	mcconvip
column scouting with methanol	20	249:48	<input type="checkbox"/>	Run Only	Continue on Fault	column scouting with methanol	mcconvip

ACQ-BSMB1
ACQ-BSMA2
ACQ-BSMB2

Total Samples

Sample Set Time Remaining: 00:08:15
Total Samples Time Remaining: 08:20:17
New Sample Set Time: 00:00:00

ACQ-BSMB1
ACQ-BSMA2
ACQ-BSMB2

Total Sample

Sample Set Time Remaining: 00:08:15
Total Samples Time Remaining: 08:20:17
New Sample Set Time: 00:00:00

ACQUITY PDA Detector

ACQUITY Binary Solvent Manager



Running

5426 psi A2 97.3 %
0.750 mL/min B1 2.7 %

Instrument Method: 2 to 35%MeOH_750mL_pH3.89V3_DL

dry lab methanol

Can I transfer the method to other systems that have a standard column heater?



Develop Methods Faster with UPLC®:

Time Savings

UPLC Methods Development Protocol

2.1 x 50 mm, 1.7 µm

<u>pH 3/ acetonitrile</u>	<u>Time</u>
Flow ramp	5 min
Column conditioning (2 blank gradients)	11 min
Sample injection (2 replicates)	11 min

<u>pH 3/ methanol</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	11 min
Sample injection (2 replicates)	11 min
Column purge	6 min

<u>pH 10/ acetonitrile</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	11 min
Sample injection (2 replicates)	12 min

<u>pH 10/ methanol</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	11 min
Sample injection (2 replicates)	11 min
Column purge	6 min

120 min

SCREENING TIME **2 Hours/column**
x 4 columns

TOTAL SCREENING TIME **8 HOURS**

EQUIVALENT HPLC Methods Development Protocol

4.6 x 150 mm, 5 µm

<u>pH 3/ acetonitrile</u>	<u>Time</u>
Flow ramp	5 min
Column conditioning (2 blank gradients)	80 min
Sample injection (2 replicates)	80 min

<u>pH 3/ methanol</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	80 min
Sample injection (2 replicates)	80 min
Column purge	35 min

<u>pH 10/ acetonitrile</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	80 min
Sample injection (2 replicates)	80 min


<u>pH 10/ methanol</u>	
Flow ramp	5 min
Column conditioning (2 blank gradients)	80 min
Sample injection (2 replicates)	80 min
Column purge	35 min

730 min

SCREENING TIME **12.2 Hours/column**
x 4 columns

TOTAL SCREENING TIME **48.8 HOURS**

- Strategies and Tools For HPLC to UPLC® Method Migration
 - Column chemistry and dimensions
 - Instrument considerations
 - Proper scaling
 - Optimization
 - Redevelopment
 - Method migration/conversion calculator



- Diane Diehl
- Eric Grumbach
- Jeff Mazzeo
- Uwe Neue
- Michael Jones
- Andy Aubin
- Craig Dobbs

And:



- Waters premier Laboratory Informatics symposium where users can network with colleagues, share experiences, and gain valuable insights into how to best deploy and use Waters Informatics suite of solutions
 - May 7-10, 2007
 - Miami Marriott Biscayne Bay
- Three-day event with optional tutorials
 - Monday: (2) ½ Day Tutorial Session
 - Tuesday – Thursday: Symposium
- Download the agenda and register at:
 - Waters booth - 3641
 - www.inform2007.com

Fill out a seminar evaluation form and enter the raffle to win a Gold Package registration* for Inform 2007 which includes all tutorials and workshops - a \$1595 value

One drawing daily.
Winner will be notified.

* Prize covers conference registration only