

Simultaneous Determination of Vitamins A and E in Infant Formula by ACQUITY UPC² with PDA Detection

Jinchuan Yang

GOAL

To test the feasibility of simultaneous determination of vitamins A and E in infant formulas by UltraPerformance Convergence Chromatography™ (UPC²)

BACKGROUND

Determination of fat-soluble vitamins in food products typically involves lengthy sample preparation (such as saponification and extraction), followed by high performance liquid chromatography (HPLC) with UV or fluorescence detection.¹ Vitamins A and E are often determined in the same analysis because of their similar sample preparation requirements.¹ Recently, a simple sample preparation procedure without saponification was proposed in a simultaneous determination of vitamins A and E in infant formula (IF) method, where various forms of vitamins A and E were separated and quantified in a single injection.² The elimination of saponification greatly increased the throughput of the analysis; however, the chromatography time was still lengthy (25 minutes), and the resolution of *cis*- and *trans*- forms of vitamin A were not evaluated in the paper.

Waters® UltraPerformance Convergence Chromatography (UPC²) leverages the unique properties of supercritical carbon dioxide, including low viscosity, high diffusivity, and liquid-like solvation power. It provides an alternative approach to normal phase LC, reversed phase LC, and gas chromatography (GC) for a wide range of analytical challenges.^{3,4} To investigate the performance of UPC² for the simultaneous determination of vitamins A and E, a feasibility study on commercial IF samples was performed.

UPC² provides a faster, simpler, and “greener” approach for the simultaneous determination of vitamins A and E in infant formulas.

THE SOLUTION

In this feasibility study, vitamin A (retinyl acetate and retinyl palmitate) and vitamin E (alpha-tocopheryl acetate and alpha-tocopherol) were extracted from IF samples by a simple liquid-liquid extraction, then analyzed on an ACQUITY UPC²™ System coupled with an ACQUITY UPC² PDA Detector. The separation of the analytes was achieved on a single ACQUITY UPC² HSS C₁₈ SB 3.0 x 100 mm, 1.8 μm Column under a gradient elution of a carbon dioxide and methanol mixture (3% to 10% methanol). The chromatograms for these compounds were extracted from the PDA data at their maximum absorbance wavelength of 320 nm, 283 nm, 293 nm for retinol esters, alpha-tocopheryl acetate, and alpha-tocopherol, respectively, as shown in Figure 1.




Figure 1. Typical chromatograms of vitamins A and E by UPC² with PDA detection. (A) standards; (B) Infant formula sample. Peaks: 1 *cis*-retinyl acetate, 2 all *trans*-retinyl acetate, 3 *cis*-retinyl palmitate, 4 all *trans*-retinyl palmitate, 5 alpha-tocopheryl acetate, and 6 alpha-tocopherol.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

UPC² provides a fast and well-resolved separation of vitamins A and E. All of the compounds, including the *cis*- and *trans*- forms of the retinol esters, were separated from each other, and were eluted before sample matrix peaks. The total cycle time, including column equilibration, was eight minutes per injection, which was at least three times faster than the typical run time (25 minutes) using other methods. Details of the method's analytical performance in linearity, sensitivity, repeatability, and recovery are shown in Tables 1 and 2. Although the sample extract can be injected directly onto the system, the LOQ result indicated that evaporation or concentration of the extract may be needed for some compounds, depending on their content levels in IF samples. In this study, the sample extracts were concentrated ten-fold by evaporation for the vitamin E determination.

UPC² is an environmentally friendly, or "green" technology. The source of the primary mobile phase, carbon dioxide, is recaptured carbon dioxide that is released from other industries, so the use of carbon dioxide does not generate any new greenhouse gases. The co-solvent (methanol) consumption was only 0.9 mL for each injection by UPC² compared to 10 mL of hexane used in the method in reference 2. This corresponds to at least a 90% reduction in solvent consumption.

SUMMARY

Simultaneous determination of *cis*- and *trans*-retinyl palmitate, *cis*- and *trans*-retinyl acetate, alpha-tocopheryl acetate, and alpha-tocopherol in commercial IF samples was achieved in a single injection using a single UPC² column on Waters' ACQUITY UPC² System with PDA detection.

Compound	Range (µg/mL) ^a	R ²	Equation ^b	LOQ (µg/mL) ^a
Retinyl acetate	0.4 to 5.5	0.9977	Y=9186x + 1107	0.5
Retinyl palmitate	0.2 to 3.5	0.9953	Y=4212x + 98	0.3
Alpha-tocopheryl acetate	17.5 to 245.0	0.9988	Y=199x + 289	13.0
Alpha-tocopherol	15.8 to 221.0	0.9987	Y=315.7x + 80.7	12.0

Table 1. Linearity of the method and estimated LOQ by UPC²/PDA.

^a Expressed in µg of analyte per mL in standard solution.

^b Y, peak area; x, concentration (µg/mL).

Compound	Repeatability (n=6)		Recovery (%) (n=3)
	Mean ± SD (µg/g) ^a	RSD (%)	Mean ± SD
Retinyl acetate	5.34 ± 0.04	0.7	91 ± 0.8
Retinyl palmitate	13.6 ± 0.3	2.4	–
Alpha-tocopheryl acetate	130.0 ± 1.0	1.1	94 ± 1.4
Alpha-tocopherol	82.3 ± 3.1	3.8	–

Table 2. Repeatability and recovery result obtained on spiked infant formula samples.

^a The values are expressed in µg vitamin per gram of infant formula powder.

The sample analysis time took eight minutes, three times faster than a typical analysis time, and the solvent consumption for each injection was 0.9 mL, one tenth of that in a normal phase LC method. This approach provides promising results in resolution, linearity, sensitivity, precision, and accuracy. UPC² has great potential to become a practical solution for the routine determination of vitamins A and E in infant formula products.

References

- DeVries JW, Silvera KR. Determination of vitamins A (retinol) and E (alpha-tocopherol) in foods by liquid chromatography: collaborative study. *J. AOAC International*. 2002; 85: 424.
- Chavez-Servin JL, Castellote AI, Lopez-Sabater MC. Simultaneous analysis of Vitamins A and E in infant milk-based formulae by normal-phase high-performance liquid chromatography-diode array detection using a short narrow-bore column. *J. Chromatogr. A*. 2006; 1122: 138.
- Aubin A. [Analysis of Fat-Soluble Vitamin Capsules Using UltraPerformance Convergence Chromatography UPC² Waters Application Note 720004394en](#). 2012 June.
- [ACQUITY UPC² System brochure. 2012. p/n 720004225en](#). 2012 March.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters is a registered trademark of Waters Corporation. ACQUITY UPC², UPC², UltraPerformance Convergence Chromatography, and The Science of What's Possible are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2013 Waters Corporation. Produced in the U.S.A.
January 2013 720004538EN TC-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com

