

**High performance liquid
chromatography
for the isolation and purification of
oligosaccharides from glycoconjugates.**

Waters Glyco-Pak™ DEAE and N

Obtain high purity carbohydrates for NMR and MS.
Prepare carbohydrates for use as substrates or to
determine biological activity.

**Fast and easy separation of acidic and neutral
carbohydrates**


Polymeric columns permit recovery of carbohydrates
without interference and provide high stability for
consistent chromatography.

Non-destructive purification

Methods provide unaltered and intact carbohydrates
ideal for structural analysis or biological activity
elucidation.

Fetuin oligosaccharides

1. Monosialylated
2. Disialylated
3. Trisialylated
4. Tetrasialylated




Glyco-Pak chemistries for the isolation and purification of oligosaccharides

The oligosaccharides (polymers of 2–50 monosaccharides) derived from glycoconjugates are one of the most challenging groups of biological molecules to separate and purify due to the large number of possible structures and linkages. High purity oligosaccharides are required for structure determination by NMR or MS as well as for use as substrates to determine their biological activity (i.e., enzymology of oligosaccharide biosynthesis, glycoconjugate modification, etc.). Single mode chromatography frequently is inadequate to purify these compounds to homogeneity.

Waters has developed a multiple column approach for oligosaccharide purification using two new column chemistries. The Waters Glyco-Pak DEAE and Glyco-Pak N are highly stable polymeric columns. They have been specifically developed and tested to provide consistent chromatography and yield carbohydrates without interferences due to the stationary phase. The Glyco-Pak DEAE column is designed for the purification of acidic oligosaccharides. The Glyco-Pak N column is a hydrophilic polymer column with distinctly different chemistry from the traditional silica based amino columns. This column is for the purification of neutral oligosaccharides.

A multi-column separation approach for the purification of oligosaccharides.




Oligosaccharides can be released chemically or enzymatically from glycoconjugates. After the removal of the sugars, the mixture of oligosaccharides is initially fractionated by gradient elution on the Glyco-Pak DEAE column. On this column, the neutral oligosaccharides are not retained and elute in the column void volume. Acidic oligosaccharides are eluted from the column according to the number of sialic acids or degree of sulphation. The acidic fractions can be further purified by an isocratic run on the Glyco-Pak DEAE column. The neutral fraction is then run on the Glyco-Pak N column. For more complex neutral mixtures, a subfractionation may be performed using a variety of Waters column chemistries. In the approach outlined above, the Waters Carbohydrate Analysis column is used to fractionate the neutral carbohydrates. Fractions from the Waters Carbohydrate Analysis column are then further purified on the Glyco-Pak N column. Used with the Waters 600E System Controller, which features a multisolvent deliv-

ery system and unique dual piston design, highly reproducible gradients are achieved ensuring consistent chromatography.


When isolating oligosaccharides for analysis by NMR, MS or for use as substrates, a minimum of 100µg of material is generally required. It is important that the mobile phase and detection be non-destructive. Ultraviolet and Refractive Index detection meet these requirements (detection limits of about 750ng for UV and 100ng for RI). Waters offers a full range of detection options from single or multi-channel UV/Visible to full spectral analysis as well as high sensitivity RI. Many of these detectors come with additional capabilities such as baseline correction for subtracting gradient shifts or spectral libraries for directly comparing peaks to previously acquired spectra.

Ovomucoid oligosaccharides




The initial fraction of oligosaccharides derived from ovomucoid is performed on the Waters Glyco-Pak DEAE column. Fraction 1 contains neutral carbohydrates and fractions 2 and 3 contain acidic carbohydrates.

Purification of neutral oligosaccharides



The neutral oligosaccharides obtained from fraction 1 of the DEAE run are now fractionated on the Waters Carbohydrate Analysis column.

Final purification of the neutral oligosaccharides



Fraction C collected from the Waters Carbohydrate Analysis column is now fractionated on the Glyco-Pak N column. This column provides an additional mode of separation with selectivity that is distinctly different from traditional NH₂ columns.


The Waters HPLC systems for complex carbohydrate analysis are complete gradient HPLC systems, featuring multisolvent inlet capability which permits the blending of up to four solvents or buffers. The micro-processor based control unit of the 600E can be programmed to store up to 15 complete gradient protocols allowing for quick conversion between modes of chromatography when performing multi-column carbohydrate separations.

High resolution separation of acidic oligosaccharides with anion exchange HPLC.

Fractionation of fetuin oligosaccharides on the Waters Glyco-Pak DEAE column

Sample: Hydrazine released fetuin oligosaccharides
Column: Glyco-Pak™ DEAE
Eluent A: H₂O
Eluent B: 100mM NaCl
Gradient: 0-50% B, 30 minutes, curve 6
Flow Rate: 0.6ml/min
Detection: UV at 205nm with baseline subtraction

1. Monosialylated
2. Disialylated
3. Trisialylated
4. Tetrasialylated



Acidic oligosaccharides are purified using the Glyco-Pak DEAE column in an anion exchange mode. This separation is based on degree of sialylation. Fractions can be collected and further purified isocratically on the Glyco-Pak DEAE column.


Oligosaccharides released chemically or enzymatically can be purified on the Glyco-Pak columns. Samples such as the Genzyme/Waters high mannose oligosaccharide standard can be purified in a single step.

Neutral oligosaccharide separations through partition HPLC.

High mannose oligosaccharides

Sample: Genzyme Endo H released high mannose oligosaccharide
Column: Glyco-Pak™ N
Eluent: CH₃CN/H₂O 70:30 (v/v)
Flow Rate: 1.0ml/min
Detection: UV at 197nm
Injection: 100mcg of high mannose standard

1. Man₄Gn
2. Man₅Gn
3. Man₆Gn
4. Man₇Gn
5. Man₈Gn



The Genzyme/Waters high mannose oligosaccharide standard separated on the Glyco-Pak N column.

Neutral oligosaccharide standards


Sample: Mixed neutral oligosaccharide standards

Column: Glyco-Pak™ N

Eluent: CH₃CN/H₂O (68:32) (v/v)

Flow Rate: 1.0ml/min


Detection: UV at 200nm



Neutral oligosaccharides are separated on the Glyco-Pak N column with a mobile phase consisting of water and 65 to 75% acetonitrile. Structures as diverse as monosaccharides, Man 9 high mannose and multiantennary complex carbohydrates can be separated isocratically free of contamination that is found with silica based amino columns. These columns provide consistent chromatography and because of nondestructive detection and mobile phase, yields carbohydrates without altering their structure.


Separation of mixed neutral oligosaccharide standard on the Glyco-Pak N column.

Neutral oligosaccharide standard structures†



The above structures are of the neutral oligosaccharide standard separated on the Glyco-Pak N column.


NMR of peak 5 of neutral oligosaccharide standard



Oligosaccharides purified on the Glyco-Pak N column can be taken directly to NMR or mass spectrometry without additional purification steps to remove contaminants. The Glyco-Pak N column has been found to yield cleaner NMRs without loss of the carbohydrate due to additional sample clean up.


Structural reporter regions of 500-MHz proton NMR of peak 5 of neutral oligosaccharide standard, a biantennary structure with terminal galactose (from fibrinogen) after purification on the Glyco-Pak N column.

NMR of peak 6 of neutral oligosaccharide standard



Structural reporter regions of 500-MHz proton NMR of peak 6 of neutral oligosaccharide standard, (Man 9 from soybean agglutinin) after purification on the Glyco-Pak N column.

Complex carbohydrate separation systems




The Waters family of complex carbohydrate separation systems are complete gradient HPLC systems, featuring precise gradient control, variable volume injection, sensitive and stable absorbance detection, data reduction, and full automation. They are modular in design and can be configured for either high pressure or low pressure chromatography.

Ordering information

Glyco-Pak DEAE Column	7.5mm × 7.5cm	PN 26730
Glyco-Pak N Column	7.8mm × 30cm	PN 26725
Waters Carbohydrate Analysis Column	3.9mm × 30cm	PN 84038

Options

Waters 600E Gradient Controller	PN 62710
Waters 484 Tunable UV/Visible Absorbance Detector	PN 80560
Waters 490 Programmable Multiwavelength Detector	PN 46000
Waters 410 Differential Refractometer	PN 70390



Waters HPLC systems for complex carbohydrates are modular in design and can be optimally configured to the needs of your laboratory. As your needs expand, your Waters complex carbohydrate system can grow with you.

Waters Chromatography Division
Millipore Corporation
 34 Maple Street
 Milford, MA 01757
 Telephone (508) 478-2000

Waters and Glyco-Pak are trademarks of Millipore Corporation. N-Glycanase is a registered trademark of Genzyme Corporation.

U.S.A. Sales Offices

New England Waters Chromatography Division / Millipore Corporation / 34 Maple Street / Milford, MA 01757 / Tel. 508-478-2000

Chicago Waters Chromatography Division 1555 Mittel Boulevard / Wood Dale, IL 60191 / Tel. 312-766-6060

New Jersey Waters Chromatography Division 35 Airport Park / Suite 230 / Morristown, NJ 07960 / Tel. 201-285-1404

North Carolina Waters Chromatography Division / 117 Edinburgh South / Suite 206 Cary, NC 27511 / Tel. 919-469-2501

San Francisco Waters Chromatography Division c/o Millipore Corporation / 448 Grandview Drive / South San Francisco, CA 94080 Tel. 415-952-9200

Southwest Waters Chromatography Division 2626 South Loop West / Suite 170 Houston, TX 77054 / Tel. 713-665-7310

Washington Waters Chromatography Division 3702 Pender Drive / Fairfax, VA 22030 Tel. 703-273-4435

Overseas Offices

Australia Millipore Pty. Ltd. / Waters Chromatography Division / Locked Bag 18 87-89 Mars Road / Lane Cove, N.S.W. 2066 Tel. Inwats (008) 222-111 / In Sydney Area (02) 427-0511 / Telex 25934

Austria Millipore Ges.m.b.h. / Hietzinger Hauptstrasse 145 / A-1130 Wien Tel. (0222) 82 89 26 / Telex 131464 Miliv A

Belgium and Luxembourg Millipore-N.V. Waters Chromatography Division / Raketstraat 60 / Rue de la Fusée 60 / 1130 Bruxelles-Brussels Tel. (02) 242.1740 / Telex 22798 milwat b

Brazil Millipore Industria e Comercio Ltda. Waters Chromatography Division / Rua Prof. Campos de Oliveira, 430 / CEP 04675-Caixa Postal 19065 / Sao Paulo / Tel. (011) 548-7011 Telex (11) 25078

Canada Waters Chromatography Division 3688 Nashua Drive / Mississauga / Ontario L4V-1M5 / Tel. 416-677-0580 Toll-free 1-800-268-4881 / Telex 06-968513

Denmark Millipore A/S / Waters Chromatography Division / Roskildevej 342 2630 Tåstrup / Tel. (02) 52 88 11 Telex 33544

Finland Millipore Oy / Waters Chromatography Division / Ruukinkuja 1 / 02320 Espoo Tel. 90-801 90 77 / Telex 12-2123 milpor sf

France Millipore S.A. / Waters Chromatography Division / Boite Postale 307 / 78054 Saint Quentin / Yvelines Cédex / Tel. (1) 30.45.70.00 Telex 698 371 F

India Waters Instruments (India) Pvt. Ltd. 1st Floor Shankar House / 1 Rajmahal Vilas Extn. / Mekhri Circle / Bangalore 560 080 Tel. (812) 34-9203 / Telex 8452028 SWPL IN

Italy Millipore S.p.A. / Waters Chromatography Division / Via Achille Grande, 23/20090 Vimodrone MI/Tel. (2)250781 Telex 312284 milspai/Telefax (2)2650324

Italy Millipore S.p.A. / Waters Chromatography Division / Viale del Caravaggio 107 / 00147 Roma RM / Tel. (06) 51.32.452 / (06) 51.31.593 (06) 51.31.512 / Telefax (06) 51.26.267

Japan Nihon Millipore Ltd. / Waters Chromatography Division / No. 5 Koike Bldg. 3-12, Kita-Shinagawa 1-chome / Shinagawa-Ku, Tokyo 140 / Tel. 81-3-471-7191 Telex 781-02324161

Mexico Millipore S.A. de C.V. / Waters Chromatography Division / Avenida Ingenieros Militares 85 / Mexico 11230 D.F. Tel. 576-9688 / Telex 1777 442

The Netherlands Millipore B.V. / Waters Chromatography Division / Penningweg 33 Postbus 166 / 4870 AD Etten-Leur (N.B.) Tel. (0 16 08) 22000 / Telex 54 840 milwa nl

Norway Millipore AB / Waters Chromatography Division / Enebakkeveien 133 / 0680 Oslo 6 Tel. (02) 67 82 53

Peoples Republic of China (Beijing) Nihon Waters Service Center Rm. 638, Xiyuan Hotel Erligou Xijiao, Beijing Tel. 8022018 / Telex 22685

Peoples Republic of China Millipore China Ltd. Waters Chromatography Division / Hope Sea Industrial Centre / 26 Lam Hing Street / Suite 105 / 1st Floor / Kowloon Bay / Hong Kong / Tel. 3-796 3736 Telex 3-7967374

Puerto Rico Millipore Corporation Puerto Rico Branch / Waters Chromatography Division KM 7.7 / State Road 172 / Cidra, Puerto Rico 00639 / Tel. 809-747-8444 Telex 385 72 76

Spain Millipore Ibérica S.A. / Waters Chromatography Division / Avda del Llano Castellano, 13 / 28034 Madrid Tel. (91) 729.03.00 / Telex 23545 milli e

Spain Millipore Ibérica S.A. / Waters Chromatography Division / Entenza 28 08015 Barcelona / Tel. (93) 325.96.16 Telex 50524 wtrs e

Sweden Millipore AB / Waters Chromatography Division / Box 233 / Långebergsgatan 30 421 23 Västra Frölunda / Tel. (031) 28 98 60 Telex 21064 milpor s

Sweden Millipore AB / Waters Chromatography Division / Grahmsvägen 15 172 46 Sundbyberg / Tel. (08) 98 89 60

Switzerland Millipore AG / Steinackerstrasse 11 8302 Kloten / Tel. 41-181-41-363 / Telex 56067

The United Kingdom and Ireland Millipore (U.K.) Limited / Waters Chromatography Division 11-15 Peterborough Road / Harrow, Middlesex HA1 2YH / Tel. (01) 864-5499 Telex 24 191 milpor g

USSR, Middle East, Eastern Europe and Africa Millipore Ges.m.b.h. / Hietzinger Hauptstrasse 145 / A-1130 Wien / Tel. (0222) 82 89 26 Telex 131464 Miliv A

West Germany Millipore GmbH / Waters Chromatography Division / Hauptstrasse 71-79 6236 Eschborn / Tel. 06196-494-0 Telex 4072704 mil d

For all other countries: Waters Chromatography Division / Millipore Corporation / 34 Maple Street Milford, MA 01757 / USA / Tel. 508-478-2000 Telex 174166

Waters
 Division of MILLIPORE