

# Performance of an ultra low elution volume 96-well plate


Claude R. Mallet, Ziling Lu, Jeff R. Mazzeo, Uwe D. Neue  
Waters Corporation

PittCon 2003  
March 10-14 2003  
Orlando, Florida

- **Sample Process**
  - Faster methods development
  - Generic and simple extraction protocol
  - Increased sensitivity and selectivity
- **Instrumentation**
  - Faster analysis (1000-analysis-per-day barrier)

**Sample Preparation**

- Raw sample:**
- CACO<sub>2</sub>, microsomes, P450, hepatocytes ... etc
  - tissue, CSF, plasma, serum urine, tears ... etc
  - water, sediment, food ... etc

**Extracted sample  
For LC/MS/MS**

Ideally, the final sample should be as clean as a non extracted standard.

**Chromatography****Polarity:**

Silica- C<sub>18</sub>, C<sub>8</sub>, C<sub>4</sub>, C<sub>2</sub>  
Hybrid- C<sub>18</sub>, C<sub>8</sub>, C<sub>4</sub>, C<sub>2</sub>  
Polymer- C<sub>18</sub>, C<sub>8</sub>, C<sub>4</sub>, C<sub>2</sub>  
Embedded polar group  
Cyano, Phenyl

**Particle size:**

2.5, 3.5, 5 or 7  $\mu$ m

**Internal diameter:**

4.6, 3.9, 2.1, 1.0,  
0.32 mm and 75  $\mu$ m

**Length:**

150, 100, 50, 30, 20 mm

**Mass Spectrometry****Source:**

ESI  
APCI  
Nano-ESI

**Mass analyzers:**

magnetic sectors  
electric sectors  
time of flight  
quadrupole  
ion trap  
FT-ICR

“Many extractions practices are based on classical methodologies of liquid-liquid or liquid-solid extraction using practices which have not changed for the last one hundred years”

## Classical

Protein precipitation

Liquid-liquid extraction

Membrane extraction

Soxhlet

## New technologies

Solid phase extraction

Solid phase micro extraction

Accelerated fluid extraction

Supercritical fluid extraction

Microwave-assisted extraction

## DRUG DISCOVERY

→ Target ID, Lead Generation, Lead Optimization

**Matrix:** CaCO<sub>2</sub>, microsomes P450, hepatocytes, rat plasma

**Sample volume:** 50 – 100  $\mu$ L

**Linearity range:** 5 – 250 ng/mL

## DRUG DEVELOPMENT

→ Pre-Clinical, Clinical Trials, Pilot Plant

**Matrix:** plasma/serum from rat, rabbit, dog, monkey, human

**Sample volume:** 50  $\mu$ L up to 1 mL

**Linearity range:** 0.1 ng/mL – 500 ng/mL

## DRUG MANUFACTURE

→ QA/QC, Patent Protection

# 96-Well Plates and Cartridges



96-well plate

barrel

# Generic reversed phase SPE Method

**Prepare Sample Solution**

**Condition/Equilibrate**  
500  $\mu\text{L}$  methanol / 500  $\mu\text{L}$  water

**Typical values for a  
10 mg bed packing**

**Load**  
500  $\mu\text{L}$  spiked sample solution

**Wash**  
500  $\mu\text{L}$  5% methanol in water

**Elute**  
500  $\mu\text{L}$  methanol

**Note: For larger bed  
packing (e.g. 30 or 60  
mg), increase the  
condition, load, wash and  
elution volumes**

**Pre concentration step**

**Evaporate and Reconstitute**

# Choice of Sorbent Weight Based on Sample Size

Sorbent per Well	Maximum Mass Capacity	Typical Sample Volume	Typical Elution Volume
5 mg	0.15 to 1 mg	10 to 100 $\mu\text{L}$	$\leq 150 \mu\text{L}$
10 mg	0.35 to 2 mg	50 to 400 $\mu\text{L}$	$\leq 250 \mu\text{L}$
30 mg	1 to 5 mg	100 $\mu\text{L}$ to 1 mL	$> 400 \mu\text{L}$
60 mg	2 to 10 mg	200 $\mu\text{L}$ to 2 mL	$> 800 \mu\text{L}$

**With decreasing limit of quantification (LOQ's), it will be necessary to include an evaporation and reconstitution step, which is extremely time consuming.**




# Waters two-Stage Well Design in 96-Well Extraction Plate


Enabling technology for 96-well plates - Oasis<sup>®</sup> sorbent in the amount required to meet your capacity and elution volume needs.

Designed for Optimal Flow Property


5 mg




10 mg



30 mg




60 mg




Waters

# Oasis<sup>®</sup> $\mu$ Elution 96-well plate





# Hydrophilic-Lipophilic Balanced copolymer



N-Vinyl-Pyrrolidone




Di-Vinyl-Benzene

**Water Loving**  
Hydrophilic monomer

- Provides wetting properties
- No impact of sorbent drying

**Fat Loving**  
Lipophilic monomer

- Provides reversed phase properties for analyte retention



Strong cation-exchange mode

Best retention at least 2 pH units below  $\text{pK}_a$

Reversed-phase mode

Propranolol  
Basic  
drug

# Generic micro-elution reversed phase SPE Method

**Prepare Sample Solution**

**Condition/Equilibrate**  
200  $\mu$ L methanol / 200  $\mu$ L water

**Smaller condition  
volumes**

**Dilute plasma  
1:1 ratio with water**

**Load**  
100  $\mu$ L spiked sample solution

**Wash**  
200  $\mu$ L 5% methanol in water

**Elute**  
25  $\mu$ L ACN:IPA 40:60 + 2% FA

**HILIC**

**Pre-concentration  
step is avoided**

**Dilute with 50  $\mu$ l water**

**C<sub>18</sub>**

# Generic micro-elution mixed mode SPE Method

**Prepare Sample Solution**

**Condition/Equilibrate**  
200  $\mu$ L methanol / 200  $\mu$ L water

**Smaller condition  
volume**

**Dilute plasma  
1:1 ratio with water**

**Load**  
100  $\mu$ L spiked sample solution

**Wash 1**  
200  $\mu$ L Water + 2 % FA

**Locks basic drug  
on ion exchanger**


**Removes polar  
interferences**

**Wash 2**  
200  $\mu$ L MeOH


**Elute**  
25  $\mu$ L ACN:IPA 40:60 + 2%  $\text{NH}_4\text{OH}$

**Pre-concentration  
step is avoided**

**Dilute with 50  $\mu$ l water**



**Propranolol**  
**Antihypertensive**




**Amitriptyline**  
**Antidepressant**




MS: Micromass Quattro Ultima	LC: Waters Alliance® 2795
Ion source: ESI (+)	Flow rate: 0.2 mL/min
Source temperature: 150 °C	Mobile phase A: Water + 0.5 % NH <sub>4</sub> OH
Gas cell: 2.0 e-3 bar Argon	Mobile phase B: ACN + 0.5 % NH <sub>4</sub> OH
Desolvation temperature: 350 °C	Column: XTerra MS C18
Capillary voltage: 3.5 kV	2.1 x 30 mm, 3.5 µm
Drying gas flow: 500 L/hr	LC conditions: 5% to 95% in 1 min.
Cone gas flow: 50 L/hr	Column temperature: ambient
Cone voltage: 35 volts	




MRM transition: Metoclopramide (IS)	m/z 299.8 → 226.7
Propranolol	m/z 259.9 → 154.9
Amitriptyline	m/z 278.1 → 232.9

# Comparison of Ppt vs HLB vs MCX




# Comparison of Ppt vs HLB vs MCX






# LC/MS/MS analysis of Terfenadine and Metabolites on mixed mode




# Calibration curve of Terfenadine on mixed mode

<b>Conc. ng/mL N=6</b>	0.5	1.0	5.0	10.0	20.0	100.0	200.0
<b>Average</b>	0.50	0.96	4.99	10.11	19.26	103.25	197.18
<b>Standard Deviation</b>	0.008	0.048	0.25	0.51	0.38	3.68	1.92
<b>RSD %</b>	2	5	5	5	2	4	1




# Calibration curve of Terfenadine-Carboxylate on mixed mode

<b>Conc. ng/mL N=6</b>	0.5	1.0	5.0	10.0	20.0	100.0	200.0
<b>Average</b>	0.49	1.02	5.05	10.12	19.71	101.69	196.39
<b>Standard Deviation</b>	0.013	0.052	0.26	0.49	0.46	3.25	4.41
<b>RSD %</b>	3	5	5	5	2	3	2



# Calibration curve of Terfenadine-Alcohol on mixed mode

<b>Conc. ng/mL N=6</b>	0.5	1.0	5.0	10.0	20.0	100.0	200.0
<b>Average</b>	0.51	0.97	5.19	10.24	19.80	102.05	196.55
<b>Standard Deviation</b>	0.008	0.041	0.27	0.44	0.77	3.15	7.92
<b>RSD %</b>	2	4	5	4	4	3	4





# Long term stability of Terfenadine and Metabolite on mixed mode

Conc. (ng/mL)	0.5	1.0	5	10	20	100	200
<b>Terfenadine</b>							
<b>Day 1 (N=6) (C.V.)</b>	0.50 (3)	0.97 (8)	4.96 (4)	9.93 (2)	19.89 (4)	105.29 (4)	195.87 (2)
<b>Day 2 (N=6)</b>	0.51 (2)	0.92 (2)	5.29 (2)	10.51 (4)	19.41 (3)	100.56 (4)	200.70 (1)
<b>Day 3 (N=6)</b>	0.50 (3)	0.98 (5)	5.09 (6)	9.64 (5)	19.22 (3)	106.38 (3)	194.47 (2)
<b>Day 4 (N=6)</b>	0.51 (2)	0.89 (1)	5.21 (3)	9.94 (5)	19.20 (3)	107.04 (2)	197.1 (1)
<b>Terfenadine-carboxylate</b>							
<b>Day 1 (N=6) (C.V.)</b>	0.50 (2)	0.96 (5)	4.88 (6)	10.45 (5)	20.74 (6)	na	na
<b>Day 2 (N=6)</b>	0.50 (3)	1.01 (6)	4.79 (1)	9.68 (4)	21.11 (3)	97.12 (5)	204.06 (4)
<b>Day 3 (N=6)</b>	0.45 (3)	0.99 (5)	5.21 (4)	10.08 (5)	20.79 (4)	98.36 (4)	194.85 (4)
<b>Day 4 (N=6)</b>	0.50 (2)	1.00 (4)	5.10 (6)	10.35 (3)	20.34 (3)	96.02 (5)	192.70 (3)
<b>Terfenadine-alcohol</b>							
<b>Day 1 (N=6) (C.V.)</b>	0.52 (1)	0.92 (2)	4.99 (2)	10.30 (4)	20.53 (2)	96.74 (4)	195.23 (5)
<b>Day 2 (N=6)</b>	0.51 (1)	0.94 (2)	5.01 (4)	10.03 (4)	20.55 (3)	100.83 (3)	191.65 (3)
<b>Day 3 (N=6)</b>	0.50 (2)	0.97 (5)	4.99 (5)	9.88 (2)	20.27 (5)	99.34 (5)	204.10 (3)
<b>Day 4 (N=6)</b>	0.50 (2)	1.00 (5)	5.03 (4)	10.36 (4)	20.23 (6)	97.43 (7)	196.04 (3)

Wide volume range for loading (50 to 800  $\mu\text{L}$ )

Low elution volume (25 – 50  $\mu\text{L}$ )

Generic protocols for reversed phase and mixed mode

Pre-concentration factor up to 10x

Sub ng/mL quantitation limits

96-well plate format for high throughput

No evaporation and reconstitution step needed