

Quad LCMS QC Reference Material

CONTENTS

- I. INTRODUCTION
- II. STORAGE AND STABILITY
- III. QC REFERENCE MATERIAL USAGE
- IV. QC REFERENCE MATERIAL USAGE EXAMPLE
- V. SUMMARY
- VI. ORDERING INFORMATION

I. INTRODUCTION

Quality Control (QC) Reference Materials contain mixtures of standards specifically chosen to provide an easy and reliable way to monitor the performance of chromatographic systems.¹ Using a QC Reference Material, you can be assured that your instrument and columns are ready to analyze your samples. Regular use of QC Reference Material's provides an opportunity to benchmark your chromatographic systems and track performance over time, making it easier to proactively identify problems and resolve them faster.

The Quad LCMS QC Reference Material is an 8 component mix used to provide a comprehensive reference standard for use with Waters Quad detectors, like Xevo® TQD and SQD, with a wide variety of conditions and methods.

Individual Components in the Quad LCMS QC Reference Material Mix

Component	Empirical formula	Exact mass (as [M+H] +)	Exact mass (as [M+H] -)	Concentration for analysis (µg/mL)
Acetaminophen	C ₈ H ₉ NO ₂	152.0712		100
Sulfaguanidine	C ₇ H ₁₀ N ₄ O ₂ S	215.0603	213.0446	50
Sulfadimethoxine	C ₁₂ H ₁₄ N ₄ O ₄ S	311.0814	309.0658	10
Val-Tyr-Val	C ₁₉ H ₂₉ N ₃ O ₅	380.2185	378.2029	25
Verapamil	C ₂₇ H ₃₈ N ₂ O ₄	455.2910		6
Terfenadine	C ₃₂ H ₄₁ NO ₂	472.3216		6
Leucine-Enkephalin	C ₂₈ H ₃₇ N ₅ O ₇	556.2771	554.2615	25
Reserpine	C ₃₃ H ₄₀ N ₂ O ₉	609.2812		18

The compounds in this mix give a mixture of responses in:

- ESI (+/-)
- Covers a wide range of *m/z*
- Optimized concentration to provide a more equal response by component in ESI+ mode
- Provides a separation in a range of chromatographic conditions used to benchmark instrument performance

II. STORAGE AND STABILITY

The standard comes in a flame sealed amber ampoule. It contains 1 mL in a matrix of 19% acetonitrile/81% water/0.008% formic acid. Due to the higher organic concentration, it is recommended, upon opening the ampoule, to immediately transfer the contents into a vial and then close the vial with a solid septum cap. Use Waters TruView™ LCMS Certified Vials with solid septum cap for best results.

The standard mix is shipped at ambient temperatures. It is highly recommended that upon receipt, the standard should be refrigerated at 4 °C for short term and for long term stored frozen (-15 °C).

Note: If any undissolved material is visible inside the ampoule/ vial, sonicate the unopened ampoule/vial until the material is completely dissolved.

III. QC REFERENCE MATERIAL USAGE

A QC Reference Material is used to benchmark, qualify and troubleshoot a chromatographic system:

- Benchmarking is the acquisition of data, using a mixture of standards such as a QC Reference Material, on a properly functioning chromatographic system with one or more new columns to determine normal performance parameter values (e.g., peak retention times, widths, intensities, areas, and tailing/asymmetry). Before benchmarking, a chromatographic system is confirmed to be in a properly functioning state by successfully completing [1] component level calibrations for all instrument modules and then [2] system level calibrations for the entire chromatographic system.² System benchmarking is done by performing replicate analyses of the QC Reference Material and calculating the mean and standard deviation for the desired performance parameter values. Upper and lower control limits (UCL, LCL) are then set. After benchmarking, these control limits are used to determine success or failure of system qualification. Benchmarking is generally completed on each new chromatographic system after it is installed and then repeated on each existing chromatographic system after it is calibrated.⁴
- Qualification is the acquisition of data, using a mixture of standards such as a QC Reference Material, on a benchmarked chromatographic system to determine if that system is still functioning properly and therefore is still “qualified” to receive analytical samples. When the QC Reference Material performance parameter values measured during qualification are within the control limits determined during benchmarking, the system passes qualification and is deemed ready for use. Qualification is performed on each benchmarked chromatographic system at regular intervals and after any maintenance or repair. Setting a minimum qualification interval of “every workday morning” or “the first workday morning of the week” is strongly recommended. Labs that have a high sample load are encouraged to also do bracketing qualification runs. This is where qualification is done before and after each block of samples, especially during unmonitored runs (e.g., nights, weekends). When the “before” and “after” qualifications pass, there is high confidence that data is trustworthy for the bracketed block of samples. The data from benchmarking and subsequent qualification runs can also be entered into a control chart, allowing the analyst to evaluate the system performance over time.³
- Troubleshooting is the sequence of activities performed to ascertain and correct deviations in chromatographic system performance. A qualification failure should trigger troubleshooting. The changes in performance parameter values (e.g., retention time shifts, and peak shape/area changes), observed using a QC Reference Material can facilitate more rapid determination of the root cause of qualification failure. Use of a QC Reference Material can thereby guide and accelerate the troubleshooting effort. After successful troubleshooting, a repeat of the qualification should then pass unless significant changes were made during troubleshooting. In such cases, re-calibration and re-benchmarking are likely needed.

IV. QC REFERENCE MATERIAL USAGE EXAMPLE

System: ACQUITY UPLC® with TQD Mass Spectrometer
 Column: CORTECS® UPLC® C₁₈, 1.6 µm, 2.1 x 50 mm
 Column temp.: 40 °C
 Mobile phase A: 0.1% acid in water
 Mobile phase B: Acetonitrile
 Gradient:


Time/min	Flow rate mL/min	%A	%B
Initial	0.8	90	10
2.5	0.8	25	75

Injection volume: 1 µL
 Sample: Quad LCMS QC Reference Material
 (p/n: 186007362) Lot: W23121301
 MS conditions: Full scan (Continuum) between 145-615 m/z with a scan time of 0.1 seconds
 Capillary voltage: 3.8 kV
 Cone voltage: 38 V
 Extractor: 3
 RF lens: 0.0
 Source temp.: 120 °C
 Desolvation temp.: 450 °C
 Desolvation flow: 1000 L/hr
 Cone flow: 30 L/hr

CV-Cone Voltage, CE-Collision Energy, PI-Product Ions

Compound	Precursor	CV	PI 1	CE	PI 2	CE
Acetaminophen	152.06	35	110	15	93	25
Sulfaguanidine	211.1	25	156	15	92	25
Sulfadimethoxine	311	25	156	25	92	25
Val-Tyr-Val	380.2	25	136.1	25	235.2	15
Verapamil	455.4	35	165.2	35	303.3	25
Terfenadine	472.2	35	436.2	25	454.3	25
Leucine-Enkephalin	556.2	35	397.2	15	425.2	15
Reserpine	609.3	35	195	35	397.2	15

Standard Chromatogram of Quad LCMS Quality Control Reference Material



V. SUMMARY

Quality Control Reference Materials are specifically designed to provide a controlled, consistent, and reliable way of monitoring system performance. Regular use of QC Reference Materials to benchmark and qualify systems gives assurance that chromatographic results are high quality. Troubleshooting with QC Reference Materials also minimize system downtimes.

VI. ORDERING INFORMATION

Description	Part Number
Quad LCMS QC Reference Material	186007362
LCMS QC Reference Material	186006963
Neutrals QC Reference Material	186006360
Reversed-Phase QC Reference Material	186006363
Preparative Chromatography Mix Standard	186006703
HILIC QC Reference Material	186007226
QDa QC Reference Material	186007345
UPC ² QC Reference Material	186007950
AutoPurification™ System Standard	716000765

References

1. A "chromatographic system" consists of the installed instrument components including all detectors (PDA, ELS, MS, etc.) plus the installed columns.
2. All calibrations are done as specified in the instrument manufacturer's protocols. In some cases, there may be no system level calibration whereupon Benchmarking occurs directly after component level calibrations.
3. Benchmarking control limits are determined by each lab according to their data precision and reproducibility needs. Setting control limits of + 3 standard deviations about the mean value is a common practice. For more information on control limits and control charts, see Mullins, E.; *Introduction to Control Charts in the Analytical Laboratory, Analyst* 1994, 119, 369–375.
4. Consult the chromatographic system manufacturer's recommendations regarding the frequency of calibration.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Waters, The Science of What's Possible, Xevo, ACQUITY UPLC, UPLC, and CORTECS are registered trademarks of Waters Corporation. TruView and AutoPurification are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2015 Waters Corporation. Produced in the U.S.A. April 2015 720005288EN IH-PDF

Waters Corporation
 34 Maple Street
 Milford, MA 01757 U.S.A.
 T: 1 508 478 2000
 F: 1 508 872 1990
www.waters.com