

# **Simple Techniques for Improving the Isolation of Synthetic Peptides**

**Jo-Ann Jablonski  
Principal Scientist  
Waters Corporation**

# Agenda

- Background
- Techniques
  - Scaling a separation
  - Focusing the gradient
  - At-column dilution
  - Temperature control
- Practical example
- Summary


# Why peptides?

- Unique role
- Attractive as potential drug candidates
- Sources

# Synthetic Peptides

- Synthetic peptides
  - Amino acids linked using organic chemical reactions
  - Average: 10-30 residues
  - Large: 70-120 residues or more
- Impurities
  - Failure Sequences
  - Deletion Sequences
  - Damaged Sequences
  - Cleavage Adducts
- Purification requirements
  - 75-98% purity of final product

# Peptide Isolation Workflow



## Most Common Operating Conditions

- Purification using reversed-phase chromatography
- C18 columns
- Elution with gradient
- Monitor with UV detection at 214-220nm
- Mass-directed isolation

# Agenda

- Background
- Techniques
  - **Scaling a separation**
  - Focusing the gradient
  - At-column dilution
  - Temperature control
- Practical example
- Summary

# Scaling Equations

## Mass

$$M_2 = M_1 \times \frac{L_2}{L_1} \times \frac{d_2^2}{d_1^2} \quad \text{where:}$$

$M_1$  is Mass on Column1

$M_2$  is Mass on Column2

$L_1$  is Length of Column1

$L_2$  is Length of Column2

$d_1$  is Diameter of Column1

$d_2$  is Diameter of Column2

## Gradient

$$t_2 = t_1 \times \frac{L_2}{L_1} \times \frac{d_2^2}{d_1^2} \times \frac{F_1}{F_2} \quad \text{where:}$$

$t_1$  is Gradient Duration for Column1

$t_2$  is Gradient Duration for Column2

$L_1$  is Length of Column1

$L_2$  is Length of Column2

$d_1$  is Diameter of Column1

$d_2$  is Diameter of Column2

$F_1$  is Flow rate for column1

$F_2$  is Flow rate for column2

## Flow rate

$$F_2 = F_1 \times \frac{d_2^2}{d_1^2} \quad \text{where:}$$

$F_1$  is Flow Rate for Column1

$F_2$  is Flow Rate Column2

$d_1$  is Diameter of Column1

$d_2$  is Diameter of Column2

# Scaling a Method

## Analytical Method

Time	Flow	%A	%B
0.00	1.46	95	5
5.00	1.46	50	50
5.50	1.46	5	95
6.50	1.46	5	95
7.00	1.46	95	5
10.00	1.46	95	5

Column: 4.6 x 50 mm

System dwell volume: 0.71 mL

Column volume: 0.698 mL

System dwell volume (in column volumes): 1.02

## Preparative Method

Time	Flow	%A	%B
0.00	25	95	5
0.30	25	95	5
5.30	25	50	50
5.80	25	5	95
6.80	25	5	95
7.30	25	95	5
10.30	25	95	5

Column: 19 x 50 mm


System dwell volume: 4.75 mL

Column volume: 11.9 mL

System dwell volume (in column volumes): 0.40

Gradient offset: 0.30

# Gradient Scaling



# Estimated Peptide Mass Loading Capacity (mg)

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.®

Length (mm)	Diameter (mm)					
	2.1	4.6	10*	19*	30**	50
50	0.04-0.11	0.3-0.6	1.5-3.0	4-9	11-22	31-62
100	0.11-0.22	0.5-1.0	2.5-5.0	9-18	22-45	62-125
150	0.15-0.33	0.8-1.6	4.0-8.0	13-27	34-68	93-186
250	0.26-0.55	1.3-2.6	6.0-12.0	22-45	56-112	155-310

\*5 µm and 10 µm OBD Prep Columns; \*\* 10 µm OBD Prep Columns

\* [OBD™], Optimum Bed Density; US Patent # 7,399,410, UK Patent # GB2408469

Estimates of mass load are broad because capacity depends on solubility of the peptide in the mobile phase.

# Flow rate (mL/min) Injection Volume ( $\mu$ L)

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.<sup>®</sup>




	Diameter (mm)					
	2.1	4.6	10*	19*	30**	50
Flow rate	0.19-0.39	0.9-1.8	4.5-9.0	16-32	40-80	111-222
Inj Vol	4.3	20	100	350	880	2450

\*5  $\mu$ m and 10  $\mu$ m OBD Prep Columns; \*\* 10  $\mu$ m OBD Prep Columns

# Agenda

- Background
- Techniques
  - Scaling a separation
  - **Focusing the gradient**
  - At-column dilution
  - Temperature control
- Practical example
- Summary


# Gradient Types



# Developing a Focused Gradient

- Estimate volume of system and column
- Make pilot run of sample
- From elution time, system volume and gradient table, estimate % that elutes the peptide
- From gradient and column volume, estimate gradient slope in pilot run in % change per column volume
- Create shallow gradient segment from 5% below to 3% above estimated elution percentage
  
- Full Description in Application Note:  
Developing Focused Gradients for Isolation and Purification  
720002955en


# Crude Peptide Analysis



XBridge Peptide BEH C<sub>18</sub> 5µm Column, 4.6 x 100 mm

Gradient: 5-50% B in 10 minutes, 3.4% change per column volume, 5 µL  
Peptide eluted at 33% B.

# Focused Gradients



# Agenda

- Background
- Techniques
  - Scaling a separation
  - Focusing the gradient
  - **At-column dilution**
  - Temperature control
- Practical example
- Summary

# Dissolving Synthetic Peptides

- No universal technique or solvent
- Apparent insolubility may include synthesis by-products
- Useful choices
  - Water with 0.1-2.0% TFA
  - Wet with small amount of DMF, dilute with water to ~10% organic
  - 100% DMSO
  - 100% Hexafluoroisopropanol (HFIP)
  - 6-12M Guanidine-HCl (not for use with mass spectrometers)
- Always results in relatively large volume of relatively strong solvent

# Standard Gradient System

## Conventional Separations Simplified Flow Diagram



**System Gradient Pump**  
(Aqueous flow delivers sample plug to column.)

**20 ml/min**  
**95% Aqueous**  
**5% Organic**



**Auto Sampler Valve**


**C  
o  
l  
u  
m  
n**

*Configurations vary.*

Fraction Collector  
Detector(s)

# Standard Prep Separation With DMSO As Sample Solvent

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



# At-Column Dilution

## At-Column-Dilution Separations *Simplified Flow Diagram*



**Gradient Pump**



**Loading Pump  
(Organic flow)**

**19 ml/min  
100% Aqueous**

**Tee**

**20 ml/min  
95% Aqueous  
5% Organic**



**Auto Sampler Valve**

**1.0 ml/min  
100% Organic**

Waters Application Note : 71500078010rA

US Patent (Waters Corporation)  
6,790,361 B2






*Configurations  
vary.*


**Fraction Collector**

**Detector(s)**

# At-Column Dilution Prep Separation



# Conventional vs. At-Column Dilution



## At-Column Dilution



- Permits injection of large volumes of strong solvents
- Improves both mass capacity and resolution
- Increases system ruggedness
- Extends column life
  - Reduces introduction of particles
  - Reduces pressure shock from viscous solvent
- Waters Application Note :71500078010rA
- US Patent (Waters Corporation)
  - 6,790,361 B2

# Agenda

- Background
- Techniques
  - Scaling a separation
  - Focusing the gradient
  - At-column dilution
  - **Temperature control**
- Practical example
- Summary


# Effect of Temperature Crude Synthetic Peptide

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®




# Effect of Temperature Peptide Separation

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



# Effect of Temperature Peptide Separation

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



# Controlling Temperature in Purification

## Practical Considerations

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.®

- Large diameter columns cannot be effectively heated from the outside
- High flow rate separations occur at the temperature of the incoming solvent
- Temperature gradient within column
- Add a coil of stainless steel tubing to the inlet end of the column for pre-heating (e.g. 5-mL sample loop)
- Immerse column and pre-heating coil into a water bath

# Column and Pre-heater Loop Immersed in Water Bath at 60°


Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



# At-Column-Dilution: Sample Loading with Temperature Control

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®

Sample in 100% DMSO carried to  
Tee by organic solvent




# Agenda

- Background
- Techniques
  - Scaling a separation
  - Focusing the gradient
  - At-column dilution
  - Temperature control
- **Practical example**
- Summary

# Temperature Evaluation Crude Peptide

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



XBridge Peptide BEH C<sub>18</sub> 4.6 x 50 mm

Gradient: 0-30% isopropanol in 5 minutes @ 40 °C; 5.13 mg/mL, 7 µL

# Focused Gradient

## Screening Gradient

Time	Flow	%A	%B
0.00	1.46	100	0
0.50	1.46	100	0
5.28	1.46	70	30
6.00	1.46	10	90
7.00	1.46	10	90
7.50	1.46	100	0
10.50	1.46	100	0

System volume: 0.77 mL

Column volume: 0.698 mL

Calculated peptide elution %: 10.85%

Screening gradient slope: 3.0%/column volume

Focused gradient segment range: 5-13% B

Calculated focused gradient slope: 0.6%/column volume

Calculated time for the focused gradient segment: 6.37 min

## Focused Gradient

Time	Flow	%A	%B
0.00	1.46	100	0
0.50	1.46	95	5
6.87	1.46	87	13
7.00	1.46	10	90
8.00	1.46	10	90
8.10	1.46	100	0
11.10	1.46	100	0


Gradient segment range: 5-13% B

Gradient slope: 0.6%/column volume

Time for the focused gradient segment: 6.37 min

# Loading Study Focused Gradient

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.<sup>®</sup>



5-13% isopropanol with 0.1% TFA @ 40 °C, 214 nm

# Scaling

## Analytical Method

Time	Flow	%A	%B
0.00	1.46	100	0
0.50	1.46	95	5
6.87	1.46	87	13
7.00	1.46	10	90
8.00	1.46	10	90
8.10	1.46	100	0
11.10	1.46	100	0

Column: 4.6 x 50 mm

System dwell volume: 0.77 mL

Column volume: 0.698 mL

System dwell volume (in column volumes): 1.10

## Preparative Method

Time	Flow	%A	%B
0.00	25	100	0
0.66	25	100	0
1.66	25	95	5
14.40	25	87	13
14.66	25	10	90
16.66	25	10	90
16.86	25	100	0
22.86	25	100	0

Column: 19 x 100 mm

System dwell volume: 9.75 mL

Column volume: 23.804 mL


System dwell volume (in column volumes): 0.41

Gradient offset: 0.66

# Preparative Chromatography


## Conventional Injection

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®






XBridge Peptide BEH C<sub>18</sub> 5 µm OBD Prep, 19 x 100 mm, 2.4 mg

# Fraction Analysis



# Fraction Analysis



# Comparison of Prep Gradient Methods

## Conventional Injection

Time	Flow	%A	%B
0.00	25	100	0
0.66	25	100	0
1.66	25	95	5
14.40	25	87	13
14.66	25	10	90
16.66	25	100	90
16.86	25	100	0
22.86	25	100	0

Gradient time: 14.40 - 1.66 = 12.74 min

Gradient volume: 12.74 min x 25 mL/min = 318.5 mL

Column volumes: 318.5 mL x 1 cv/23.804 mL = 13.38 cv

Gradient slope: 8%/13.38 cv = 0.6%/col vol

## At-Column Dilution Injection

Time	Flow	%A	%B
0.00	23.75	100	0
2.26	23.75	100	0
3.26	23.75	100	0
16.00	23.75	92	8
16.26	23.75	15	85
18.26	23.75	15	85
18.46	23.75	100	0
24.46	23.75	100	0

Gradient time: 16.00 - 3.26 = 12.74 min

Gradient volume: 12.74 min x 25 mL/min = 318.5 mL

Column volumes: 318.5 mL x 1 cv/23.804 mL = 13.38 cv

Gradient slope: 8%/13.38 cv = 0.6%/col vol

Loading pump flow rate during gradient = 1.25 mL/min


1.25 mL/min = 5% of the total flow

Total gradient method flow rate = 25 mL/min


# At-Column Dilution with Temperature Control

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®

Sample in 100% DMSO carried to  
Tee by organic solvent




# Loading Method Comparison











# Preparative Chromatography At-Column Dilution

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



# Fraction Analysis



# Agenda

- Background
- Techniques
  - Scaling a separation
  - Focusing the gradient
  - At-column dilution
  - Temperature control
- Practical example
- **Summary**

## Summary

- Focusing the gradient improves resolution and load
- Temperature control is a useful tool for improving peptide detection, peak shape, and resolution
- At-Column Dilution permits injection of large volumes of strong solvents, improves mass capacity and resolution, increases system ruggedness, and extends column life

## Supplementary Information

- Developing Focused Gradients for Isolation and Purification
  - Waters Technical Note, 720002955EN
- Effective Use of Temperature Control in Compound Isolation
  - Waters Technical Note, 720002954EN
- At-Column Dilution Application Notes
  - Waters Application Notes, 71500078010
- Prep 150 LC System: Considerations for Analytical to Preparative Scaling
  - Waters Application Note, 720005458EN
- Peptide Isolation Using the Prep 150 LC System
  - Waters Application Note, 720005455EN
- Topics in Liquid Chromatography Part 2, Optimum Bed Density [OBD™] Columns
  - Waters White Paper, 720001939EN